



博登科技有限公司
B&D Technologies Co.,Ltd

Dental Burs 牙科车针、磨头

Vol. 8 : 2024

About Us

B&D Technologies Company Limited are specialized in Rotary things of dental, Jewelry, Beauty or Hobbies & DIY (Do it Yourself) .

Our products include Diamond Burs and discs, Tungsten Carbide Burs or cutters, Flexible Rubber Burs, Wheels or Points, Polishing Brushes, Mounted Stones, CAD/CAM Milling Burs or any Diamond Stones or Polishers.


Our Rotatory Instruments are widely applied to Cut, Grind, Finish, Polish any materials of Zirconia, Lithium Disilicates, Porcelain (Ceramics) , Alloy or Metal, Acrylic (Resin) , PEEK, Composites etc . They can match any Handpiece of FG (Friction Grip) , RA/CA(Right Angle), HP conveniently and precisely .

We have obtained the certificates of ISO 13485 and CE since the year of 2015 . We provide innovative and professional products and service.

Whoever you are dentist or technician, we always make you relaxing, safe and efficient. We are your best friend in dentistry, Jewellery or Hobbies.



ISO 13485



CE 0197



Established in 2007



Clinic Part
临床类

Page 3 - 28 第3至28页



Diamond burs 金刚石车针

Page 4 - 12 第4至12页

Diamond Ortho-strips
金刚石锯弓

Page 13 第13页

Diamond Finishing Strips
金刚石抛光条

Page 14 第14页

Tungsten Carbide burs
钨钢车针

Page 15 - 18 第15至18页

White Stones 白车石
Green Stones 绿车石

Page 19 第19页

Diamond Stones, Polishers
陶瓷（橡胶）金刚车石磨头

Page 20 第20页

**Finishing or Polishing Burs
after Debonding**
（脱胶后）洁牙磨头

Page 21 第21页

Twin-Color Silicon Polishers
硅胶抛光磨头

Page 22 第22页

Silicon Polishers
硅胶抛光磨头

Page 23 第23页

Polishing Brushes 抛光刷

Page 24 第24页

Mandrels, Sterile Blocks
夹针、消毒盒

Page 25 第25页

Kit 套装

Page 26 - 31 第26至31页

Samples Collection
样品册


Page 32 第32页

Dental Diamond Bur

牙科金刚砂车针

Advantages

1. Strictly comply with ISO Standard ;
2. Sharper and More Efficient ;
3. Durable and Valuable ;
4. Full Range of shapes and grits.
5. Fastly develop new items.
6. Individual Design and OEM Service



Grit

- Black / Super Coarse
- Green / Coarse
- Blue / Medium
- Red / Fine
- Yellow / Extra Fine

0.8 19.0 001(801) BR-49	1.0 19.0 001(801) BR-45	1.2 19.0 001(801) BR-46	1.4 19.0 001(801) BR-41	1.6 19.0 001(801) BR-40	1.8 19.0 001(801) BR-31	2.1 19.0 001(801) BR-28	2.3 19.0 001(801) BR-26	2.5 19.0 001(801) BR-25	2.9 19.0 001(801) BR-23
001-008M	001-010M	001-012M	001-014M	001-016M	001-018M	001-021M	001-023M	001-025M	001-029M
圆球型 Ball, Round /BR									

6mm 1.2 22.0 697(801L) BR-L46	1.4 22.0 697(801L) BR-L41	1.6 22.0 697(801L) BR-L40	1.8 22.0 697(801L) BR-L31	1.2/2.3 19.0 002(802) BC-43	1.4/2.3 19.0 002(802) BC-42	1.6/3.5 19.0 002(802) BC-32	1.8/2.3 19.0 002(802) BC-31	1.2/10.0 21.0 494(802L) ENDO-012	1.4/10.0 21.0 494(802L) ENDO-014
697-012M	697-014M	697-016M	697-018M	002-012M	002-014M	002-016M	002-018M	494-012M	494-014M
圆球型 Round with Long Neck /BR					紧口圆球型 Round with Collor /BC				

1.6/10.0 21.0 494(802L) ENDO-016	1.8/10.0 21.0 494(802L) ENDO-018	1.0(0.8)/1.0 19.0 010(805) SI-45	1.2(0.9)/1.2 19.0 010(805) SI-46	1.4(1.0)/1.6 19.0 010(805) SI-47	1.6(1.1)/1.8 19.0 010(805) SI-48	1.4(1.1)/4.0 19.0 225(807) SI-50	1.6(1.2)/4.0 19.0 225(807) ---	1.8(1.3)/5.0 19.0 225(807) SI-51
494-016M	494-018M	010-010M	010-012M	010-014M	010-016M	225-014M	225-016M	225-018M
倒锥型 Inverted Cone / SI								

1.0/2.3 19.0 019(806) DI-41	1.4/3.2 19.0 019(806) DI-42	3.2/5.0 19.0 038(811) EX-11	3.5/7.0 20.0 039(811L) EX-12	4.5/0.7 19.0 044(818) WR-44	4.2/0.9 19.0 067(909) WR-11	4.2/1.8 19.0 068(909) WR-13
019-010M	019-014M	038-032M	039-035M	044-045M	067-042M	068-042M
紧口倒锥型 Inverted Cone with collar/DI		双圆锥型 Double cone/EX		轮型 Wheel/WR		

1.0/4.0 19.0 109(835) SF-41	1.2/4.0 19.0 109(835) SF-31	1.4/4.0 19.0 109(835) SF-42	1.6/4.0 19.0 109(835) ---	1.2/6.0 21.0 110(836) SF-20	1.4/6.0 21.0 110(836) SF-21	1.6/6.0 21.0 110(836) SF-15	1.2/8.0 22.0 111(837) SF-11	1.4/8.0 22.0 111(837) SF-12	1.6/8.0 22.0 111(837) SF-13
109-010M	109-012M	109-014M	109-016M	110-012M	110-014M	110-016M	111-012M	111-014M	111-016M
平头圆柱型 Flat End Cylinder/SF									

1.2/8.0 22.0 130(885) CP-11	1.4/8.0 22.0 130(885) CP-12	1.6/8.0 22.0 130(885) CP-13	1.2/10.0 24.0 131(886) ---	1.4/10.0 24.0 131(886) ---	1.6/10.0 24.0 131(886) ---	1.4(0.7)/10.0 24.0 173(848) TF-11
130-012M	130-014M	130-016M	131-012M	131-014M	131-016M	173-014M
尖头圆柱型 Ogival end Cylinder/CP						平头锥型
						规格
						型号

金刚石粒度

- 黑色/极粗粒度
- 绿色/粗粒度
- 蓝色/标准粒度
- 红色/较细粒度
- 黄色/极细粒度

	1.0/4.0 19.0 139(838) SR-41		1.2/4.0 19.0 139(838) SR-31		1.4/4.0 19.0 139(838) SR-42		1.4/6.0 21.0 140(880) SR-21		1.2/8.0 22.0 141(881) SR-11		1.4/8.0 22.0 141(881) SR-12		1.6/8.0 22.0 141(881) SR-13		1.2/0.2 21.0 150(839) EX-18
139-010M 139-012M 139-014M 140-014M 141-012M 141-014M 141-016M 150-012M															
圆头柱型 Round End, Cylinder/SR														肩台车针 End Cut	

	1.2/6.0 21.0 164(852) TC-26		1.4/8.0 22.0 165(858) TC-21		1.0/10.0 22.0 166(859) TC-27		1.2/10.0 22.0 166(859) TC-28		1.4/10.0 22.0 166(859) TC-10		1.6/10.0 22.0 166(859) TC-11		1.8/10.0 22.0 166(859) TC-16
164-012M 165-014M 166-010M 166-012M 166-014M 166-016M 166-018M													
尖头锥型 Pointed cone/TC													

	1.0(0.6)/3.0 19.0 169(845) TF-41		1.2(0.7)/4.0 19.0 170(845) TF-42		1.4(0.8)/4.0 19.0 170(845) TF-43		1.6(0.9)/4.0 19.0 170(845) TF-31		1.8(1.0)/5.0 20.0 170(845) TF-23		2.1(1.5)/5.0 20.0 170(845) TF-22		1.4(0.7)/6.5 21.0 171(846) TF-20		1.6(0.8)/6.5 21.0 171(846) TF-21		1.0(0.6)/8.0 22.0 172(847) TF-07		1.2(0.7)/8.0 22.0 172(847) TF-08
169-010M 170-012M 170-014M 170-016M 170-018M 170-021M 171-014M 171-016M 172-010M 172-012M																			
平头锥型 Flat End Cone/TF																			

	1.4(0.7)/8.0 22.0 172(847) TF-09		1.6(0.8)/8.0 22.0 172(847) TF-10		1.8(0.9)/8.0 22.0 172(847) TF-56		2.3(1.3)/9.0 22.0 172(847) TF-14		1.4(0.7)/10.0 24.0 173(848) TF-11		1.6(0.8)/10.0 22.0 173(848) TF-12		1.8(0.9)/10.0 22.0 173(848) TF-13
172-014M 172-016M 172-018M 172-023M 173-014M 173-016M 173-018M													

	1.0(0.6)/4.0 19.0 196(849) CR-41		1.2(0.7)/4.0 19.0 196(849) CR-31		1.4(0.8)/4.0 19.0 196(849) CR-42		1.6(0.9)/4.0 19.0 196(849) CR-12		2.0(1.0)/5.0 19.0 196(849) CR-11		1.2(0.7)/6.5 21.0 197(855) ---		1.4(0.7)/6.5 21.0 197(855) TR-20		1.6(0.8)/6.5 21.0 197(855) TR-21		1.8(1.0)/6.5 21.0 197(855) TR-62		2.0(1.1)/6.5 21.0 197(855) TR-63
196-010M		196-012M		196-014M		196-016M		196-020M		197-012M		197-014M		197-016M		197-018M		197-020M	
圆头锥型 Round End Cone/CR/TR																			

	1.2(0.7)/8.0 22.0 198(856) ---		1.4(0.8)/8.0 22.0 198(856) TR-25		1.6(0.9)/8.0 22.0 198(856) TR-24		1.8(1.0)/8.0 22.0 198(856) TR-23		1.2(0.6)/10.0 22.0 199(850) ---		1.4(0.7)/10.0 22.0 199(850) TR-11		1.6(0.8)/10.0 22.0 199(850) TR-12		1.8(0.9)/9.0 22.0 199(850) TR-13		2.2(1.3)/9.0 22.0 199(850) TR-14		2.2(1.3)/11.0 24.0 200(850L) TR-15
198-012M 198-014M 198-016M 198-018M 199-012M 199-014M 199-016M 199-018M 199-022M 200-022M																			

	2.4(1.4)/11.0 24.0 200(850L) TR-19		1.0(0.8)/2.8 19.0 233(808) EX-41		1.4(1.1)/4.0 19.0 234(808L) EX-31		1.8(1.3)/5.0 19.0 237(808R) EX-20		2.1(1.6)/5.0 19.0 237(808R) EX-21		3.2(1.8)/5.0 19.0 237(808R) EX-26
200-024M 233-010M 234-014M 237-018M 237-021M 237-032M											
圆头倒截锥型 Round End Inverted cone/EX											

 1.2/6.0 21.0 248(861) FO-41	 1.4/6.0 21.0 248(861) FO-20	 1.4/8.0 22.0 249(862) FO-21	 1.6/8.0 22.0 249(862) FO-22	 1.2/10.0 24.0 250(863) FO-11	 1.4/10.0 24.0 250(863) FO-12	 1.6/10.0 24.0 250(863) FO-13	 1.6/3.2 19.0 257(368,830) FO-30	 1.8/3.6 19.0 257(368,830) FO-32
248-012M	248-014M	249-014M	249-016M	250-012M	250-014M	250-016M	257-016M	257-018M
尖头圆柱型 Ogival end Cylinder/FO							花蕾芽型 Bud/FO	

 2.1/5.0 19.0 257(368,830) ---	 2.3/5.0 19.0 257(368,830) FO-23	 2.5/5.0 19.0 257(368,830) FO-24	 2.8/5.0 19.0 257(368,830) FO-25	 3.0/5.0 19.0 257(368,830) FO-26	 3.2/5.0 19.0 257(368,830) FO-27	 1.8/3.5 19.0 277(379,833) ---	 2.0/3.9 19.0 277(379,833) TR-22	 2.3/4.5 19.0 277(379,833) ---	 3.0/5.0 19.0 277(379,833) TR-28
257-021M	257-023M	257-025M	257-028M	257-030M	257-032M	277-018M	277-020M	277-023M	277-030M
蛋型 Egg Football/TR									

 1.2/6.0 21.0 288(867,877) SO-20	 1.2/8.0 22.0 289(868,878) ---	 1.4/8.0 22.0 289(868,878) SO-21	 1.6/8.0 22.0 289(868,878) ---	 1.2/10.0 24.0 290(869,879) ---	 1.4/10.0 24.0 290(869,879) ---	 1.6/10.0 24.0 290(869,879) ---
288-012M	289-012M	289-014M	289-016M	290-012M	290-014M	290-016M
鱼雷圆柱型 Torpedo/SO						

 1.6(1.0)/8.0 22.0 298(878K) ---	 1.8(1.1)/8.0 22.0 298(878K) ---	 1.6(1.0)/10.0 24.0 299(879K) ---	 1.8(1.1)/10.0 24.0 299(879K) ---	 1.4/8.0 22.0 115(842SP) BC-11	 1.6/8.0 22.0 178(848SP) BC-12	 1.6/5.0 19.0 465(820) ---
298-016M	298-018M	299-016M	299-018M	115-014M	178-016M	465-016M
鱼雷圆锥型 Torpedo,Conical				邻面去釉 Interdental		

 cutting depth 0.3mm 1.6(1.0)/6.0 21.0 552(834) VP-20	 cutting depth 0.4mm 1.8(1.0)/6.0 21.0 552(834) VP-21	 cutting depth 0.5mm 2.1(1.1)/6.0 21.0 552(834) VP-22	 1.8(1.1)/4.5 19.0 544(845KR) RS-31 or EX-38	 1.8(1.1)/7.5 21.0 545(846KR) RS-21 or EX-29	 1.6(0.8)/8.0 22.0 546(847KR) 847KR016	 1.8(0.9)/10.0 24.0 553(848KR) 848KR018	 1.9(1.2)/10.5 23.0 553(848KR) RS-11 or EX-17
552-016M	552-018M	552-021M	544-018M	545-018M	546-016M	553-018M	553-019M
定深车针 Depth Marking				平头圆柱型 Taper, Round Edge			


 1.6.5mm FG S 313	 1.0 16.0 001(801) BR-S45	 1.2 16.0 001(801) BR-S46	 1.4 16.0 001(801) BR-S41	 1.2(0.9)/1.2 16.0 010(805) SI-S46	 1.6(1.1)/1.8 16.0 010(805) SI-S48	 1.4/8.0 18.5 165(858) TC-S21	 1.0(0.6)/3.0 16.5 169(845) TF-S41	 1.6(0.9)/4.0 14.5 170(845) TF-SS31
FG Short	001-S010M	001-S012M	001-S014M	010-S012M	010-S016M	165-S014M	169-S010M	170-SS016M
短柄 Short								

 1.6(0.9)/4.0 16.5 170(845) TF-S31	 1.8(1.0)/5.0 17.0 170(845) TF-S23	 2.1(1.5)/5.0 17.0 170(845) TF-S22	 1.6(0.8)/6.5 18.5 171(846) TF-S21	 1.4(0.7)/6.5 18.5 171(846) TF-S20	 1.4(0.7)/6.5 18.5 197(855) TR-S20	 1.6(0.8)/6.5 18.5 197(855) TR-S21	 1.8(0.9)/9.0 21.0 199(850) TR-S13	 1.4/8.0 19.5 249(862) FO-S21
170-S016M	170-S018M	170-S021M	171-S016M	171-S014M	197-S014M	197-S016M	199-S018M	249-S014M

								
Super Coarse Black	1.4 19.0 001(801) BR-41S	1.6 19.0 001(801) BR-40S	1.8 19.0 001(801) BR-31S	2.1 19.0 001(801) BR-28S	2.3 19.0 001(801) BR-26S	2.5 19.0 001(801) BR-25S	1.2/10.0 21.0 494(802L) ENDO-012S	1.4/10.0 21.0 494(802L) ENDO-014S
黑色环极粗砂	001-014S	001-016S	001-018S	001-021S	001-023S	001-025S	494-012S	494-014S
圆球型 Ball , Round /BR							根管治疗车针 /Root canal therapy	

									
1.6/10.0 21.0 494(802L) ENDO-016S	1.8/10.0 21.0 494(802L) ENDO-018S	1.6(1.1)/1.8 19.0 010(805) SI-48S	3.2/5.0 19.0 038(811) EX-11S	3.5/7.0 20.0 039(811L) EX-12S	1.8(1.3)/5.0 19.0 237(808R) EX-20S	2.1(1.6)/5.0 19.0 237(808R) EX-21S	1.2/4.0 19.0 109(835) SF-31S	1.4/4.0 19.0 109(835) SF-42S	1.4/6.0 21.0 110(836) SF-21S
494-016S	494-018S	010-016S	038-032S	039-035S	237-018S	237-021S	109-012S	109-014S	110-014S
倒锥型 Inverted Cone /SI			双圆锥型 Double cone/EX		圆头倒锥型 Round End Inverted cone/EX		平头圆柱型 Flat End Cylinder/SF		

									
1.6/6.0 21.0 110(836) SF-15S	1.2/8.0 22.0 111(837) SF-11S	1.4/8.0 22.0 111(837) SF-12S	1.6/8.0 22.0 111(837) SF-13S	1.2/8.0 22.0 141(881) SR-11S	1.4/8.0 22.0 141(881) SR-12S	1.6/8.0 22.0 141(881) SR-13S	1.4/8.0 22.0 165(858) TC-21S	1.4/10.0 22.0 166(859) TC-10S	1.6/10.0 22.0 166(859) TC-11S
110-016S	111-012S	111-014S	111-016S	141-012S	141-014S	141-016S	165-014S	166-014S	166-016S
圆头柱型 Round End, Cylinder/SR							尖头锥型 Pointed cone/TC		

						
1.8/10.0 22.0 166(859) TF-16S	1.6(0.8)/6.5 21.0 171(846) TF-21S	1.4(0.7)/8.0 22.0 172(847) TF-09S	1.6(0.8)/8.0 22.0 172(847) TF-10S	1.4(0.7)/10.0 24.0 173(848) TF-11S	1.6(0.8)/10.0 22.0 173(848) TF-12S	1.8(0.9)/10.0 22.0 173(848) TF-13S
166-018S	171-016S	172-014S	172-016S	173-014S	173-016S	173-018S
平头锥型 Flat End Cone/TF						

					
1.6(0.8)/6.5 21.0 197(855) TR-21S	1.4(0.8)/8.0 22.0 198(856) TR-25S	1.6(0.9)/8.0 22.0 198(856) TR-24S	1.4(0.7)/10.0 22.0 199(850) TR-11S	1.6(0.8)/10.0 22.0 199(850) TR-12S	1.8(0.9)/9.0 22.0 199(850) TR-13S
197-016S	198-014S	198-016S	199-014S	199-016S	199-018S
圆头锥型 Round End Cone/CR/TR					

						
1.2/6.0 21.0 248(861) FO-41S	1.4/6.0 21.0 248(861) FO-20S	1.4/8.0 22.0 249(862) FO-21S	1.6/8.0 22.0 249(862) FO-22S	1.2/10.0 24.0 250(863) FO-11S	1.4/10.0 24.0 250(863) FO-12S	1.6/10.0 24.0 250(863) FO-13S
248-012S	248-014S	249-014S	249-016S	250-012S	250-014S	250-016S
尖头圆柱型 Ogival end Cylinder/FO						

								
1.8/3.6 19.0 257(368,830) FO-32S	2.1/5.0 19.0 257(368,830) FO-34S	2.3/5.0 19.0 257(368,830) FO-24S	2.5/5.0 19.0 257(368,830) FO-25S	2.8/5.0 19.0 257(368,830) FO-26S	3.0/5.0 19.0 257(368,830) FO-27S	2.0/3.9 19.0 277(379,833) TR-22S	3.0/5.0 19.0 277(379,833) TR-28S	1.4/8.0 22.0 289(868,878) SO-21S
257-018S	257-021S	257-023S	257-025S	257-028S	257-030S	257-032S	277-020S	277-030S
花蕾芽型 Bud/FO							蛋型 Egg Football/TR	
								鱼雷圆柱型 Torpedo/SO

								
Coarse Green	001(801) BR-46C	001(801) BR-41C	001(801) BR-40C	001(801) BR-31C	001(801) BR-28C	001(801) BR-26C	001(801) BR-25C	001(801) BR-23C
绿色环粗砂	001-012C	001-014C	001-016C	001-018C	001-021C	001-023C	001-025C	001-029C
圆球型 Ball, Round /BR								

								
494(802L) ENDO-012C	494(802L) ENDO-014C	494(802L) ENDO-016C	494(802L) ENDO-018C	010(805) SI-47C	010(805) SI-48C	038(811) EX-11C	039(811L) EX-12C	068(909) WR-13C
494-012C	494-014C	494-016C	494-018C	010-014C	010-016C	038-032C	039-035C	068-042C
根管治疗车针/Root canal therapy			倒锥型 Inverted Cone / SI			双圆锥型 Double cone/EX		轮型 Wheel/WR

								
109(835) SF-31C	109(835) SF-42C	109(835) ---	110(836) SF-20C	110(836) SF-21C	110(836) SF-15C	111(837) SF-11C	111(837) SF-12C	111(837) SF-13C
109-012C	109-014C	109-016C	110-012C	110-014C	110-016C	111-012C	111-014C	111-016C
平头圆柱型 Flat End Cylinder/SF								

					
130(885) CP-11C	130(885) CP-12C	130(885) CP-13C	131(886) ---	131(886) ---	131(886) ---
130-012C	130-014C	130-016C	131-012C	131-014C	131-016C
Ogival end Cylinder/CP					

					
140(880) SR-21C	141(881) SR-11C	141(881) SR-12C	141(881) SR-13C	142(882) ---	142(882) ---
140-014C	141-012C	141-014C	141-016C	142-014C	142-016C
圆头柱型 Round End, Cylinder/SR					

						
165(858) TC-21C	166(859) TC-10C	166(859) TC-11C	166(859) TC-16C	167(859L) ---	167(859L) ---	167(859L) ---
165-014C	166-014C	166-016C	166-018C	167-014C	167-016C	167-018C
尖头锥型 Pointed cone/TC						

								
170(845) TF-43C	171(846) TF-20C	172(847) TF-09C	172(847) TF-10C	172(847) TF-56C	172(847) TF-14C	173(848) TF-11C	173(848) TF-12C	173(848) TF-13C
170-014C	171-014C	172-014C	172-016C	172-018C	172-023C	173-014C	173-016C	173-018C
平头锥型 Flat End Cone/TF								

	1.4(0.7)/6.5 21.0 197(855) TR-20C		1.6(0.8)/6.5 21.0 197(855) TR-21C		1.8(1.0)/6.5 21.0 197(855) TR-62C		1.4(0.8)/8.0 22.0 198(856) TR-25C		1.6(0.9)/8.0 22.0 198(856) TR-24C		1.8(1.0)/8.0 22.0 198(856) TR-23C		1.4(0.7)/10.0 22.0 199(850) TR-11C		1.6(0.8)/10.0 22.0 199(850) TR-12C		1.8(0.9)/9.0 22.0 199(850) TR-13C		2.2(1.3)/9.0 22.0 199(850) TR-14C
197-014C 197-016C 197-018C 198-014C 198-016C 198-018C 199-014C 199-016C 199-018C 199-022C																			
圓頭錐型 Round End Cone/CR/TR																			


	1.4/8.0 22.0 249(862) FO-21C		1.6/8.0 22.0 249(862) FO-22C		1.2/10.0 24.0 250(863) FO-11C		1.4/10.0 24.0 250(863) FO-12C		1.6/10.0 24.0 250(863) FO-13C		1.8/9.5 21.0 190 FO-54C
249-014C 249-016C 250-012C 250-014C 250-016C 190-018C											
尖頭圓柱型 Ogival end Cylinder/FO											

	1.8/3.6 19.0 257(368,830) FO-32C		2.1/5.0 19.0 257(368,830) ---		2.3/5.0 19.0 257(368,830) FO-34C		2.5/5.0 19.0 257(368,830) FO-24C		2.8/5.0 19.0 257(368,830) FO-25C
257-018C 257-021C 257-023C 257-025C 257-028C									
花蕾芽型 Bud/FO									


	1.2/8.0 22.0 289(868,878) ---		1.4/8.0 22.0 289(868,878) SO-21C		1.6/8.0 22.0 289(868,878) ---		1.2/10.0 24.0 290(869,879) ---		1.4/10.0 24.0 290(869,879) ---		1.6/10.0 24.0 290(869,879) ---
289-012C 289-014C 289-016C 290-012C 290-014C 290-016C											
魚雷圓柱型 Torpedo/SO											

					
 Fine Red	1.2 19.0 001(801) BR-46F	1.4 19.0 001(801) BR-41F	1.6 19.0 001(801) BR-40F	1.8 19.0 001(801) BR-31F	2.1 19.0 001(801) BR-28F
紅色環紋細砂					
圓球型 Ball , Round /BR					

								
1.4/6.0 21.0 110(836) SF-21F	1.2/8.0 22.0 111(837) SF-11F	1.4/8.0 22.0 111(837) SF-12F	1.6/8.0 22.0 111(837) SF-13F	1.4/6.0 21.0 140(880) SR-21F	1.2/8.0 22.0 141(881) SR-11F	1.4/8.0 22.0 141(881) SR-12F	1.6/8.0 22.0 141(881) SR-13F	1.2/0.2 21.0 150(839) EX-18F
110-014F 111-012F 111-014F 111-016F 140-014F 141-012F 141-014F 141-016F 150-012F								
平頭圓柱型 Flat End Cylinder/SF			圓頭柱型 Round End, Cylinder/SR				End Cut	

							
1.2/6.0 21.0 164(852) TC-26F	1.2/8.0 22.0 165(858) ---	1.4/8.0 22.0 165(858) TC-21F	1.0/10.0 22.0 166(859) TC-27F	1.2/10.0 22.0 166(859) TC-28F	1.4/10.0 22.0 166(859) TC-10F	1.6/10.0 22.0 166(859) TC-11F	1.8/10.0 22.0 166(859) TC-16F
164-012F 165-012F 165-014F 166-010F 166-012F 166-014F 166-016F 166-018F							
尖頭錐型 Pointed cone/TC							

	1.4(0.8)/4.0 19.0 170(845) TF-43F		1.4(0.7)/6.5 21.0 171(846) TF-20F		1.6(0.8)/6.5 21.0 171(846) TF-21F		1.2(0.7)/8.0 22.0 172(847) TF-08F		1.4(0.7)/8.0 22.0 172(847) TF-09F		1.6(0.8)/8.0 22.0 172(847) TF-10F		1.4(0.7)/10.0 24.0 173(848) TF-11F		1.6(0.8)/10.0 22.0 173(848) TF-12F		1.8(0.9)/10.0 22.0 173(848) TF-13F
170-014F	171-014F	171-016F	172-012F	172-014F	172-016F	173-014F	173-016F	173-018F									
平头锥型 Flat End Cone/TF																	

	1.6(0.9)/4.0 19.0 196(849) CR-12F		2.0(1.0)/5.0 19.0 196(849) CR-11F		1.4(0.7)/6.5 21.0 197(855) TR-20F		1.6(0.8)/6.5 21.0 197(855) TR-21F		1.8(1.0)/6.5 21.0 197(855) TR-62F		1.4(0.8)/8.0 22.0 198(856) TR-25F		1.6(0.9)/8.0 22.0 198(856) TR-24F		1.8(1.0)/8.0 22.0 198(856) TR-23F		1.2(0.6)/10.0 22.0 199(850) ---		1.4(0.7)/10.0 22.0 199(850) TR-11F
196-016F	196-020F	197-014F	197-016F	197-018F	198-014F	198-016F	198-018F	199-012F	199-014F										
圆头锥型 Round End Cone/CR/TR																			

	1.6(0.8)/10.0 22.0 199(850) TR-12F		1.8(0.9)/9.0 22.0 199(850) TR-13F		2.2(1.3)/9.0 22.0 199(850) TR-14F		1.4(1.1)/4.0 19.0 234(808L) EX-31F		1.8(1.3)/5.0 19.0 237(808R) EX-20F		2.1(1.6)/5.0 19.0 237(808R) EX-21F
199-016F	199-018F	199-022F	234-014F	237-018F	237-021F						
圆头倒截锥型 Round End Inverted cone/EX											

	1.2/6.0 21.0 248(861) FO-41F		1.4/6.0 21.0 248(861) FO-20F		1.4/8.0 22.0 249(862) FO-21F		1.6/8.0 22.0 249(862) FO-22F		1.2/10.0 24.0 250(863) FO-11F		1.4/10.0 24.0 250(863) FO-12F		1.6/10.0 24.0 250(863) FO-13F
248-012F	248-014F	249-014F	249-016F	250-012F	250-014F	250-016F							
尖头圆柱型 Ogival end Cylinder/FO													

	1.6/3.2 19.0 257(368,830) FO-30F		1.8/3.6 19.0 257(368,830) FO-32F		2.1/5.0 19.0 257(368,830) ---		2.3/5.0 19.0 257(368,830) FO-34F		2.5/5.0 19.0 257(368,830) FO-24F		2.8/5.0 19.0 257(368,830) FO-25F		3.0/5.0 19.0 257(368,830) FO-26F		3.2/5.0 19.0 257(368,830) FO-27F
257-016F	257-018F	257-021F	257-023F	257-025F	257-028F	257-030F	257-032F								
花蕾芽型 Bud/FO															

	1.2/6.0 21.0 288(867,877) SO-20F		1.2/8.0 22.0 289(868,878) ---		1.4/8.0 22.0 289(868,878) SO-21F		1.6/8.0 22.0 289(868,878) ---
288-012F	289-012F	289-014F	289-016F				
魚雷圓柱型 Torpedo/SO							

	0.6 16.8 001 CD-50F		0.7/2.2 16.8 138 CD-51F		0.7/3.8 16.8 139 CD-52F		0.7/2.2 16.8 161 CD-53F		0.6/2.2 16.8 246 CD-54F		0.9/1.2 16.8 254 CD-55F		0.7/2.2 16.8 246 CD-56F		0.7/3.6 16.8 247 CD-57F		0.8/3.2 16.8 108 CD-58F		1.3/6.9 16.8 171 CD-59F
001-S006F	138-S007F	139-S007F	161-S007F	246-S006F	254-S009F	246-S007F	247-S007F	108-S008F	171-S013F										
儿童车针 FG Short Burs for Kids / CD																			



Extra Fine
Yellow



1.2
19.0
001(801)
BR-46EF



1.4
19.0
001(801)
BR-41EF



1.6
19.0
001(801)
BR-40EF



1.8
19.0
001(801)
BR-31EF



2.1
19.0
001(801)
BR-28EF



2.3
19.0
001(801)
BR-26EF

黄色环 细砂

001-012EF

001-014EF

001-016EF

001-018EF

001-021EF

001-023EF

圆球型 Ball, Round /BR



1.2/6.0
21.0
164(852)
TC-26EF



1.2/8.0
22.0
165(858)



1.4/8.0
22.0
165(858)
TC-21EF



1.0/10.0
22.0
166(859)
TC-27EF



1.2/10.0
22.0
166(859)
TC-28EF



1.4/10.0
22.0
166(859)
TC-10EF



1.6/10.0
22.0
166(859)
TC-11EF



1.8/10.0
22.0
166(859)
TC-16EF

164-012EF

165-012EF

165-014EF

166-010EF

166-012EF

166-014EF

166-016EF

166-018EF

尖头锥型 Pointed cone/TC



1.6(0.8)/6.5
21.0
171(846)
TF-21EF



1.6(0.8)/8.0
22.0
172(847)
TF-10EF



1.4(0.7)/10.0
24.0
173(848)
TF-11EF



1.6(0.8)/10.0
22.0
173(848)
TF-12EF



1.8(0.9)/10.0
22.0
173(848)
TF-13EF

171-016EF

172-016EF

173-014EF

173-016EF

173-018EF

平头锥型 Flat End Cone/TF



1.6(0.9)/4.0
19.0
196(849)
CR-12EF



2.0(1.0)/5.0
19.0
196(849)
CR-11EF



1.2(0.7)/6.5
21.0
197(855)



1.4(0.7)/6.5
21.0
197(855)
TR-20EF



1.6(0.8)/6.5
21.0
197(855)
TR-21EF



1.8(1.0)/6.5
21.0
197(855)
TR-62EF




1.2(0.7)/8.0
22.0
198(856)



1.4(0.8)/8.0
22.0
198(856)
TR-25EF



1.6(0.9)/8.0
22.0
198(856)
TR-24EF



1.8(1.0)/8.0
22.0
198(856)
TR-23EF

196-016EF

196-020EF

197-012EF

197-014EF

197-016EF

197-018EF

198-012EF

198-014EF

198-016EF

198-018EF

圆头锥型 Round End Cone/CR/TR



1.2(0.6)/10.0
22.0
199(850)



1.4(0.7)/10.0
22.0
199(850)
TR-11EF



1.6(0.8)/10.0
22.0
199(850)
TR-12EF



1.8(0.9)/9.0
22.0
199(850)
TR-13EF



1.4(1.1)/4.0
19.0
234(808L)
EX-31EF



1.8(1.3)/5.0
19.0
237(808R)
EX-20EF



2.1(1.6)/5.0
19.0
237(808R)
EX-21EF

199-012EF

199-014EF

199-016EF

199-018EF

234-014EF

237-018EF

237-021EF

圆头倒截锥型 Round End Inverted cone/EX



1.0/4.3
23.5
540(889)
FO-40EF



1.2/6.0
21.0
248(861)
FO-41EF



1.4/6.0
21.0
248(861)
FO-20EF



1.4/8.0
22.0
249(862)
FO-21EF



1.6/8.0
22.0
249(862)
FO-22EF



1.2/10.0
24.0
250(863)
FO-11EF



1.4/10.0
24.0
250(863)
FO-12EF



1.6/10.0
24.0
250(863)
FO-13EF

540-010EF

248-012EF

248-014EF

249-014EF

249-016EF

250-012EF

250-014EF

250-016EF

尖头圆柱型 Ogival end Cylinder/FO



1.6/3.2
19.0
257(368,830)
FO-30EF



1.8/3.6
19.0
257(368,830)
FO-32EF



2.1/5.0
19.0
257(368,830)



2.3/5.0
19.0
257(368,830)
FO-34EF



2.5/5.0
19.0
257(368,830)
FO-24EF



2.8/5.0
19.0
257(368,830)
FO-25EF



3.0/5.0
19.0
257(368,830)
FO-26EF



3.2/5.0
19.0
257(368,830)
FO-27EF



1.2/6.0
21.0
288(867,877)
SO-20EF



1.4/8.0
22.0
289(868,878)
SO-21EF

257-016EF

257-018EF

257-021EF

257-023EF

257-025EF

257-028EF

257-030EF


257-032EF

288-012EF

289-014EF

花蕾芽型 Bud/FO

鱼雷圆柱型 Torpedo/SO

 25mm FG XL 316									
	1.0 25.0 001(801) BR-X45	1.2 25.0 001(801) BR-X46	1.4 25.0 001(801) BR-X41	1.6 25.0 001(801) BR-X40	1.8 25.0 001(801) BR-X31	1.2(0.9)/1.2 25.0 010(805) SI-X46	1.2/8.0 28.0 111(837) SF-X11	1.4/8.0 28.0 111(837) SF-X12	
FG Surgical 加长柄	001-X010M	001-X012M	001-X014M	001-X016M	001-X018M	010-X012M	111-X012M	111-X014M	
1.6/8.0 28.0 111(837) SF-X13	1.6/10.0 29.0 166(859) TC-X11	1.4(0.7)/10.0 29.0 173(848) TF-X11	1.6(0.8)/10.0 29.0 173(848) TF-X12	1.8(0.9)/10.0 29.0 173(848) TF-X13	1.6(0.8)/10.0 29.0 199(850) TR-X12	2.2(1.3)/11.0 29.0 200(850L) TR-X15	1.4/8.0 28.0 249(862) FO-X21	1.6/8.0 28.0 249(862) FO-X22	
111-X016M	166-X016M	173-X014M	173-X016M	173-X018M	199-X016M	200-X022M	249-X014M	249-X016M	
Super Coarse 1.6 25.0 001(801) BR-X40S	1.8 25.0 001(801) BR-X31S	Coarse 1.4 25.0 001(801) BR-X41C	1.6 25.0 001(801) BR-X40C	1.8 25.0 001(801) BR-X31C					
001-X016S	001-X018S	001-X014C	001-X016C	001-X018C					
 22.0mm RA, CA 204									
	1.2 22.0 001(801) BR-R46	1.4 22.0 001(801) BR-R41	1.6 22.0 001(801) BR-R40	1.8 22.0 001(801) BR-R31	2.3 22.0 001(801) BR-R26	1.0(0.8)/1.0 22.0 010(805) SI-R45	1.2(0.9)/1.2 22.0 010(805) SI-R46	1.4(1.0)/1.6 22.0 010(805) SI-R47	
RA, CA 弯机柄	001-R012M	001-R014M	001-R016M	001-R018M	001-R023M	010-R010M	010-R012M	010-R014M	
1.6(1.1)/1.8 22.0 010(805) SI-R48	1.0/4.0 22.0 109(835) SF-R41	1.2/4.0 22.0 109(835) SF-R31	1.4/4.0 22.0 109(835) SF-R42	1.2/8.0 24.0 111(837) SF-R11	1.4/8.0 24.0 111(837) SF-R12	1.6/8.0 24.0 111(837) SF-R13	1.2/8.0 24.0 141(881) SR-R11	1.4/8.0 24.0 141(881) SR-R12	1.4/8.0 24.0 165(858) TC-R21
010-R016M	109-R010M	109-R012M	109-R014M	111-R012M	111-R014M	111-R016M	141-R012M	141-R014M	165-R014M
1.6/10.0 27.0 166(859) TC-R11	1.0(0.6)/3.0 22.0 169(845) TF-R41	1.2(0.7)/4.0 22.0 170(845) TF-R42	1.8(1.0)/8.0 25.0 198(856) TR-R23	1.6(0.8)/10.0 27.0 199(850) TR-R12	1.8/3.6 22.0 257(368,830) FO-R32	1.6 22.0 001(801) BR-R40F	1.8 22.0 001(801) BR-R31F	2.3 22.0 001(801) BR-R26F	1.4/8.0 24.0 165(858) TC-R21F
166-R016M	169-R010M	170-R012M	198-R018M	199-R016M	257-R018M	001-R016F	001-R018F	001-R023F	165-R014F
1.8/3.6 22.0 257(368,830) FO-R32F	Extra Fine 1.6 22.0 001(801) BR-R40EF	1.8 22.0 001(801) BR-R31EF	2.3 22.0 001(801) BR-R26EF	去除肉芽组织 车针套装 Degranulation Kit					
257-R018F	001-R016EF	001-R018EF	001-R023EF	001-R010	001-R025	001-R030	001-R035		

用于制备烧结氧化锆的金刚石车针

Diamond burs for preparing sintered Zirconia

氧化锆是一种具有良好发展前景的假牙材料。氧化锆的特殊性能对制备提出了极高的要求。本系列金刚石钻头专门针对氧化锆陶瓷的极高硬度而设计。它们具有比标准金刚石钻头更长的使用寿命和更大的切削能力。它们是根管治疗中切割氧化锆冠的理想工具。

Zirconia is a material with a promising future in prosthetic dentistry. The special properties of zirconia place extreme demands on the preparation. This series of our diamond burs are tailored to the extreme hardness of zirconia ceramic. They have a long service life and much greater cutting capacity than standard diamond burs. They are ideal for sectioning zirconia crowns during endodontic treatment.

									
2.1/10.0 24.0	3.2/5.0 19.0	1.2/8.0 22.0	1.4/8.0 22.0	2.3(1.3)/9.0 22.0	1.6(0.8)/10.0 22.0	1.4/10.0 24.0	1.6(0.8)/8.0 22.0	1.8(0.9)/10.0 24.0	1.6(0.9)/8.0 22.0
494 (802L)	038 (811)	141 (881)	141 (881)	172 (847)	173 (848)	290 (869,879)	546 (847KR)	553 (848KR)	198 (856)
Z494-021M	Z038-032M	Z141-012M	Z141-014M	Z172-023M	Z173-016M	Z290-014M	Z546-016M	Z553-018M	Z198-016M
—	ZEX-11	ZSR-11	ZSR-12	ZTF-14	ZTF-12	—	—	—	ZTR-24
									
1.4 22.0	1.2/8.0 22.0	1.4/8.0 22.0	1.4/10.0 24.0	1.6(0.8)/8.0 22.0	1.8(0.9)/10.0 24.0	1.6(0.9)/8.0 22.0	2.5/5.0 19.0	1.6/3.2 19.0	1.0/11.5 24.0
697 (801L)	141 (881)	141 (881)	290 (869,879)	546 (847KR)	553 (848KR)	198 (856)	257 (368,830)	277 (379,833)	167 (859L)
Z697-014F	Z141-012F	Z141-014F	Z290-014F	Z546-016F	Z553-018F	Z198-016F	Z257-025F	Z277-016F	Z167-010F
ZBR-L41F	ZSR-11F	ZSR-12F	—	—	—	ZTR-24F	ZFO-24F	—	—
									
1.4 22.0	1.2/8.0 22.0	1.4/8.0 22.0	1.4/10.0 24.0	1.6(0.8)/8.0 22.0	1.8(0.9)/10.0 24.0	1.6(0.9)/8.0 22.0	2.5/5.0 19.0	1.6/3.2 19.0	1.0/11.5 24.0
697 (801L)	141 (881)	141 (881)	290 (869,879)	546 (847KR)	553 (848KR)	198 (856)	257 (368,830)	277 (379,833)	167 (859L)
Z697-014EF	Z141-012EF	Z141-014EF	Z290-014EF	Z546-016EF	Z553-018EF	Z198-016EF	Z257-025EF	Z277-016EF	Z167-010EF
ZBR-L41EF	ZSR-11EF	ZSR-12EF	—	—	—	ZTR-24EF	ZFO-24EF	—	—

快速切削全陶瓷，烧结氧化锆冠的黑色金刚石车针（超粗）

Fast reducing the All-Ceramics, Zirconia crown with the Black burs (Super coarse)

									
1.2/8.0 22.0	1.4/8.0 22.0	1.6(0.8)/10.0 22.0	1.6/10.0 22.0	1.8/10.0 22.0	1.8(0.9)/9.0 22.0	2.2(1.3)/9.0 22.0	1.6(0.8)/8.0 22.0	1.8(0.9)/10.0 24.0	2.3/5.0 19.0
141 (881)	141 (881)	173 (848)	166 (859)	166 (859)	199 (850)	199 (850)	546 (847KR)	553 (848KR)	257 (368,830)
Z141-012S	Z141-014S	Z173-016S	Z166-016S	Z166-018S	Z199-018S	Z199-022S	Z546-016S	Z553-018S	Z257-023S
ZSR-11S	ZSR-12S	ZTF-12S	ZTC-11S	ZTC-16S	ZTR-13S	ZTR-14S	—	—	ZFO-34S

牙体制备专用车针

1. 高品质专业备牙车针，有前后牙之分。
2. 制备肩台时不会出现羽状肩台，避免了因备牙而引起的修复体无法就位和颈缘崩裂现象。
3. 保证肩台宽度一致。

Especially for All-ceramic restoration. Pin guide not Diamond coated avoids injury of the gingival edge;
508T-014XC wearing-off the contact points mesial-distal;
001-014M determination of prep-bounds
198C-016XC 021XC 508T-016XC 021XC for circular wearing-off without damaging the prep-stage;
508T-016F;021F for precise grinding of the pre-stage with the help of a precision guiding-pin.

									
1.4/8.0 21.5	1.4 19.0	1.6/8.0 20.5	1.6/9.0 21.5	1.6/9.0 21.5	1.4/8.0 21.5	1.6 16.0	2.1/8.0 19.0	2.1/9.0 20.0	2.1/9.0 20.0
508T-014XC	001-014M	198C-016XC	508T-016XC	508T-016F	508T-014XC	001-016M	198C-021XC	508T-021XC	508T-021F



品名	FG 十支黄色塑料盒
Name	FG 10 pcs yellow plastic box
尺寸/Size	44mm*30.5mm
重量/Weight	6.6g



品名	FG 十支红色塑料盒
Name	FG 10 pcs red plastic box
尺寸/Size	44mm*30.5mm
重量/Weight	6.6g



品名	FG 十支蓝色塑料盒
Name	FG 10 pcs blue plastic box
尺寸/Size	44mm*30.5mm
重量/Weight	6.6g



品名	FG 十支绿色塑料盒
Name	FG 10 pcs green plastic box
尺寸/Size	44mm*30.5mm
重量/Weight	6.6g



品名	FG 十支黑色塑料盒
Name	FG 10 pcs black plastic box
尺寸/Size	44mm*30.5mm
重量/Weight	6.6g



品名	FG 五支黄色塑料盒
Name	FG 5 pcs yellow plastic box
尺寸/Size	35mm*29mm
重量/Weight	6.4g



品名	FG 五支红色塑料盒
Name	FG 5 pcs red plastic box
尺寸/Size	35mm*29mm
重量/Weight	6.4g



品名	FG 五支蓝色塑料盒
Name	FG 5 pcs blue plastic box
尺寸/Size	35mm*29mm
重量/Weight	6.4g



品名	FG 五支绿色塑料盒
Name	FG 5 pcs green plastic box
尺寸/Size	35mm*29mm
重量/Weight	6.4g



品名	FG 五支黑色塑料盒
Name	FG 5 pcs black plastic box
尺寸/Size	35mm*29mm
重量/Weight	6.4g



品名	FG 十支推拉盒
Name	FG 10 pcs plastic box
尺寸/Size	35mm*33mm
重量/Weight	6.8g



品名	HP 30支塑料盒
Name	HP 30 pcs plastic box
尺寸/Size	144mm*71mm
重量/Weight	77.2g



品名	HP 十支塑料盒
Name	HP 10 pcs plastic box
尺寸/Size	108mm*69mm
重量/Weight	30.5g



品名	HP 十支吸塑
Name	HP 10pcs blister card
尺寸/Size	61mm*165mm
重量/Weight	22.4g



品名	FG 五支吸塑
Name	FG 5pcs blister card
尺寸/Size	95.5mm*56mm
重量/Weight	4.8g



品名	FG 五支吸塑
Name	FG 5pcs blister card
尺寸/Size	81mm*47mm
重量/Weight	4.2g



品名	HP 五支吸塑
Name	HP 5pcs blister card
尺寸/Size	98.5mm*71mm
重量/Weight	10.5g



品名	FG 十支吸塑
Name	FG 10pcs blister card
尺寸/Size	80mm*69mm
重量/Weight	6.7g



品名	FG 十支吸塑
Name	FG 10pcs blister card
尺寸/Size	105.5mm*80mm
重量/Weight	6.9g



品名	FG 五连包
Name	FG Five packs
尺寸/Size	149.5mm*50mm
重量/Weight	3.2g

产品套装 Product Kit



备冠/前牙套装
Crown Preparation Anterior Teeth Kit



备冠/后牙套装
Crown Preparation Molar Teeth Kit



备洞套装
Cavity Preparation Kit



(高) 嵌体预备套装
Inlay Preparation Kit



烤瓷牙抛光套装
Ceramics/Composite Polishing Kit



儿童微创修复套装
Kiddies Kit



烤瓷制备套装
Ceramic Restoration Kit



备肩台套装
Porcelain Shouldered Abutment Polishing Kit



烤瓷制备/修复套装
Porcelain Preparation / Repair Kit



制备/修复套装
Preparation Repair Kit



护龈综合套装
Gingival Bevel Burs Kit



前牙全冠牙体预备套装
Crown Preparation Anterior Teeth Kit



磨牙全冠牙体预备套装
Crown Preparation Molar Teeth Kit



复合树脂金刚砂车针套装
Composite Resin Kit



微创窝洞制备套装 (粗砂)
Minimally Invasive Cavity Preparation Coarse Kit



微创窝洞制备套装 (细砂)
Minimally Invasive Cavity Preparation Fine Kit



树脂修形套装
Composite Repair Kit



陶瓷复合抛光套装
Porcelain Ceramics Composite Polishing Kit



瓷冠制备套装
Porcelain Crown Preparation Kit



常规粗砂套装
Regular Coarse Sand Kit



磨牙修复套装
Molars Restoration Kit



备肩台套装
Shoulder preparation Kit




备冠/前牙套装
Crown Preparation Anterior Teeth Kit



邻面去釉套装
Interproximal enamel kit



氧化锆金刚砂车针套装
Zirconia Diamond Kit



氧化锆金刚砂车针套装
Diamond Zirconia Finishing Kit



氧化锆金刚砂车针套装
Zirconia Kit



破除氧化锆金刚砂车针套装
Diamond Zirconia Cutting Kit



去除肉芽组织车针套装
Degranulation Kit

Diamond Ortho-strips

金刚石锯弓




Creating additional space by interproximal enamel reduction (IPR) :

When opening the interdental space and reducing, finishing and polishing enamel in orthodontics, care must be taken to avoid grooves and scratches. it's important to proceed step by step, following a coarse-to-fine sequence.



Item No. 1

Thickness (mm)



Separating

OS40

0.40 mm



Reducing

OS30

0.30 mm



Contouring

OS25


0.25 mm



Finishing

OS20

0.20 mm



Pre-polishing

OS15

0.15 mm



Polishing

OS10


0.10 mm



No-Diamond Coated
with serrated edge

OS05

0.05 mm



Diamond finishing strips

金刚石抛光条



Perforated Diamond finishing strips with and without serrated edge ;

Double side Coated / Single side Coated ;

with or without honeycomb design .

	45 μ m	Length x Width x Thickness : 148 x 3.75 x 0.16mm	No. : SWS37		Double side Coated
	45 μ m	Length x Width x Thickness : 148 x 3.75 x 0.13mm	No. : WS37		Single side Coated
	30 μ m	Length x Width x Thickness : 148 x 3.75 x 0.12mm	No. : SWS37F		Double side Coated
	30 μ m	Length x Width x Thickness : 148 x 3.75 x 0.10mm	No. : WS37F		Single side Coated
	15 μ m	Length x Width x Thickness : 148 x 3.75 x 0.10mm	No. : SWS37EF		Double side Coated
	15 μ m	Length x Width x Thickness : 148 x 3.75 x 0.08mm	No. : WS37EF		Single side Coated
	45 μ m	Length x Width x Thickness : 148 x 3.75 x 0.16mm	No. : SDS37		Double side Coated
	45 μ m	Length x Width x Thickness : 148 x 3.75 x 0.13mm	No. : DS37		Single side Coated
	30 μ m	Length x Width x Thickness : 148 x 3.75 x 0.12mm	No. : SDS37F		Double side Coated
	30 μ m	Length x Width x Thickness : 148 x 3.75 x 0.10mm	No. : DS37F		Single side Coated
	15 μ m	Length x Width x Thickness : 148 x 3.75 x 0.10mm	No. : SDS37EF		Double side Coated
	15 μ m	Length x Width x Thickness : 148 x 3.75 x 0.08mm	No. : DS37EF		Single side Coated



10 pcs of each per pack
or 10 pcs mixed 2 • + 4 • + 4 •



Tungsten Carbide Bur

钨钢车针



Advantages

1. Strictly comply with ISO Standard
2. Sharper and More Efficient
3. Durable and Valuable
4. FG, FG Short, FG Long, RA, RA Long etc



Ball, Plain Cut

ISO 500 ... 001	001	005	006	008	010	012	014	016	018	021	023	025
FG 314	FG 1/4	FG 1/2	FG 1	FG 2	FG 3	FG 4	FG 5	FG 6	FG 7	FG 8	FG 9	
FG S 313			FG S1	FG S2	FG S3	FG S4	FG S5					
FG X 316			FG X1	FG X2	FG X3	FG X4	FG X5	FG X6	FG X7	FG X8		
RA 204		RA 1/2	RA 1	RA 2	RA 3	RA 4	RA 5	RA 6	RA 7	RA 8		
RA X 205				RA X2	RA X3	RA X4	RA X5	RA X6	RA X7	RA X8		



Inverted Cone

ISO 500 ... 010	001	008	010	012	014	016	018	021	023
FG 314	FG 34	FG 35	FG 36	FG 37	FG 38	FG 39			
RA 204	RA 34	RA 35	RA 36	RA 37	RA 38	RA 39	RA 40	RA 41	



Cylinder, Plain Cut

ISO 500 ... 107	006	008	009	010	012	014	016
Working Length/mm	4.0	4.0	4.0	4.0	4.0	4.0	4.0
FG 314	FG 55	FG 56	FG 57	FG 58	FG 59	FG 60	
FG X 316			FG X57	FG X58			



Flat End Taper, Plain Cut

ISO 500 ... 168	006	008	009	010	012	016	009	010	012	016
Working Length/mm	4.0	4.0	4.0	4.0	4.0	4.0	6.0	6.0	6.0	6.0
FG 314	FG 168	FG 169	FG 170	FG 171	FG 172	FG 169L	FG 170L	FG 171L	FG 172L	
FG X 316			FG X170	FG X171						



Pear Plain Cut




ISO 500 ... 233 001	009	012	006	008	010	012	010	012
Working Length/mm	2.5	2.9	1.7	1.8	2.0	2.2	4.1	4.1
FG 314	FG 245	FG 246	FG 329	FG 330	FG 331	FG 332	FG 331L	FG 332L
FG S 313	FG S245							
RA 204	RA 245							




Cylinder, Cross Cut



ISO 500 ... 107 007	008	009	010	012	014	016	018	021	023
Working Length/mm	4.0	4.0	4.0	4.0	4.0	4.0	4.0	4.0	4.0
FG 314	FG 555	FG 556	FG 557	FG 558	FG 559	FG 560			
FG S 313		FG S556	FG S557						
FG X 316			FG X557	FG X558	FG X559				
RA 204		RA 556	RA 557	RA 558	RA 559	RA 560	RA 561	RA 562	RA 563
RA X 205			RA X557						



Cylinder, Cross Cut



ISO 500 ... 110 007	009	010	012	014
Working Length/mm	6.0	6.0	6.0	6.0
FG 314	FG 556L	FG 557L	FG 558L	FG 559L



Flat End Taper,
Cross Cut



ISO 500 ... 168 007	008	009	010	012	016	021	023	009	010	012	016
Working Length/mm	4.0	4.0	4.0	4.0	4.0	4.0	4.0	6.0	6.0	6.0	6.0
FG 314	FG 698	FG 699	FG 700	FG 701	FG 702			FG 699L	FG 700L	FG 701L	FG 702L
FG X 316			FG X700	FG X701	FG X702						
RA 204		RA 669	RA 700	RA 701	RA 702	RA 703	RA 704				
RA X 205				RA X701							



							
ISO 500 ... 137 007	008	009	010	012	014	016	
Working Length/mm	4.0	4.0	4.0	4.0	4.0	4.0	
FG 314	FG 1155	FG 1156	FG 1157	FG 1158	FG 1159	FG 1160	
FG S 313	FG S1157						

							
ISO 500 ... 194 006	009	010	012	016	ISO 500 ... 194 007	012	016
Working Length/mm	4.0	4.0	4.0	4.0	Working Length	4.0	4.0
FG 314	FG 1169	FG 1170	FG 1171	FG 1172	FG 314	FG 1701	FG 1702




						
ISO 500 ... 137 007	009	010	012	014		
Working Length/mm	4.0	4.0	4.0	4.0		
FG 314	FG 1556	FG 1557	FG 1558	FG 1559		
FG X 316	FG X1557	FG X1558				

							
ISO 500 ...	009	012	ISO 500 ... 468 ...	014	ISO 500 ... 150 001	010	012
Working Length/mm	3.0	5.5	Working Length/mm	4.0	Working Length/mm	0.2	0.2
FG 314	ET3	ET5	FG 314	FG 468	FG 314	FG 957	FG 958



								
ISO 500 ... 237.008	010	012	ISO 500 ... 137.008	010	012	ISO 500 ... 107 008	010	012
Working Length/mm	2.0	2.2	Working Length/mm	4.0	4.0	Working Length/mm	4.0	4.0
FG 314	FG 1931	FG 1932		FG 1957	FG 1958		FG 2057	FG 2058





Trimming & Finishing Carbide Burs




White Stone 白车石 Finish,Polish Filling Composite 去树脂,填充物及抛光。



Green Stone 绿车石 Grind Ceramics,Precious Metal 研磨烤瓷,贵金属等。






Diamond Stones 陶瓷金刚车石 Finish Zirconia, Lithium Disilicate, Ceramics . 研磨氧化锆、玻璃陶瓷、烤瓷。



Diamond Rubber Polishers 橡胶金刚石抛光磨头

Universally used to polish Zirconia, Lithium Disilicate, Ceramics, Composite etc.
广泛的用于抛光氧化锆、玻璃陶瓷、烤瓷、树脂等。



Dental Restorative Finisher & Polisher

牙科修复体修整与抛光磨头

For contouring, finishing and polishing cured composite restorations, and for polishing enamel surfaces after dental cleaning.

Aluminum Oxide Restorative Finisher

For contouring, finishing and smoothing cured composite restorations.

氧化铝复合树脂修整磨头

用于已固化复合树脂修复体的轮廓修整、表面精修和平滑。

Al₂O₃ • Polyurethane • Medium 200#
10 pieces / box



REF-100WB

10 × 1 mm



REF-060CB

6 × 7 mm



REF-050PB

5 × 11 mm

Diamond Restorative & Enamel Polisher - Medium

For polishing cured composite restorations and enamel surfaces after dental cleaning.

金刚石复合树脂及牙釉质中粗抛光磨头

用于已固化复合树脂修复体抛光及洁牙后的牙釉质表面抛光。

Diamond • Resin bond • Medium
5 pieces / box



REP-100DB

10 × 2 mm



REP-055TB

5.5 × 9.5 mm

Diamond Restorative & Enamel Polisher - Fine

For fine polishing of cured composite restorations and enamel surfaces after dental cleaning.

金刚石复合树脂及牙釉质精细抛光磨头

用于已固化复合树脂修复体精细抛光及洁牙后的牙釉质表面抛光。

Diamond • Resin bond • Fine
5 pieces / box



REP-100DW

10 × 2 mm



REP-055TW

5.5 × 9.5 mm

Cleaning Burs after debonding (脱胶后)洁牙磨头

Remove Calculus or Stains , No damage to the enamel.

 3 pieces mixed



2.3
9.0
Coarse



2.3
7.0



2.3
7.5



2.3
9.0
Medium



2.3
7.0



2.3
7.5

PC00C

PC02C

PC03C

PC00

PC02

PC03

 4 Medium



3.0
7.5



5.0
13.0



10.0
0.7



6.0
9.0

BR101

BR103

BR104

BR106

Silicon
矽粒子

 25



6x10 mm



10x0.7 mm



6x9 mm



SM103

SM104

SM106

53RS-2

Twin-Color Polishers for Composite

双色·树脂专用抛光磨头



Composite, Acrylic teeth's Polishing,
Polishing after remove Calculus or Stains , No damage to the enamel.
Polishing after remove the debonding



Item No.  12

Tr103C

Tr104C

Tr105C

Tr106C

Thickness (mm)

5 x 10.5mm

10 x 1 mm

10 x 2.0 (0.7) mm

6 x 8 mm



Item No.  12

Tr103M

Tr104M

Tr105M

Tr106M

Thickness (mm)

5 x 10.5mm

10 x 1 mm

10 x 2.0 (0.7) mm

6 x 8 mm



Item No.  12

Tr103F

Tr104F

Tr105F

Tr106F

Thickness (mm)

5 x 10.5mm

10 x 1 mm

10 x 2.0 (0.7) mm

6 x 8 mm

Silicon Polisher

For Ceramics

硅胶磨头
烤瓷抛光



White , Coarse

Cf101C

Cf103C

Cr101C

Cr103C

Cr104C

Cr105C

Cr106C

Cr108C

Silicon Polisher

For Ceramics

硅胶磨头
烤瓷抛光



Pink , Medium

Cf101M

Cf103M

Cr101M

Cr103M

Cr104M

Cr105M

Cr106M

Cr108M

Silicon Polisher

For Ceramics, Composite

硅胶磨头
烤瓷、树脂抛光



Blue , Fine

Cf101F

Cf103F

Cr101F

Cr103F

Cr104F

Cr105F

Cr106F

Cr108F

Silicon Polisher

For Metal, Alloy

硅胶磨头
金属、合金抛光



Black , Coarse

Af101C

Af103C

Ar101C

Ar103C

Ar104C

Ar105C

Ar106C

Ar108C

Silicon Polisher

For Metal, Alloy

硅胶磨头
金属、合金抛光



Brown , Medium

Af101M

Af103M

Ar101M

Ar103M

Ar104M

Ar105M

Ar106M

Ar108M

Silicon Polisher

For Metal, Alloy, Composite

硅胶磨头
金属、合金、树脂抛光



Fine , Light Grey

Af101F

Af103F

Ar101F

Ar103F

Ar104F

Ar105F

Ar106F

Ar108F

Silicon Polisher

For Composite

硅胶磨头
树脂抛光



Fine , Yellow

Rf101F

Rf103F

Rr101F

Rr103F

Rr104F

Rr105F

Rr106F

Rr108F

Polishing Brushes

Nature Material, No hurt to teeth
Final polishing composite,
nature teeth or Ceramics
天然毛质、不伤牙
适用于树脂、真牙、烤瓷等



Felt
毛毡

R1Ns



Leather
羊皮

R1Ry



Goat Hair
羊毛

R1G



Thread
白棉线

R1T

Polishing Brushes

Abrasive Fiber
磨料丝 (含碳化硅
或氧化铝)

For Composite,
Ceramics

With Al₂O₃ 氧化铝



R1Al



R4Al



R5Al



R7Al

With Silicon Carbide 碳化硅



R1SiC



R4SiC



R5SiC



R7SiC



White Bristle

R4Dw



White Nylon

R4Cw



Clear Nylon

R4Cc



Blue Nylon

R4Cb



Pink Nylon

R4Cr



Purple Nylon

R4Cp



Yellow Nylon

R4Cy



White Goat Hair

R4G



White Bristle

R7Dw



Clear Nylon

R7Cc



White Goat Hair

R7G



White Bristle

R5Dw



White Nylon

R5Cw



Clear Nylon

R5Cc



Blue Nylon

R5Cb



Pink Nylon

R5Cr



Purple Nylon

R5Cp



Yellow Nylon

R5Cy



White Goat Hair

R5G



4 Angles, Soft

R01W



6 Angles, Hard

R02W

Mandrels

夹 针



31RS



31RSL



35RS



52RS



53RS



26RS



51R



No. 201
for FG burs * 15



No. 211
for 10*FG + 5*RA



No. 221
for 5*FG + 10*RA



No. 241
for RA burs * 12



No. 202
for Root canal files * 18



No. 212
For HP burs * 18



No. 015A
for 18*FG



No. 015C
for 12*FG + 6*RA



No. 015E
for 18*RA



No. 015D
for 20*FG + 10*RA



No. 013E
for 16*FG + 9*RA



No. 022D
for 18*RA



No. 016B
for 8*HP



No. 023
132 * 107 * 65 mm
for FG, RA, HP, Files Mixed

Plastic Blocks, Sterilisable

塑料消毒盒

KIT 套装

Provide kinds of kits for clinic use , which offer the smartest operation ;
Whatever for Ceramics, Zirconia , Lithium Disilicate, Composite,Metal or Alloy ;
From Preparation to Final polishing ;
High Efficiency to shining, No need Glaze after polishing .



BR1901

BR-49, BR-45, BR-46, BR-41, BR-40,
BR-31, BR-28, BR-26, BR-25, BR-23 .

Diamond burs of Ball shape Assorted in FG shank
圆球状金刚石车针混装(高速)



BR1902

BR-X45, BR-X46, BR-X41, BR-X40, BR-X31,
BR-X40S, BR-X31S, BR-X41C, BR-X40C, BR-X31C


Diamond burs of Ball shape Assorted in FG LONG
圆球状金刚石车针混装(高速加长)



BR1903

BR-R46, BR-R41, BR-R40, BR-R31, BR-R26,
BR-R40F, BR-R31F, BR-R26F, BR-R40EF, BR-R26EF


Diamond burs of Ball shape Assorted in RA shank
圆球状金刚石车针混装(低速)



4100

508T-014S, BR-41, 198-016S, 508T-016S, 508T-016F ;
508T-014S, BR-540, 198-021S, 508T-021S, 508T-021F


Safe End Diamond Preparation
前后牙预备(安全端金刚石车针)



4101

TF-21, TF-10, TF-12, TC-10, WR-13
SR-11, TR-11, EX-18, TR-12F, TC-11F

Crown & Bridge Preparation for Anterior
前牙桥体/冠体预备



4102

EX-20S, SR-21, FO-34, TR-62C, TR-21
TR-24, SO-21, EX-18, TR-21F, EX-21F


Crown & Bridge Preparation for Posterior
后牙桥体/冠体预备



18041V

SI-45 x 2 ; SI-46 x 3,
SI-47 x 3, SI-48 x 2 .

Diamond burs of Inverted Cone Assorted in FG shank
常用倒锥金刚石车针混装(高速)



1805TA

TR-11 x 2, TR-12 x 2, TF-11 x 2 ;
TF-12, TC-10, TC-27, TC-28 x 1

Regular Diamond burs of Taper Assorted in FG shank
常用锥形金刚石车针混装 (高速)



1806

TC-11EF, 820-016EF, 820-031F,
820-016F, TC-10F

Interproximal Enamel Reduction (IPR) for Ortho
邻面去釉套装



4103

BR-45, BR-40, BC-32, EX-31
TF-20, SF-41, SF-21, SI-48

Cavity Preparation, Restoratives
备洞



4104

SI-48, DI-42, TF-31, TF-22, SF-21
EX-18, TF-21F, SF-21F, FO-22EF, EX-21EF

Inlay & Onlay Preparation
嵌体/高嵌体预备



4105

BR-49, SF-41, SR-41, TC-26, CR-41
EX-41, FO-40EF / FG1/4, FG1/2, FG 1

Micro Invasive Surgery
微创备牙



6000

BR-45, VP-20, VP-21, VP-22,
SR-41, TR-25, SO-20, TR-13F, CF101F

Porcelain Veneers' Preparation
瓷贴面备牙套装



6002

CD-50F, CD-51F, CD-52F, CD-53F, CD-54F
CD-55F, CD-56F, CD-57F, CD-58F, CD-59F

Children's Teeth Preparation
儿童牙预备



4106

BR-40EF, TC-11EF, TF-11EF, TF-12EF, TR-11EF,
EX-21EF, FO-40EF, FO-21EF, FO-22EF, FO-30EF

Ceramics / Composite Polishing
烤瓷牙/复合树脂抛光



FG008D

BR-49, BR-49, SI-49, SF-49, TF-49

Student or Beginner's Diamond FG 0.8mm
新手入门金刚石车针, 高速



FG009D

001-009M, 001-009M, 010-009M,
109-009M, 170-009M

Student or Beginner's Diamond FG 0.9mm
新手入门金刚石车针, 高速



RA008D

BR-R49, BR-R49, SI-R49, SF-R49, TF-R49

Student or Beginner's Diamond RA 0.8mm
新手入门金刚石车针, 慢速



RA009D

001-R009M, 001-R009M, 010-R009M,
109-R009M, 170-R009M


Student or Beginner's Diamond RA 0.9mm
新手入门金刚石车针, 慢速



TF7100

7004, 7404, 7214, 7714, 7903L

TC Trimming & Finishing burs, FG
钨钢精修车针, FG



1900 T

FG 1, FG 2, FG 3, FG 4, FG 5


Tungsten carbide burs of Ball shape Assorted in FG shank
圆球状钨钢车针混装(高速)



FG008T

FG 1, FG 1, FG 34, FG555, FG698

Tungsten Carbide of FG shank 0.8mm for Students
初学者用的钨钢车针 (头径0.8mm, 高速)



RA008T

RA 1, RA 1, RA 34, RA 555, RA 698

Tungsten Carbide of RA shank 0.8mm for Students
初学者用的钨钢车针 (头径0.8mm, 低速)



7109FW

FL2-Fw, CN1-Fw, RD1-Fw, RD2-Fw, TC1-Fw,

FG White Stone Burs, to Finish Composite
白车石, 高速



7109RW

FL2-Rw, CN1-Rw, RD1-Rw, RD2-Rw, TC1-Rw,

RA White Stone Burs, to Finish Composite
白车石, 低速



7109FG

FL2-Fg, CN1-Fg, RD1-Fg, RD2-Fg, TC1-Fg,

FG Green Stone burs, Finish Ceramics
绿车石, 高速



7109RG

FL2-Rg, CN1-Rg, RD1-Rg, RD2-Rg, TC1-Rg,

RA Green Stone burs, Finish Ceramics
绿车石, 低速



1906

Cr103C x 2, Cr104C, Cr106C, Cr103M, Cr104M, R4AL, R5AL; R4Cb, R5Cb, R4Cr, R5Cr

Composite Finish and Polish RA
树脂成型和抛光,弯机柄



1907

Cr103C, Cr104C, Cr106C, Cr103M, Cr104M, Cr106M, R4AL, R1AL, R5AL; R1Ry, R1G, R1T

Composite Finish and Polish RA
树脂成型和抛光,弯机柄



QR2101

2 step: Sr003F: Pro-Polishing
Sr003UF: Polishing

Composite Polish, RA
树脂抛光,弯机柄



QR2102

2 step: Sr304F: Pro-Polishing
Sr304UF, Sr003UF: Polishing

Composite Polish, RA
树脂抛光,弯机柄



QR2103

3 step: Sr304CM: Smooth;
Sr003F: Pro-Polishing; Sr003UF: Polishing

Composite Polish, RA
树脂抛光,弯机柄



QR2104

3 step: Sr003M: Smooth;
Sr003F: Pro-Polishing; Sr003UF x 2: Polishing

Composite Polish, RA
树脂抛光,弯机柄



CP 6011

2 Step: Sr 204 M x 3: Pro-polish;
Sr 204 F x 3: Polish

Composite, Ceramic, Zirconia Polishers . RA
树脂、烤瓷、氧化锆抛光



CP 6001

3 Step: Sr 204 C x 2: Smooth;
Sr 204 M x 2: Pro-polish; Fine: Sr 204 F x 2: Polish

Composite, Ceramic, Zirconia Polishers . RA
树脂、烤瓷、氧化锆抛光



7115

Smooth: DCR101, DCR102, DCR122, DCR105
Pro-polish: Sr101F, Sr103F, Sr105F, Sr106F
Polish: Sr101UF, Sr103UF, Sr105UF, Sr106UF

Diamond stones & polishers for Zirconia / Lithium Disilicate, 氧化锆, 玻璃陶瓷研磨抛光



1904

X151, X1151, EC25, EC28, FG X6,
EX-28, BR-X31, TR-X12, TR-X15, TF-X13

Extraction Burs
拔阻生齿车针



1905

FG X2, FG X4, FG X6, FG X558, FG X701,
BR-X41, BR-X40, BR-X31, TF-X12, TR-X12

Oral Surgery burs FG Long
手术车针(高速加长)



4110V

Smooth: Rw01, Rw02 x 2
Pro-polish: Cr101M, Cr103M, Cr104M, Cr106M,
Polish: Cr101F, Cr103F, Cr104F, Cr106F

Composite Finish and Polish RA
树脂成型和抛光,弯机柄



4112

Smooth: Cr101C, Cr103C, Cr104C, Cr106C
Pro-polish: Cr101M, Cr103M, Cr104M, Cr106M
Polish: Cr101F, Cr103F, Cr104F, Cr106F

Porcelain / Natural Teeth Polishing RA
烤瓷牙,自然牙抛光,弯机柄



4114

Smooth: Ar101C, Ar103C, Ar104C, Ar106C
Pro-polish: Ar101M, Ar103M, Ar104M, Ar106M
Polish: Ar101F, Ar103F, Ar104F, Ar106F

Amalgam / Metal Polishing RA
银汞合金,金属抛光,弯机柄



4115

Smooth: Sr103M, Sr105M, Sr106M; Pro-polish: Sr103F, Sr105F, Sr106F
Polish: Sr101EF, Sr103EF, Sr105EF, Sr106EF,
Final Polish: Sr101UF, Sr103UF

Diamond Polishers for Zirconia / Lithium Disilicate
氧化锆,玻璃陶瓷抛光,弯机柄

Crown & Bridge Preparation for Anterior
前牙桥体/冠体预备



4101A

TF-21, TF-10, TF-12 x 2, WR-13, TC-21,
TR-11, TC-10/EX-18, SR-11/ TR-12F, TC-11F




Crown & Bridge Preparation for Posterior
后牙桥体/冠体预备



4102A

EX-11, EX-20S, TR-22S, BR-31C, TR-62C, FO-34/
TR-24, TR-21/EX-18, SO-21/EX-21F, TR-21F




Cavity Preparation, Restoratives
备洞



4103A

BR-45, BR-40, BR-46, BR-31, BR-41/BC-32,
SI-47, SI-48/SF-41, SF-21/EX-31, TF-20




Inlay/Onlay Preparation
嵌体/高嵌体预备



4104A

DI-42, SI-48, TF-22, TF-31 / TF-21F, SF-21F /
TF-21, SF-21, EX-18, BC-32 / EX-21EF, FO-22EF




Porcelain Veneer's Preparation
瓷贴面备牙




6000A

SR-41, BR-45 / VP-20 x 2, VP-21 x 2, VP-22 x 2/
SO-20, TR-25, TR-13F, TC-28EF



Trimming & Finishing Tungsten
Carbide burs FG 钨钢精修车针



6003A

7004, 7005, 7006, 7404, 7406, 7114,
7214, 7714, 7901, 7902, 7903, 7903L




Zirconia, Lithium Disilicate Polishing kit, RA
氧化锆, 玻璃陶瓷抛光(低速)




7115A

Smooth : DCR101 , DCR102, DCR122, DCR105 ; Pro-polish: Sr101F, Sr103F, Sr105F, Sr106F
Polish: Sr101UF, Sr103UF, Sr105UF, Sr106UF




Zirconia, Lithium Disilicate High Gloss polishing kit, RA
氧化锆抛光(低速)



4115A

Smooth : Sr103M, Sr105M, Sr106M ; Pro-polish : Sr103F, Sr105F, Sr106F
Polish: Sr101EF, Sr103EF, Sr105EF, Sr106EF, Final Polish: Sr101UF, Sr103UF



Composite or filling material Polisher RA
美学树脂、填充材料抛光



CP6001A

3 Step: Sr 204 C x 2 : Smooth;
Sr 204 M x 2 : Pro-polish ; Fine : Sr 204 F x 2 : Polish



Composite Polishing kit , RA
树脂抛光套装(低速)



4110A2

Smooth: FL2-Rwx 2 , CN1-Rw x 2 ; Pro-polish: Cr101M, Cr103M, Cr104M, Cr106M,
Polish: Cr101F, Cr103F, Cr104F, Cr106F



Ceramics / Porcelain Teeth Polish , RA
烤瓷牙抛光(低速)



4112A

Smooth: Cr101C, Cr103C, Cr104C, Cr106C
Pro-polish : Cr101M, Cr103M, Cr104M, Cr106M
Polish: Cr101F, Cr103F, Cr104F, Cr106F



Amalgam, Alloy, Metal Polish , RA
银汞合金抛光(低速)




4114A

Smooth : Ar101C, Ar103C, Ar104C, Ar106C
Pro-polish : Ar101M, Ar103M, Ar104M, Ar106M
Polish: Ar101F, Ar103F, Ar104F, Ar106F



Micro Invasive Surgery Burs
微创手术用车针



4105A

BR-49 x 2, SR-41 x 2, SF-41, EX-41, TC-26,
CR-41, FG1/4, FG1/2, FG 1, FO-40EF

Extra Fine Diamond burs for polishing
特细粒度金刚砂抛光车针,高速



4106A

CR-11EF, BR-40EF, TC-27EF, TC-11EF, TF-12EF, TF-11EF,
FO-40EF, TR-11EF, FO-22EF, FO-21EF, FO-30EF, EX-21EF

Composite, Ceramics Finishing kit, FG
树脂成型(高速)



0302A

FL2-Fw x 3 ; CN1-Fw x 3
FL2-Fg x 3 ; CN1-Fg x 3

Composite, Ceramics Finishing kit, RA
树脂成型(低速)



0306A

FL2-Rw x 3 ; CN1-Rw x 3
FL2-Rg x 3 ; CN1-Rg x 3

White Stone Burs kit, FG
白车石(高速)



302F1

FL2-Fw x 2, CN1-Fw x 2, RD1-Fw x 2
RD2-Fw x 2, TC1-Fw x 2, IC2-Fw x 2

White Stone Burs kit, RA
白车石(低速)



306R1

FL2-Rw x 2, FL3-Rw x 2, CN1-Rw x 2,
RD1-Rw x 2, RD2-Rw x 2, TC1-Rw x 2.

HP Kit to adjust Zirconia
氧化锆修整椅旁套装(HP)



Z31

ZF33, DY001EF to Smoothing
D203F for Pre-polishing ; D203EF for Polishing

HP Kit to adjust Lithium Disilicate
玻璃陶瓷修整椅旁套装(HP)



L31

ZF33, DY001EF to Smoothing
D203F for Pre-polishing ; D203UF for Polishing

Samples Collection For Clinic Items

样品册



Flip-Cart Plastic Pages ;
Smart, Easily fill all burs ;
Economic and Recyclable ;
Nice-looking for showing to Customers;
Convenient to take to anywhere .



No. 170B5

Empty Collection Wallet
样品本(不带车针)



4 Colors Available ,
100 holes per page,
5 pages per set,
Smart to FG or RA Burs



No. 1701D

Diamond Burs Collection Wallet
金刚石车针样品本



More than 300 specification ,
Diamond burs ,
FG, FG XL, FG 2XL,
FG Short, RA, RA Long



No. 1702T

Tungsten Carbide Collection Wallet
钨钢车针样品本



Nearly 200 specification ,
TC burs
FG, FG XL, FG Short, RA, RA Long
T & F Gold, Endo z etc.



No. 1703P

Polishing burs Collection Wallet
抛光器械样品本



Nearly 150 specifications ,
Silicon Rubber burs,
Diamond Abrasives & polishers,
Polishing brushes etc




Four Colors Available

Brown

Black

Blue

Yellow





Laboratory Part 技工类

Page 33-66 第33-66页



Diamond bur 金刚石磨头

Page 34-35 第34-35页

Fully-sintered Diamond bur 全砂烧结磨头

Page 36-37 第36-37页

Sintered Diamond Disc Diamond Disc

烧结切片、金刚石切片

Page 38-40 第38-40页

CAD / CAM Milling Bur 铣刀

Page 41-43 第41-43页

Tungsten Carbide bur 钨钢车针

Page 44-47 第44-47页

Diamond Stones & Polisher 陶瓷(橡胶) 结合剂金刚石磨头 (全瓷、全锆磨头)

Page 48-52 第48-52页

Silicon Polisher 硅胶抛光磨头

Page 53-54 第53-54页

Mounted stone, Sandpaper 车石、砂纸棒

Page 55 第55页

Polishing Brushes, Felt burs Mandrels

抛光刷、毛毡磨头、夹针

Page 56-58 第56-58页

Hard Z Cutter etc . 硬锆磨头等

Page 59 第59页

Kit 套装

Page 60-62 第60-62页

Samples Collection 样品册

Page 63 第63页

FlexSoft Partial Denture 隐形义齿胶

Page 64 第64页

Diamond Bur HP (Electroplated) 金刚石磨头 (电镀)



Sharp and Durable,
Grind or polish Ceramics
锋利、耐用
磨削平稳、不偏摆
不留黑印、适合车瓷、



Grit


- Black / Super Coarse
- Green / Coarse
- Blue / Medium
- Red / Fine

Fine
Red
细砂
红圈




编号 No. This grit can grind Zirconia 可用来车氧化锆

Medium
Blue
标准砂
蓝圈




Round / 圆球

Inverted Cone




Wheel / 平轮

Cylinder / 平柱



Round Cylinder / 圆头柱


Needle / 尖针



Taper / 平锥

Round Taper / 圆头锥

Inverted Cone




Transparent Showing plates 展示台 120x120x8mm



No. 1216

7 x 7 holes
1.60mm, for FG burs



No. 1223

7 x 7 holes
2.35mm, for RA or HP burs.


Fully-Sintered Diamond Bur . HP

全砂烧结磨头



Extremely Durable,
Grind or polish smoothly
Especially for Metal

从里到外都是砂
平稳、耐用，适合车金




Grit

- Black/Super Coarse
- Green/Coarse
- Blue/Medium
- Red/Fine
- Yellow/Extra Fine

Super
Coarse
Black
特粗砂
黑圈



编号 No. D060S D070S K080S K100S J070S J100S A050S A060S C050S C060S



C070S M050S M070S B050S B070S L050S L060S L070S G060S F050S F070S

Coarse
Green
粗砂
绿圈



R050S R070S D050C D060C J060C J070C A037C A050C A060C C037C






C050C C060C B037C B050C B060C L037C L050C L060C N050C G050C F037C

Medium
Blue
标准砂
蓝圈



F050C R037C R050C D031 D037 D050 D060 K031 K037 K050




Sintered Diamond Disc

烧结切割片

Grinding or Trimming of sintered Zirconia or any other Extra hard workpiece

用于切割烧结后的氧化锆或者各种超硬材质工件。



The galvanically sintered edge of the diamond discs offers excellent cutting performance.

提供三种不同粒度的砂粒，供您试用。



70 μm

Item No.



S394/190/0.25mm

S394/220/0.25mm

S394/190/0.30mm

S394/220/0.30mm

Size (mm)

19 x 0.25mm

22 x 0.25mm

19 x 0.30mm

22 x 0.30mm



60 μm

Item No.



S395/190/0.25mm

S395/220/0.25mm

S395/190/0.30mm

S395/220/0.30mm


Size (mm)

19 x 0.25mm

22 x 0.25mm

19 x 0.30mm

22 x 0.30mm



40 μm

Item No.



S396/190/0.25mm

S396/220/0.25mm

S396/190/0.30mm

S396/220/0.30mm

Size (mm)

19 x 0.25mm

22 x 0.25mm

19 x 0.30mm


22 x 0.30mm

Diamond Disc

金刚石切片



How to order these Discs:



Extra Thin, 0.15mm
Double sided
Extra Fine
超薄, 双面砂



编号 No.	A19D15	A22D15	Am13D15	Am16D15	Am19D15	Am22D15
--------	--------	--------	---------	---------	---------	---------

Thin, 0.20 x 19mm
Double sided
Fine Grit
较薄, 双面砂



编号 No.	A19D20	B19D20	C19D20	N19D20	Am19D20	Bm19D20	Em19D20	Hm19D20
--------	--------	--------	--------	--------	---------	---------	---------	---------

Thin, 0.20 x 22mm
Double sided
Fine Grit
较薄, 双面砂




编号 No.	Am22D20	Bm22D20	Cm22D20	Em22D20	Im22D20	Jm22D20	Hm22D20
--------	---------	---------	---------	---------	---------	---------	---------




A22D20	B22D20	C22D20	E22D20	F22D20	H22D20	J22D20	N22D20
--------	--------	--------	--------	--------	--------	--------	--------

Hard, 0.3 - 0.4mm
Double sided
Medium Grit
厚硬, 双面砂



编号 No.	Am19D30	Am22D30	Am22D40	A19D40	A22D40
	19.0 x 0.3 mm	22.0 x 0.3 mm	22.0 x 0.4 mm	19.0 x 0.4 mm	22.0 x 0.4 mm




Single-sided
仅单面镀砂



编号 No.	A22S20	Am22S20	Am22S30
	22.0 x 0.2 mm, Fine	22.0 x 0.2 mm	22.0 x 0.3 mm, Medium


Fast Cutting
Plaster Disc ,
0.3mm,
Double sided
石膏模切片
平整、无偏摆



30mm



35mm



35mm



35mm


编号 No.

C30D30

A35D30

C35D30

J35D30



40mm

A40D30



40mm

C40D30




40mm

J40D30




45mm

Cb45D30



45mm

Ab45D30



20mm

85mm

Ch85D30

Tungsten Carbide Burs . HP 钨钢车针



Shape

- S1, ISO 001, Round
- S2, ISO 002, Inverted cone
- S21, S31, ISO 107, Cylinder
- S33, ISO 168, Taper
- S31R, ISO 137, Round, Cylinder

										
S1	001 0.5	001 0.6	001 0.7	001 0.8	001 0.9	001 1.0	001 1.2	001 1.4	001 1.6	001 1.8
球型/Round										
编号 No.	1005	1006	1007	1008	1009	1010	1012	1014	1016	1018
参考号 Ref.	HP1/4	HP1/2	HP1s	HP 1	HP 1b	HP 2	HP 3	HP 4	HP 5	HP 6

										
001 2.1	001 2.3	001 2.5	001 2.7	001 3.1	S2	010 0.8	010 1.0	010 1.2	010 1.4	010 1.6
					倒锥/Inv. Cone					
1021	1023	1025	1027	1031		2008	2010	2012	2014	2016
HP 7	HP 8	HP 9	HP 10	HP 12		HP 34	HP 35	HP 36	HP 37	HP 38


										
010 1.8	010 2.1	010 2.3	S21	107 0.8	107 0.9	107 1.0	107 1.2	107 1.4	107 1.6	107 1.8
			圆柱/Cylinder							
2018	2021	2023		21008	21009	21010	21012	21014	21016	21018
HP 39	HP 40	HP 41		HP 55	HP 56	HP 57	HP 58	HP 59	HP 60	HP 61

										
107 2.1	107 2.3	S31	107 0.8	107 0.9	107 1.0	107 1.2	107 1.4	107 1.6	107 1.8	107 2.1
		圆柱/Cylinder								
21021	21023		31008	31009	31010	31012	31014	31016	31018	31021
HP 62	HP 63		HP 555	HP 556	HP 557	HP 558	HP 559	HP 560	HP 561	HP 562

										
107 2.3	S33	168 0.8	168 0.9	168 1.0	168 1.2	168 1.6	168 2.1	168 2.3	S33L	171 0.9
	平锥/Taper								平锥加长/Taper L	
31023		33008	33009	33010	33012	33016	33021	33023		33L009
HP 563		HP 698	HP 699	HP 700	HP 701	HP 702	HP 703	HP 704		HP 699 L

										
171 1.0	171 1.2	171 1.6	171 2.1	S31R						
				圆头柱 Round Cyl.						
33L010	33L012	33L016	33L021		31R010	31R012	31R014			
HP 700 L	HP 701 L	HP 702 L	HP 703 L		HP 1557	HP 1558	HP 1559			

Tungsten Carbide Cutter . HP 钨钢磨头



Super Fine Cross Cut (SFE) 特细双纹 黄圈 ISO 500.104.110.

145 023 16.0	274 040 11.5	274 060 14.0	257 023 7.0	198 016 8.0	201 023 14.0	201 040 14.0	201 060 15.0	161 060 14.0	237 023 5.0
编 号 No. C023SFE	F040SFE	F060SFE	G023SFE	L016SFE	L023SFE	L040SFE	L060SFE	M060SFE	P023SFE
参考号 Ref. 0230.023	7230.040	5430.060	1030.023	1230.016	0130.023	5730.040	5530.060	5930.060	7730.023

Fine Cross Cut (FE) 较细双纹 红圈 ISO 500.104.140.

289 023 8.0	292 023 16.0	145 023 16.0	143 060 13.0	274 060 14.0	257 040 9.0	257 060 12.0	198 016 8.0	201 023 14.0
Q023SFE	Q023SFE	C023FE	C060FE	F060FE	G040FE	G060FE	L016FE	L023FE
0830.023	0430.023	0220.023	5820.060	5420.060	6120.040	5320.060	1220.016	0120.023

Medium Cross Cut (CE) 标准双纹 蓝圈 ISO 500.104.190.

201 040 14.0	201 060 14.0	161 060 15.0	237 023 5.0	237 060 12.0	292 023 16.0	116 023 16.0	145 023 16.0	144 060 12.0
L040FE	L060FE	M060FE	P023FE	P060FE	Q023FE	A023CE	C023CE	C060CE
5720.040	5520.060	5920.060	7720.023	5120.060	0420.023	0310.023	0210.023	5810.060


001 023 2.3	001 030 3.0	001 040 4.0	001 050 5.0	274 040 11.5	274 060 14.0	257 023 7.0	257 040 9.0	257 060 12.0	201 023 14.0
D023CE	D030CE	D040CE	D050CE	F040CE	F060CE	G023CE	G040CE	G060CE	L023CE
7110.023	7110.030	7110.040	7110.050	7210.040	5410.060	1010.023	6110.040	5310.060	0110.023

								Medium Plain Cut (CN) 标准单纹 蓝圈 ISO 500.104.175.	
201 040 14.0	201 060 15.0	167 023 14.50	167 040 14.5	161 060 15.0	237 023 5.0	237 060 11.0	292 023 16.0	001 030 3.0	
L040CE	L060CE	M023CE	M040CE	M060CE	P023CE	P060CE	Q023CE	D030CN	
5710.040	5510.060	—	—	5910.060	7710.023	5110.060	0410.023	7170.030	

							Coarse Cross Cut (SCE) 粗糙双纹 绿圈 ISO 500.104.220			
001 040 4.0	001 050 5.0	201 040 14.0	201 060 15.0	161 060 14.0	257 060 12.0	274 060 14.0	201 060 14.0		274 060 14.0	257 060 12.0
D040CN	D050CN	L040CN	L060CN	M060CN	G060CN	F060CN	C060SCE		F060SCE	G060SCE
7170.040	7170.050	5770.040	5570.060	5970.060	5370.060	5470.060	5850.060		5450.060	5350.060

			Super Coarse Cross Cut (XCE) 特粗双纹 黑圈 ISO 500.104.223			Straight Edge with Spiral Cut (QSCN) 直刃螺旋纹 蓝圈 ISO 500.104.176					
201 040 14.0	201 060 14.0	237 060 11.0			263 070 14.0	274 060 15.0		201 023 14.5	200 045 13.0	263 060 12.0	
L040SCE	L060SCE	P060SCE	L070XCE			F060QSCN		L023QSCN	L045QSCN	LB060QSCN	
5750.040	5550.060	5150.060	6055.070			7275.060		0175.023	5675.045	5275.060	

Straight Edge with Cross Cut (GTi) 直刃双纹 黑圈 ISO 500.104.194				Spiral Cut Fine (QF) 较细螺旋纹 红圈 ISO 500.104.134			
							
141 023 8.0	194 023 8.0	201 040 14.0	289 023 8.0	274 040 11.5	274 060 15.0	201 023 14.5	201 040 13.0
CS023GTi	LS023GTi	L040GTi	QS023GTi	F040QF	F060QF	L023QF	L040QF
1726.023	0726.023	5726.040	0826.023	7224.040	5424.060	0124.023	5724.040

Spiral Cut Extra Fine (QSF) 特细螺旋纹 黄圈 ISO 500.104.137									
		141 023 8.0	274 060 14.0	198 023 8.0	201 040 14.0	289 023 8.0	118 060 3.5		467 023 9.0
		CS023QSF	F060QSF	LS023QSF	L040QSF	QS023QSF	ES801		AL023S4
		1725.023	5425.060	0725.023	5725.040	0825.023	118.060		467.023

Tungsten Carbide Burs for Milling Machines


钨钢车针(精密研磨仪用)



Precisely match Various System

Extra - long working life & high cutting capacity ;
Used for trimming, polishing implant teeth shoulders ;
Milling precious or Non-precious metals or Titanium .

Medium
Cross Cut
(CM, LM)
标准双纹
ISO 500.103.190




Fine
Cross Cut
(CF, LF)
较细单纹
ISO 500.103.135



Coarse
Cross Cut
(CC, LC)
粗糙双纹
ISO 500.103.220



Max Mill
(CW, LW)
蜡型针
ISO 500.103.217



GZH

PRE-SINTERED ZIRCONIA POLISHING BURS

预烧结氧化锆
抛光磨头

- Polishing pre-sintered zirconia simplifies finishing after final sintering
抛光预烧结氧化锆可简化最终烧结后的加工步骤
- No contamination of zirconia, as no colour pigments are used
由于未使用颜色颜料，氧化锆无污染
- Developed for the lower strength of pre-sintered zirconia
为预烧结氧化锆的较低强度而开发的



Set GZHS-01

Pre-sintered Zirconia Finishing

预烧结氧化锆
打磨

NO.	GZH243-5C	GZH275-6C	GZH040-15C
D x H mm	5.5 x 13	6.5 x 13	15 x 0.4

HP Shank, Coarse
HP 长柄, 粗糙

Pre-sintered Zirconia Polishing 2-in-1 Burs

预烧结氧化锆
抛光
二合一磨头

NO.	GZH507-10MF	GZH243-10MF	GZH243-6MF	GZH507-7MF	GZH204-10MF
D x H mm	10 x 22.5	10 x 20	6 x 16	7 x 19	10 x 24

HP Shank, Medium + Fine
HP 长柄, 中等 + 精细

1 STEP
MEDIUM



2 STEP
FINE



10 Pieces

DuCoBur[®] 陶瓷结合剂金刚石磨头 (全瓷、全锆磨头)

Diamond Stones for Zirconia or AllCeramics
(Ceramic Pregnant Diamond Grinding Burs)



Widely Used for Zirconia, Lithium Disilicate, Ceramics etc .
COMFORTABLE at the minimal pressure;
HEATLESS no need water cooling ;
MORE DURABLE under our suggested PRM

不产热、无需水冷 更耐磨、经济实惠 更锋利、手感好 不绷瓷、更轻松 型号齐、用途广、啥都磨得动

Medium Blue
蓝色、中等粒度
1

Item No.	DC009	DC010	DC010B	DC011	DC014	DC015	DC016	DC017	DC003A
Size (mm)	4 x 4	4 x 7	4 x 11	3 x 7	2.5 x 4.5	3.5 x 10	2.5 x 7	3.5 x 11	6.5 x 2
DC001	DC002	DC002B	DC012	DC012B	DC013	DC004A			
4 x 13	5 x 13	6.5 x 13	13 x 2	13 x 3.5	13 x 2	16 x 2/0.7			
DC005A	DC004	DC005	DC006	DC007					
16 x 2	22 x 2/0.7	22 x 2	22 x 4.5	15 x 3.5					

Medium Black
黑色、中等粒度
1

Item No.	DX001	DX002	DX002B	DX012	DX012B	DX013	DX004A
Size (mm)	4 x 13	5 x 13	6.5 x 13	13 x 2	13 x 3.5	13 x 2	16 x 2/0.7
DX005A	DX004	DX005	DX006	DX007			
16 x 2	22 x 2/0.7	22 x 2	22 x 4.5	15 x 3.5			



Medium
Yellow

黄色、中等粒度



1



Medium
Yellow



1

黄色、中等粒度















Item No.	DY001	DY002	DY002B	DY012	DY012B	DY013	DY004A
Size (mm)	4 x 13	5 x 13	6.5 x 13	13 x 2	13 x 3.5	13 x 2	16 x 2/0.7



DY005A

16 x 2



DY004

22 x 2/0.7



DY005

22 x 2



DY006

22 x 4.5



DY007

15 x 3.5



Medium
Green

绿色、中等粒度



1

 Medium Green 绿色、中等粒度									
 1									
Item No.	DG009	DG010	DG010B	DG011	DG014	DG015	DG016	DG017	DG003A
Size (mm)	4 x 4	4 x 7	4 x 11	3 x 7	2.5 x 4.5	3.5 x 10	2.5 x 7	3.5 x 11	6.5 x 2



DG001

4 x 13



DG002

5 x 13



DG002B

6.5 x 13



DG012

13 x 2



DG012B

13 x 3.5



DG013

13 x 2



DG004A

16 x 2/0.7



DG005A

16 x 2



DG004

22 x 2/0.7



DG005

22 x 2



DG006

22 x 4.5




DG007

15 x 3.5

Extra Fine
Yellow
Gold shank


Suitable to adjusting Lithium Disilicate, All Ceramics, Zirconia of Occlusal, Lingual, Edge, Intracoronar, Proximal etc. Dentists will feel very convenient if they can prepare a few pieces beside the chairs.

专门针对诊所的全瓷、玻璃陶瓷打磨抛光、氧化锆牙调咬合、发育沟、窝沟、边缘、内冠、邻面、舌侧等。黄色、特细粒度。粒度细、不绷瓷。



DuCoBur® 豆酷巴

Diamond Polishers for Zirconia or AllCeramics
(Rubber Pregnant Diamond Polishing Burs)



Widely Used for Zirconia, Lithium Disilicate, Ceramics etc.

COMFORTABLE at the minimal pressure;
HEATLESS no need water cooling ;
MORE DURABLE under our suggested PRM

不产热、无需水冷 更耐磨、经济实惠 手感好、很服帖 抛就亮、更轻松 型号齐、用途广、啥都抛得亮

Super Coarse
Black
黑色、特粗粒度



Item No.

D1S

D2AS

D2S

D3S

D3BS

Size (mm)

26 x 2

20 x 2

26 x 2

22 x 3

26 x 2

Coarse
Green
绿色、粗糙粒度



Item No.

D1C

D2AC

D2C

D3C

D3BC

Size (mm)

26 x 2

20 x 2

26 x 2

22 x 3

26 x 2

 Medium
Blue
蓝色、中等粒度
 1



Item No.

D1M

D2AM

D2M

D3M

D3BM

Size (mm)

26 x 2

20 x 2

26 x 2

22 x 3

26 x 2



D201M

5 x 16



D203M

6 x 16



D205M

12 x 2



D205BM

12 x 4



D207M

12 x 2.5



D13M

14 x 2



D17M

17 x 2.2



D26M

26 x 2.2

 Fine
Red
红色、较细粒度
 1



D1F

26 x 2



D2AF

20 x 2



D2F

26 x 2



D3F

22 x 3



D3BF

26 x 2



D201F

5 x 16



D203F

6 x 16



D205F

12 x 2



D205BF

12 x 4



D207F

12 x 2.5



D13F

14 x 2



D17F

17 x 2.2



D26F

26 x 2.2

Extra Fine
Yellow
黄色、特细粒度



1



D1EF



D2AEF



D2EF



D3EF



D3BEF

Item No.

Size (mm)

26 x 2

20 x 2

26 x 2

22 x 3

26 x 2



D201EF

5 x 16



D203EF

6 x 16



D205EF

12 x 2



D205BEF

12 x 4



D207EF

12 x 2.5



D13EF

14 x 2



D17EF

17 x 2.2



D26EF

26 x 2.2

Ultra Fine
Grey
灰白色、特细粒度



1



D1UF



D2AUF



D2UF



D3UF



D3BUF

Item No.

Size (mm)

26 x 2

20 x 2

26 x 2

22 x 3

26 x 2



D201UF

5 x 16



D203UF

6 x 16



D205UF

12 x 2



D205BUF

12 x 4



D207UF

12 x 2.5




D13UF

14 x 2



D17UF

17 x 2.2



D26UF

26 x 2.2

**Silicon Rubber Wheels ,
Points & Burs
For Polish Ceramics
胶轮、胶粒、硅胶磨头
(烤瓷牙抛光)**




White, Coarse , to adjust and remove Porcelain
白色, 粗糙, 用于烤瓷工件的修整。
Pink, Medium , to Smooth or polish Porcelain,
粉色, 中等, 用于烤瓷毛面平整和预抛光。
Blue, Fine , to high-lustre polishing Porcelain,
蓝色, 较细, 用于烤瓷工件的细抛光。


Coarse
White
白色、粗糙



Item No.



C 104 C



C 107 C



C 101 C



C203C



C201C



C205C



C207C

Size (mm)

22 x 3

22 x 3

6 x 22

6 x 13

5 x 13

12 x 2

13 x 2

Medium
Pink
粉色、中等粒度



Item No.



C 104 M



C 107 M



C 101 M



C203M



C201M



C205M



C207M

Size (mm)

22 x 3

22 x 3

6 x 22

6 x 13

5 x 13


12 x 2

13 x 2

Fine
Blue
蓝色、较细



Item No.



C 104 F



C 107 F



C 101 F



C203F



C201F



C205F



C207F

Size (mm)

22 x 3

22 x 3

6 x 22

6 x 13

5 x 13

12 x 2

13 x 2

**Rubber Wheels , Points & Burs
For Polish Metal or Alloy
胶轮、胶粒、橡胶磨头
(金属、合金抛光)**



Black, Coarse , to adjust Cobalt-Chromium
Alloy and hard metal
黑色, 粗糙, 用于磨削钴铬合金和硬金属。
Brown, Medium , to remove metal and
semi-precious metal
棕色, 中等, 用于烤瓷毛面平整和预抛光。
Green, Fine , to fast smooth scratch or polish metal,
绿色, 较细, 用于快速去除划痕和抛光金属。

Coarse
Black
黑色、较粗



Item No.



A102 C



A103C



A104 C



A203C



A201C



A206C



A207C

Size (mm)

22 x 3

22 x 3

6 x 22


4 x 12

5 x 13

12 x 2

13 x 2

Medium
Brown
棕色、中等粒度
100



Item No.	A102 M	A103 M	A104 M	A203 M	A201 M	A206 M	A207 M
----------	--------	--------	--------	--------	--------	--------	--------

Size (mm)	22 x 3	22 x 3	6 x 22	4 x 12	5 x 13	12 x 2	13 x 2
-----------	--------	--------	--------	--------	--------	--------	--------

Fine
Green
绿色、较细
100



Item No.	A102 F	A103 F	A104 F	A203 F	A201 F	A206 F	A207 F
----------	--------	--------	--------	--------	--------	--------	--------

Size (mm)	22 x 3	22 x 3	6 x 22	4 x 12	5 x 13	12 x 2	13 x 2
-----------	--------	--------	--------	--------	--------	--------	--------

Silicon Polishing Burs
For Composites, Resins
硅胶磨头
(树脂胶托抛光)



Green, Coarse, to reduce, finish Acrylic, Resin.
绿色, 粗糙, 用于修整和成型树脂胶托。
Gray, Medium, to Smooth Acrylic, Resin.
灰黑色, 中等, 用于平整和预抛光树脂胶托。
Yellow, Fine, to final polish Acrylic, Resin.
黄色, 较细, 用于树脂胶托的细抛光。

Coarse
Green
绿色、粗糙
100



Item No.	Rh 101 C	Rh 103 C	Rh 106 C	R103 C	R106 C	R209 C
----------	----------	----------	----------	--------	--------	--------

Size (mm)	13 x 20	10 x 24	10 x 20	10 x 24	10 x 20	6 x 13
-----------	---------	---------	---------	---------	---------	--------

Medium
Gray
灰黑、中等粒度
100



Item No.	Rh 101 M	Rh 103 M	Rh 106 M	R103 M	R106 M	R209 M
----------	----------	----------	----------	--------	--------	--------

Size (mm)	13 x 20	10 x 24	10 x 20	10 x 24	10 x 20	6 x 13
-----------	---------	---------	---------	---------	---------	--------

Fine
Yellow
黄色、较细
100



Item No.	Rh 101 F	Rh 103 F	Rh 106 F	R103 F	R106 F	R209 F	R202 F
----------	----------	----------	----------	--------	--------	--------	--------

Size (mm)	13 x 20	10 x 24	10 x 20	10 x 24	10 x 20	6 x 13	4 x 12
-----------	---------	---------	---------	---------	---------	--------	--------

Green Stones
Contour Ceramic



绿车石
调烤瓷



Pink Stones
for Metal
粉车石
调金属,合金



White Stone
for Metal
Composite,
白车石
调金属,树脂



Mounted Sandpaper 砂纸棒



Polishing Brushes 抛光刷

				
百洁布 暗红,中粗	百洁布 红,硬	百洁布 棕,中粗	百洁布 紫,软	百洁布 普通
01MT25	01Mp25	01My25	01Ms25	01M22
Red.Soft	Red.coarse	Brown. Medium	Purple.fine	Blue

						
碳化硅 磨料丝	氧化铝 磨料丝	曲钢丝	曲铜丝	黑猪鬃	灰马毛	灰羊毛
01SiC22	01Al22	01A22	01B22	01D22	01E22	01F22
Silicon Carbide	Al. Oxide	Crimped steel wire	Crimped brass wire	Chungking Bristle	Grey Pony hair	Grey goat hair

					
白羊毛	灰羊毛 +超纤皮	黑猪鬃 19mm	灰马毛	灰羊毛	白羊毛
01G22	01Fr22	01D19	01E19	01F19	01G19
White goat hair	Goat hair / chamois	Chungking Bristle	Grey Pony hair	Grey goat hair	White goat hair

							
黄棉布	白棉布	橘红 羊皮	双面绒 羊皮	咖啡 超纤	白厚 超纤	羊绒	白棉线
01P22	01Q22	01R22	01Ry22	01R22.2	01Ry22.2	01S22	01T22
Yellow cloth	White cloth	Chamois leather	Sheep Leather	Brown Micro-Fiber	White Microfiber	Cashmere	Cotton yarn



						
黑猪鬃	灰马毛	白羊毛	曲钢丝	曲铜丝	黑猪鬃	灰羊毛
02D22	02E22	02G22	04A12	04B12	04D12	04F12
Chungking Bristle	Pony hair	White goat hair	Steel	Brass	Bristle	Goat Hair

						
碳化硅 磨料丝	氧化铝 磨料丝	曲钢丝	曲铜丝	黑猪鬃	灰马毛	灰羊毛
05SiC05	05Al05	05A05	05B05	05D05	05E05	05F05
SiC	Al. Oxide	Steel	Brass	Bristle	Pony Hair	Goat Hair

Polishing brush & Wheel in Bristle, Cotton Cloth 抛光刷



These brushes are made of bristle hair or White, yellow pure cotton cloth, assembled on polishing latch to finish or polish Plaster Moulds, Acrylic
台式抛光机配用的抛光刷/布轮, 用于抛光塑料, 树脂, 模型等。



Item	2A	2B	2C	6A	6B	6C	J03	J02	J01	P80	03B48
D	79	73	68	64	60	54	70	58	51	80	48
d	48	48	48	35	35	35	35	30	25	40	25
e	6.5	6.5	6.5	6.5	6.5	6.5	8.5/5	8/5	7/5	6.5	6.35
B	16	16	16	16	16	16	24	16.5	15	36	9
Rows	2	2	2	2	2	2	3	2	1	4	1



Converging




Straight


Felt Bobs & Burs 毛毡磨头




03N25	03Ns25	3P	2P	1P	1R
25 x 5 mm, 硬质	25 x 5 mm, 软质	15 x 25 mm	20 x 40 mm	22 x 50 mm	22 x 50 mm



Mandrels Stainless Reinforced 不锈钢夹针



Mandrels Steel 普通夹针



Hard Z Cutter 硬锆磨头


Made of Nano-grade Zirconia, used to grinding, finishing or polishing Ceramics and Acrylic, Resin, Zirconia etc. Also applied to Soft materials such as wood, Ceramic or stones etc

用于打磨,成型,抛光烤瓷,树脂,氧化锆(烧结前)。特殊的纳米级氧化锆胚体精制而成。



Sharpening Stones 磨刀石

Adjust to minimum speed;
Touch the stones softly.
Using water will be better



Cleaning Brushes By Hand

清洁刷(手用)



Stainless Steel Disc 不锈钢小锯片

for Soft Plastic, Acrylic, Wood etc.




Stone Disc 刚玉车片



3Z09.1**Zirconia Luster Polishing Kit**
氧化锆抛光套装


- + 2 pieces of Diamond Stones DG001, DG012
- + 2 pieces of Diamond Polishers, Blue, D203M, D3M
- + 2 pieces of Diamond Polishers, Yellow, D203EF, D3EF

3Z00.2**Grinding and polishing Kit of Zirconia**
氧化锆打磨抛光套装


- + 2 pieces of Tungsten carbide Burs HP 699, F040SFE
- + 1 set of Diamond Disc Em19D20 Mounted
- + 4 pieces of HP Diamond Burs: MK37, FM27, FD30, MD14
- + 4 pieces of Diamond Stones: DG017, DG002, DG004, DG005
- + 4 pieces of Diamond Polishers: D205M, D205F, D205EF, D201EF

3Z00.3**Grinding and Polishing Kit for Zirconia**
氧化锆打磨抛光套装

- + 1 piece HP Tungsten Carbide Burs HP 699
- + 5 pieces Diamond Burs: MG50, MD14, FD30, FM27, MK37;
- + 1 set of Diamond Disc Hm19D20 Mounted
- + 1 piece of Diamond Polisher D1M
- + 1 set of Silicon polisher C104C Mounted

3L04.2**Grinding and Polishing Kit for Lithium Disilicate**
玻璃陶瓷打磨抛光套装

- + 1 set of Diamond Disc HM19D20, Mounted,
- + 3 pieces of Diamond Stones DG017, DG004, DY001EF
- + 3 pieces of HP Diamond Burs MD10, MB12, MR23
- + 3 pieces of Diamond Polisher D1M, D13F, D13UF
- + 1 piece of Silicon Polisher C205C

3L04.3**Grinding and Polishing Kit for Lithium Disilicate**
玻璃陶瓷打磨抛光套装

- + 1 set of Diamond Disc HM19D20, Mounted,
- + 2 pieces of Diamond Stones DG004, DY001EF,
- + 3 pieces of HP Diamond Burs MR23, ML16, MM27,
- + 1 piece of Diamond Polisher D1M
- + 1 piece of Silicon Polisher C205C

3MC02.3**Adjusting and Polishing Kit for Ceramics**
烤瓷工件打磨抛光套装

- + 1 piece of HP Diamond Burs MR23
- + 1 piece of HP Silicon Polisher C203C
- + 1 piece of Diamond Polisher D26f
- + 2 pieces of Polishing Brushes 01G22, N22

3MR09.1**Adjust & Polish Invisible Braces**
隐形牙套调整抛光套装

- + 1 piece of Tungsten Carbide Burs AL02354
- + 1 set of Diamond Disc C19D20, Mounted
- + 2 pieces of Polishing Brushes 01Ms25, 01S22

3A00.1

Popular HP Kit for grinding or polish Metal or Alloy
常用打磨抛光金属、合金的套装



* Universally grind, remove, polish and finish of Alloy, metal, Jewellery, stones etc.
It widely applied to Laboratory, DIY or Hobbies.
+ 5 pieces of HP sintered diamond burs ,
+ 4 pieces of HP Tungsten Carbide Cutters
+ 10 pieces of HP Tungsten Carbide burs

+ 4 pieces of White stones
+ 2 pieces of Pink stones
+ 10 pieces of HP Silicon Polishers
+ 3 pieces of Rubber wheel mounted
+ 1 pieces of Woolen Felt wheel mounted ,
+ 5 pieces of HP Polishing brushes

3MC00.1

Popular HP Kit for grinding or polish Ceramics / Porcelain
常用打磨抛光陶瓷的套装



* Universally grind, remove, polish and finish of Ceramics or Porcelains etc.
It widely applied to Laboratory, DIY or Hobbies.
+ 10 pieces of HP Diamond burs ,
+ 4 pieces of Diamond Discs mounted
+ 6 pieces of HP Green Stones

+ 3 pieces of Silicon Rubber wheel mounted
+ 8 pieces of Silicon Rubber Polishers
+ 4 pieces of Polishing brushes

3LR00.1

Popular HP Kit for grinding or polish Acrylic, Resin
常用打磨抛光树脂胶托的套装



* Universally grind, remove, polish and finish of Acrylic, Resin etc.
It widely applied to Laboratory, DIY or Hobbies.
+ 2 pieces of HP Tungsten Carbide Cutters
+ 32 pieces of HP Silicon rubber Polishers

+ 10 pieces of Polishing brushes
+ 3 pieces of Woolen Felt polishing brushes

3Z00.1

Popular HP Kit for grinding or polish Zirconia, Lithium Disilicate
常用打磨抛光氧化锆、玻璃陶瓷的套装



* Universally grind, polish and finish of Zirconia and Lithium Disilicate etc.
It widely applied to Laboratory, DIY or Hobbies.
+ 4 pieces of HP Diamond Burs
+ 2 pieces of Diamond Discs Mounted

+ 6 pieces of Diamond Stones
+ 4 pieces of Diamond polishing wheels
+ 2 pieces of Diamond polishing burs
+ 5 pieces of Diamond Twist polishers
+ 1 pieces of Polishing brush .

ZS100

Hot HP Mounted stones mixed
车石套装, 常用。




* 100 Pcs, Hot Mounted stones in Green, Pink, White , applied to Alloy, metal, ceramics, Jewellery, stones, hard wood etc.
+ 30 pieces of HP Green Stones
+ 30 pieces of HP Pink Stones

+ 30 pieces of White Stones
+ 10 pieces of Brown Stones

ZD 101

Popular HP diamond burs KIT
常用金刚石磨头套装




MD23, MK47, MJ60, MA23, MC18, MM27, MB21, ML40, MG16, MR23

10 Hot Spec. of HP Diamond, universally grind, remove, polish and finish of Alloy, metal, Jewellery, stones etc. It widely applied to Laboratory, DIY or Hobbies.

用途: 选取10支最常用的金刚石磨头, 形状全, 大小适中, 可以用来处理陶瓷工件。

ZD 101F

Popular HP diamond burs KIT
常用金刚石磨头套装




FD23, FK47, FJ60, FA23, FC18, FM27, FB21, FL40, FG16, FR23

10 Hot Spec. of HP Diamond of Fine grit, universally grind, remove, polish and finish of Zirconia, Lithium Disilicate, Porcelain etc. It widely applied to Laboratory, DIY or Hobbies.

用途: 选取10支最常用的金刚石磨头, 形状全, 大小适中, 可以用来处理氧化锆、玻璃陶瓷, 烤瓷等工件。

ZD 102

Small & Tiny Cut, HP diamond burs KIT
金刚石磨头, 小号套装



MD10, MD14, MD23, MK16, MA16, MA23, MC25, MM12, MM18, MM27, MB12, MB16, MB21, ML16, ML23, MN18, MG16, MR23, MQ23, MU50

20 Spec. small or sharp burs. To create tiny hole, slim curves or tiny, deep grind smartly. Used to Alloy, metal, ceramics, Jewellery, stones, hard wood etc. It widely applied to Laboratory, DIY or Hobbies.

用途: 20支金刚石磨头, 用来开小洞, 细边, 细槽, 或者小的磨削。很轻巧, 灵活。适用于各种金属, 合金, 陶瓷, 珠宝, 宝石, 木头等。

ZD 103

Fast Trimming, HP diamond burs KIT
金刚石磨头套装, 快速磨削




CQ23, CL40, SD35, SA50, SC40, SP70, SG50, SV70, SV60, SF60

20 Spec. HP Diamond for fast trimming. Grinding of Plastic, Jewellery, stones, Hard wood, etc. Save time when removing and work efficiently

用途: 10支常用的金刚石磨头, 头部相对比较大, 可以快速磨削塑料, 木头, 陶瓷和石头等等, 适用于牙科, 珠宝加工和DIY爱好者。

ZD 104

Hot Fully-Sintered diamond KITs
全砂烧结磨头套装, 常用




C 050S, F050S, L050S, B050C, A037, G050, C037F, F050F, M050F, G050F

10 spec. Hot Fully-sintered Diamond HP burs, universal grind, remove, polish and finish of Alloy, metal, Jewellery, stones etc. It widely applied to Lab, DIY or Hobbies.

用途: 10个常用全砂烧结磨头, 5种不同粗细。广泛适用于研磨, 抛光金属, 合金, 珠宝, 宝石。经久耐用, 方便实惠。

ZD 105

Fast Trimming Fully-Sintered diamond KIT
全砂烧结磨头套装, 大磨削




K080S, K100S, J100S, C070S, M070S, B070S, L070S, G060S, F070S, R070S

If workpieces are extremely hard and waste too much time when grinding, this fast trimming kit may help to work efficiently. Applied to Alloy, metal, Jewellery, stones, hard wood, etc.

用途: 10个大头的全砂烧结磨头, 用于粗磨, 去块, 效率高, 经济实惠。

ZT 101

Popular HP Tungsten carbide burs KIT
常用低速钨钢针



HP 1/4, HP 1/2, HP 1, HP 2, HP 3, HP 4, HP 5, HP 8, HP 34, HP 36, HP 38, HP 41, HP 556, HP 557, HP 560, HP 563, HP 699, HP 701, HP 702, HP 704

20 Spec. HP Tungsten Carbide burs are used to create tiny hole, slim curves or tiny, deep grind very smartly. It is applied to Alloy, metal, ceramics, Jewellery, stones, hard wood etc.

用途: 20个型号钨钢车针, 用于雕牙花, 钻小孔, 洞, 开细槽, 细线等, 金属, 合金, 宝石, 坚木都可以用。

ZT 102

Hot Tungsten carbide cutters KIT
钨钢磨头套装, 常用



D050SP, L040SC, L060SC, F060SC, B023SC, L023SC, L040FC, C023FC, L040EFC, F060CC

10 spec. Tungsten carbide Cutter, universal grind, remove, polish and finish of Acrylic, Resin, stones, soft metal or wood.

用途: 10个常用钨钢磨头, 形状全, 大小适中。主要用来研磨, 抛光塑料, 树脂, 软金属, 石头等, 方便实惠。

ZF 101

Acrylic Polishing Kit
树脂基托抛光套装



F060SC, L040SC, U12 x 2, U13 x 2, U16 x 2

Simple and straight way to adjust Acrylic without the problem associated with pumice.

用途: 2个型号的钨钢磨头快速磨削, 三个不同粒度的硅胶磨头依次进行抛光。

ZA 101

Non-precise metal, Alloy Adjust and polish
非贵金属, 铬钴调整与抛光



W-01, W-09, W-08, W-16, W-03, W-06
A201C, A203C, A206C, A201F, A203F

To adjust and polish Cobalt-chromium, nickel chromium Alloy. White stones for adjusting or grinding, Rubber burs for polishing.

用途: 试用于铬钴、镍铬合金等非贵金属的研磨和抛光 6支白玉石用来研磨和调整, 5支硅胶磨头用来抛光。

ZP 100

Sandpaper Polishing burs Kit
砂纸棒套装




240 #, 400 #, 600 #, 800 #, 1000 #, 1500 #, 2000 #, 2500 #, 3000 #, 5000 #

10 different Grit Sandpaper burs are used for polishing Alloy, metal, Jewellery stones, hard woods etc. Working step by step, mirror surface in coming soon, water cooling or dry are acceptable.

用途: 10个不同粒度的砂纸棒, 从240目到5000目不等, 请根据需要, 选取合适的粒度, 耐心细致地使用。

Z 2001

Popular Felt Polishing bur Kit
毛毡抛光磨头套装



12 pieces of High quality Felt polishing burs are widely used to mould, Alloy, metal, Jewellery, stones, Hard woods etc.

12 pieces of High quality Felt polishing burs are widely used to mould, Alloy, metal, Jewellery, stones, Hard woods etc.

用途: 12支常用的毛毡磨头, 广泛适用于各种模具, 金属, 陶瓷, 珠宝, 宝石和木头等。

Samples Collection

样品册

- Flip-Cart Plastic Pages ;
- Smart, Easily fill all burs ;
- Economic and Recyclable ;
- Nice-looking for showing to Customers ;
- Convenient to take to anywhere .



No. 171B5

Empty Collection Wallet
样品本(不带车针)



4 Colors Available ,
100 holes per page,
5 pages per set,
Smart to HP Burs



No. 1704D

Sintered Diamond Burs Collection
全砂烧结磨头样品本



115 specifications ,
Fully Sintered Diamond,
HP shanks



No. 1705D

HP Diamond Burs, Mounted Stones
金刚石磨头, 车石样品本



More than 100 specifications ,
Electroplating Diamond burs
HP Shanks



No. 1707T

Tungsten carbide Burs Collection
高品质钨钢磨头样品本



More than 100 specifications ,
High Quality TC burs,
HP Shanks



No. 1708P

Polishing Burs, Mandrels Collection
抛光磨头, 夹针样品本



More than 100 specifications ,
Silicon, Rubber, Diamond burs
HP Shanks



隐形义齿胶

FlexSoft Partial Dentures



The alternative to conventional Partial Dentures ;
Flexible . Soft. Easy-to-finish.Easy-to-polish
弹性好，易成型，易抛亮，美观大方。



编号 No.	规格 Spec.
D2	Crystal ; 18 gs
C1	Ligth Pink ; 9 gs
C2	Ligth Pink ; 18 gs
B1	Standard Pink ; 9 gs
B2	Standard Pink ; 18 gs
A1	Deep Pink ; 9 gs
A2	Deep Pink ; 18 gs
For bulk pak of 1 kilo.	
C3	Ligth Pink ; 1000 gs
B3	Standard Pink ; 1000 gs
A3	Deep Pink ; 1000 gs



Instructions for Use

Suitable for: Partial and complete dentures, occlusal splints and nightguards.

Not Suitable for: Small unilateral restorations.

Warnings: Use eye protection and heat-resistant gloves when injecting.

Use proper ventilation, vacuum system and masks when grinding prosthodontic resins.

Precautions: Cartridge must be injected within 1 hour of opening the foil pouch.

I. PREPARING THE MASTER MODEL FOR DUPLICATION

All undercuts that will affect the path of insertion should be blocked-out with wax prior to duplication. Do not relieve necklines of teeth or areas desirable for retention. For free-end saddles, it is not necessary to remove the distal undercut of the last natural tooth, unless the area is too deep. For free-end saddle cases where anterior teeth are also being replaced, undercuts in the anterior should also be blocked out.

II. SET-UP

Duplicate blocked-out master model. Articulate the duplicate model and the master model (you will need the master model articulated at finishing steps to check occlusion). Set up the case on the duplicate model using only wax (no base plate).

Allow 1 mm of thickness under the teeth to flow to the labial side and into the retention diatorics in the teeth. Use a small round carbide bur to create mechanical retention in the teeth, from mesial to distal and connected to the center diatoric hole.

III. WAXING

Important: Wax the case with anatomical features to reduce grinding during the finishing stage.

Recommended thickness of wax-up:

1. Palatal thickness 1.2 - 1.5 mm
2. Buccal & Labial Flange 1.2 - 2.0 mm
3. Clasps 1.0 - 1.5 mm
4. Lingual Flange 1.2 - 2.0 mm
5. Lingual Major Connector (Mandibular) 1.7 - 2.0 mm

IV. INVESTING AND SPRUING

Invest the case in the center of the flask, aiming the posterior portion toward the sprue opening. After the stone (or plaster) has set, use two sprues from the flask opening to the most posterior portion of the wax on both sides (see diagram). Sprues should be approximately 7 mm. width. Do not create sharp angles with the sprues, as the heated resin must flow in a straight or slightly curved path. Apply a thin layer of petroleum jelly on the entire stone surface. Then invest top half of flask.

V. REMOVING WAX

After the stone has set, place flask in boiling water for 8 minutes. Open the flask and wash out the wax, making sure all wax residues are eliminated. Allow the flask to cool. Apply thin layers of separator on the entire surface of the stone and allow it to dry completely. Do not apply separator on teeth. When dry, close flask.

VI. INJECTING

Wear safety goggles and heat-protective gloves for this procedure.

Spray inside of the cartridge sleeve with mold release spray and place it into furnace. Turn on the furnace and allow it to reach approximately 550°F (290°C), preheating cartridge sleeve for approximately 15 minutes.

Remove cartridge from foil-lined bag (do not remove cartridge from bag until this step). Spray the cartridge and bronze disk with mold release spray. Insert the cartridge into the cartridge sleeve with the cap end first. Insert the bronze disc into the cartridge sleeve against the end of the cartridge. Use the thick bronze disc for medium tubes and the thin bronze disc for the large tubes. Heat cartridge for 12 minutes.

Position the flask with the sprue end up, directly under the shaft of the press. After 12 minutes, remove the cartridge sleeve from the furnace. Keep the cartridge sleeve in a horizontal position until resting on top of the flask, and then turn the cartridge sleeve into a vertical position. Immediately use the press to inject the heated cartridge into the flask using a steady fast motion. Continue applying pressure until tube bursts and continue with maximum downward pressure for 3 minutes.

Remove the cartridge sleeve from the flask with a slight twisting motion. After removing the cartridge sleeve, let the flask cool for at least 30 minutes. Do not put flask into a cold water bath. Use a plaster nipper to de-flask. Do not use a hammer. Remove the compressed cartridge from the cartridge sleeve with a knock-out punch. Important: Cartridge sleeve must be pre-heated for 10 minutes BEFORE the next injection.

VII. FINISHING AND POLISHING

Remove appliance from the model. Cut off the sprues with a cut-off wheel. Trim the edges of appliance with a coarse grinding wheel. Avoid overheating. Use rubber wheels to smooth scratches. Place the appliance on the articulated master model to check the fit, retention, thickness and occlusion. Clean the tissue surface with a small bristle brush.

To polish, use a coarse dental abrasive, such as pumice, with a wet rag wheel. Apply final shine with a hi-shine compound on low speed with a final buff wheel. Be careful not to overheat appliance.




铝桶

Empty Aluminum Cartridge

编号 No.	规格 D x H, mm
25085	25 x 85 mm
25100	25 x 100 mm
25CAP	24.2 x 16 mm, Cap
22065	22 x 65 mm
22085	22 x 85 mm
22100	22 x 100 mm
22120	22 x 120 mm
22CAP	21.2 x 16 mm, Cap



BurDental 博登科技



博登科技有限公司

B&D Technologies Co.,Ltd

地址: 广东省广州市番禺区洛浦街恒达产业园A-522

Address : A-522, Heng Da Industrial Park, Luo Pu Street, Pan Yu District, 510000,
GuangZhou, GuangDong Province, China

网站(Web): www.burdental.com ; 邮箱(Email): info@burdental.com

手机(Mobile, WhatsApp) : 0086 1866 6599 173