



博登科技有限公司
B&D Technologies Co.,Ltd

Dental Burs 牙科车针、磨头

Vol. 8 : 2024

About Us

B&D Technologies Company Limited are specialized in Rotary things of dental, Jewelry, Beauty or Hobbies & DIY (Do it Yourself) .

Our products include Diamond Burs and discs, Tungsten Carbide Burs or cutters, Flexible Rubber Burs, Wheels or Points, Polishing Brushes, Mounted Stones, CAD/CAM Milling Burs or any Diamond Stones or Polishers.

Our Rotatory Instruments are widely applied to Cut, Grind, Finish, Polish any materials of Zirconia, Lithium Disilicates, Porcelain (Ceramics) , Alloy or Metal, Acrylic (Resin) , PEEK, Composites etc . They can match any Handpiece of FG (Friction Grip) , RA/CA(Right Angle), HP conveniently and precisely .

We have obtained the certificates of ISO 13485 and CE since the year of 2015 . We provide innovative and professional products and service.

Whoever you are dentist or technician, we always make you relaxing, safe and efficient. We are your best friend in dentistry, Jewellery or Hobbies.



ISO 13485



CE 0197



Established in 2007



Clinic Part
临床类

Page 3 - 28 第3至28页



Diamond burs 金刚石车针

Page 4 - 12 第4至12页

Diamond Ortho-strips
金刚石锯弓

Page 13 第13页

Diamond Finishing Strips
金刚石抛光条

Page 14 第14页

Tungsten Carbide burs
钨钢车针

Page 15 - 18 第15至18页

White Stones 白车石
Green Stones 绿车石

Page 19 第19页

Diamond Stones, Polishers
陶瓷（橡胶）金刚车石磨头

Page 20 第20页

**Finishing or Polishing Burs
after Debonding**
（脱胶后）洁牙磨头

Page 21 第21页

Twin-Color Silicon Polishers
硅胶抛光磨头

Page 22 第22页

Silicon Polishers
硅胶抛光磨头

Page 23 第23页

Polishing Brushes 抛光刷

Page 24 第24页

Mandrels, Sterile Blocks
夹针、消毒盒

Page 25 第25页

Kit 套装

Page 26 - 31 第26至31页

Samples Collection
样品册


Page 32 第32页

Dental Diamond Bur

牙科金刚砂车针

Advantages

1. Strictly comply with ISO Standard ;
2. Sharper and More Efficient ;
3. Durable and Valuable ;
4. Full Range of shapes and grits.
5. Fastly develop new items.
6. Individual Design and OEM Service



Grit

- Black / Super Coarse
- Green / Coarse
- Blue / Medium
- Red / Fine
- Yellow / Extra Fine

0.8 19.0 001(801) BR-49	1.0 19.0 001(801) BR-45	1.2 19.0 001(801) BR-46	1.4 19.0 001(801) BR-41	1.6 19.0 001(801) BR-40	1.8 19.0 001(801) BR-31	2.1 19.0 001(801) BR-28	2.3 19.0 001(801) BR-26	2.5 19.0 001(801) BR-25	2.9 19.0 001(801) BR-23
001-008M	001-010M	001-012M	001-014M	001-016M	001-018M	001-021M	001-023M	001-025M	001-029M
圆球型 Ball, Round /BR									

6mm 1.2 22.0 697(801L) BR-L46	1.4 22.0 697(801L) BR-L41	1.6 22.0 697(801L) BR-L40	1.8 22.0 697(801L) BR-L31	1.2/2.3 19.0 002(802) BC-43	1.4/2.3 19.0 002(802) BC-42	1.6/3.5 19.0 002(802) BC-32	1.8/2.3 19.0 002(802) BC-31	1.2/10.0 21.0 494(802L) ENDO-012	1.4/10.0 21.0 494(802L) ENDO-014
697-012M	697-014M	697-016M	697-018M	002-012M	002-014M	002-016M	002-018M	494-012M	494-014M
圆球型 Round with Long Neck /BR					紧口圆球型 Round with Collor /BC				

1.6/10.0 21.0 494(802L) ENDO-016	1.8/10.0 21.0 494(802L) ENDO-018	1.0(0.8)/1.0 19.0 010(805) SI-45	1.2(0.9)/1.2 19.0 010(805) SI-46	1.4(1.0)/1.6 19.0 010(805) SI-47	1.6(1.1)/1.8 19.0 010(805) SI-48	1.4(1.1)/4.0 19.0 225(807) SI-50	1.6(1.2)/4.0 19.0 225(807) ---	1.8(1.3)/5.0 19.0 225(807) SI-51
494-016M	494-018M	010-010M	010-012M	010-014M	010-016M	225-014M	225-016M	225-018M
倒锥型 Inverted Cone / SI								

1.0/2.3 19.0 019(806) DI-41	1.4/3.2 19.0 019(806) DI-42	3.2/5.0 19.0 038(811) EX-11	3.5/7.0 20.0 039(811L) EX-12	4.5/0.7 19.0 044(818) WR-44	4.2/0.9 19.0 067(909) WR-11	4.2/1.8 19.0 068(909) WR-13
019-010M	019-014M	038-032M	039-035M	044-045M	067-042M	068-042M
紧口倒锥型 Inverted Cone with collar/DI		双圆锥型 Double cone/EX		轮型 Wheel/WR		

1.0/4.0 19.0 109(835) SF-41	1.2/4.0 19.0 109(835) SF-31	1.4/4.0 19.0 109(835) SF-42	1.6/4.0 19.0 109(835) ---	1.2/6.0 21.0 110(836) SF-20	1.4/6.0 21.0 110(836) SF-21	1.6/6.0 21.0 110(836) SF-15	1.2/8.0 22.0 111(837) SF-11	1.4/8.0 22.0 111(837) SF-12	1.6/8.0 22.0 111(837) SF-13
109-010M	109-012M	109-014M	109-016M	110-012M	110-014M	110-016M	111-012M	111-014M	111-016M
平头圆柱型 Flat End Cylinder/SF									

1.2/8.0 22.0 130(885) CP-11	1.4/8.0 22.0 130(885) CP-12	1.6/8.0 22.0 130(885) CP-13	1.2/10.0 24.0 131(886) ---	1.4/10.0 24.0 131(886) ---	1.6/10.0 24.0 131(886) ---	1.4(0.7)/10.0 24.0 173(848) TF-11
130-012M	130-014M	130-016M	131-012M	131-014M	131-016M	173-014M
尖头圆柱型 Ogival end Cylinder/CP						平头锥型
						头部直径 头部最小直径 头部长度 总长 ISO编码/欧洲编码 产品代码 规格 型号

金刚石粒度

- 黑色/极粗粒度
- 绿色/粗粒度
- 蓝色/标准粒度
- 红色/较细粒度
- 黄色/极细粒度

	1.0/4.0 19.0 139(838) SR-41		1.2/4.0 19.0 139(838) SR-31		1.4/4.0 19.0 139(838) SR-42		1.4/6.0 21.0 140(880) SR-21		1.2/8.0 22.0 141(881) SR-11		1.4/8.0 22.0 141(881) SR-12		1.6/8.0 22.0 141(881) SR-13		1.2/0.2 21.0 150(839) EX-18
139-010M		139-012M		139-014M		140-014M		141-012M		141-014M		141-016M		150-012M	
圆头柱型 Round End, Cylinder/SR														肩台车针 End Cut	

	1.2/6.0 21.0 164(852) TC-26		1.4/8.0 22.0 165(858) TC-21		1.0/10.0 22.0 166(859) TC-27		1.2/10.0 22.0 166(859) TC-28		1.4/10.0 22.0 166(859) TC-10		1.6/10.0 22.0 166(859) TC-11		1.8/10.0 22.0 166(859) TC-16
164-012M		165-014M		166-010M		166-012M		166-014M		166-016M		166-018M	
尖头锥型 Pointed cone/TC													

	1.0(0.6)/3.0 19.0 169(845) TF-41		1.2(0.7)/4.0 19.0 170(845) TF-42		1.4(0.8)/4.0 19.0 170(845) TF-43		1.6(0.9)/4.0 19.0 170(845) TF-31		1.8(1.0)/5.0 20.0 170(845) TF-23		2.1(1.5)/5.0 20.0 170(845) TF-22		1.4(0.7)/6.5 21.0 171(846) TF-20		1.6(0.8)/6.5 21.0 171(846) TF-21		1.0(0.6)/8.0 22.0 172(847) TF-07		1.2(0.7)/8.0 22.0 172(847) TF-08
169-010M		170-012M		170-014M		170-016M		170-018M		170-021M		171-014M		171-016M		172-010M		172-012M	
平头锥型 Flat End Cone/TF																			

	1.4(0.7)/8.0 22.0 172(847) TF-09		1.6(0.8)/8.0 22.0 172(847) TF-10		1.8(0.9)/8.0 22.0 172(847) TF-56		2.3(1.3)/9.0 22.0 172(847) TF-14		1.4(0.7)/10.0 24.0 173(848) TF-11		1.6(0.8)/10.0 22.0 173(848) TF-12		1.8(0.9)/10.0 22.0 173(848) TF-13
172-014M		172-016M		172-018M		172-023M		173-014M		173-016M		173-018M	

	1.0(0.6)/4.0 19.0 196(849) CR-41		1.2(0.7)/4.0 19.0 196(849) CR-31		1.4(0.8)/4.0 19.0 196(849) CR-42		1.6(0.9)/4.0 19.0 196(849) CR-12		2.0(1.0)/5.0 19.0 196(849) CR-11		1.2(0.7)/6.5 21.0 197(855) ---		1.4(0.7)/6.5 21.0 197(855) TR-20		1.6(0.8)/6.5 21.0 197(855) TR-21		1.8(1.0)/6.5 21.0 197(855) TR-62		2.0(1.1)/6.5 21.0 197(855) TR-63
196-010M		196-012M		196-014M		196-016M		196-020M		197-012M		197-014M		197-016M		197-018M		197-020M	
圆头锥型 Round End Cone/CR/TR																			

	1.2(0.7)/8.0 22.0 198(856) ---		1.4(0.8)/8.0 22.0 198(856) TR-25		1.6(0.9)/8.0 22.0 198(856) TR-24		1.8(1.0)/8.0 22.0 198(856) TR-23		1.2(0.6)/10.0 22.0 199(850) ---		1.4(0.7)/10.0 22.0 199(850) TR-11		1.6(0.8)/10.0 22.0 199(850) TR-12		1.8(0.9)/9.0 22.0 199(850) TR-13		2.2(1.3)/9.0 22.0 199(850) TR-14		2.2(1.3)/11.0 24.0 200(850L) TR-15
198-012M		198-014M		198-016M		198-018M		199-012M		199-014M		199-016M		199-018M		199-022M		200-022M	

	2.4(1.4)/11.0 24.0 200(850L) TR-19		1.0(0.8)/2.8 19.0 233(808) EX-41		1.4(1.1)/4.0 19.0 234(808L) EX-31		1.8(1.3)/5.0 19.0 237(808R) EX-20		2.1(1.6)/5.0 19.0 237(808R) EX-21		3.2(1.8)/5.0 19.0 237(808R) EX-26
200-024M		233-010M		234-014M		237-018M		237-021M		237-032M	
圆头倒截锥型 Round End Inverted cone/EX											

 1.2/6.0 21.0 248(861) FO-41	 1.4/6.0 21.0 248(861) FO-20	 1.4/8.0 22.0 249(862) FO-21	 1.6/8.0 22.0 249(862) FO-22	 1.2/10.0 24.0 250(863) FO-11	 1.4/10.0 24.0 250(863) FO-12	 1.6/10.0 24.0 250(863) FO-13	 1.6/3.2 19.0 257(368,830) FO-30	 1.8/3.6 19.0 257(368,830) FO-32
248-012M	248-014M	249-014M	249-016M	250-012M	250-014M	250-016M	257-016M	257-018M
尖头圆柱型 Ogival end Cylinder/FO							花蕾芽型 Bud/FO	

 2.1/5.0 19.0 257(368,830) ---	 2.3/5.0 19.0 257(368,830) FO-23	 2.5/5.0 19.0 257(368,830) FO-24	 2.8/5.0 19.0 257(368,830) FO-25	 3.0/5.0 19.0 257(368,830) FO-26	 3.2/5.0 19.0 257(368,830) FO-27	 1.8/3.5 19.0 277(379,833) ---	 2.0/3.9 19.0 277(379,833) TR-22	 2.3/4.5 19.0 277(379,833) ---	 3.0/5.0 19.0 277(379,833) TR-28
257-021M	257-023M	257-025M	257-028M	257-030M	257-032M	277-018M	277-020M	277-023M	277-030M
蛋型 Egg Football/TR									

 1.2/6.0 21.0 288(867,877) SO-20	 1.2/8.0 22.0 289(868,878) ---	 1.4/8.0 22.0 289(868,878) SO-21	 1.6/8.0 22.0 289(868,878) ---	 1.2/10.0 24.0 290(869,879) ---	 1.4/10.0 24.0 290(869,879) ---	 1.6/10.0 24.0 290(869,879) ---
288-012M	289-012M	289-014M	289-016M	290-012M	290-014M	290-016M
鱼雷圆柱型 Torpedo/SO						

 1.6(1.0)/8.0 22.0 298(878K) ---	 1.8(1.1)/8.0 22.0 298(878K) ---	 1.6(1.0)/10.0 24.0 299(879K) ---	 1.8(1.1)/10.0 24.0 299(879K) ---	 1.4/8.0 22.0 115(842SP) BC-11	 1.6/8.0 22.0 178(848SP) BC-12	 1.6/5.0 19.0 465(820) ---
298-016M	298-018M	299-016M	299-018M	115-014M	178-016M	465-016M
鱼雷圆锥型 Torpedo,Conical				邻面去釉 Interdental		

 cutting depth 0.3mm 1.6(1.0)/6.0 21.0 552(834) VP-20	 cutting depth 0.4mm 1.8(1.0)/6.0 21.0 552(834) VP-21	 cutting depth 0.5mm 2.1(1.1)/6.0 21.0 552(834) VP-22	 1.8(1.1)/4.5 19.0 544(845KR) RS-31 or EX-38	 1.8(1.1)/7.5 21.0 545(846KR) RS-21 or EX-29	 1.6(0.8)/8.0 22.0 546(847KR) 847KR016	 1.8(0.9)/10.0 24.0 553(848KR) 848KR018	 1.9(1.2)/10.5 23.0 553(848KR) RS-11 or EX-17
552-016M	552-018M	552-021M	544-018M	545-018M	546-016M	553-018M	553-019M
定深车针 Depth Marking				平头圆柱型 Taper, Round Edge			

 1.6.5mm FG S 313	 1.0 16.0 001(801) BR-S45	 1.2 16.0 001(801) BR-S46	 1.4 16.0 001(801) BR-S41	 1.2(0.9)/1.2 16.0 010(805) SI-S46	 1.6(1.1)/1.8 16.0 010(805) SI-S48	 1.4/8.0 18.5 165(858) TC-S21	 1.0(0.6)/3.0 16.5 169(845) TF-S41	 1.6(0.9)/4.0 14.5 170(845) TF-SS31
FG Short	001-S010M	001-S012M	001-S014M	010-S012M	010-S016M	165-S014M	169-S010M	170-S016M
短柄 Short								

 1.6(0.9)/4.0 16.5 170(845) TF-S31	 1.8(1.0)/5.0 17.0 170(845) TF-S23	 2.1(1.5)/5.0 17.0 170(845) TF-S22	 1.6(0.8)/6.5 18.5 171(846) TF-S21	 1.4(0.7)/6.5 18.5 171(846) TF-S20	 1.4(0.7)/6.5 18.5 197(855) TR-S20	 1.6(0.8)/6.5 18.5 197(855) TR-S21	 1.8(0.9)/9.0 21.0 199(850) TR-S13	 1.4/8.0 19.5 249(862) FO-S21
170-S016M	170-S018M	170-S021M	171-S016M	171-S014M	197-S014M	197-S016M	199-S018M	249-S014M

								
Super Coarse Black	1.4 19.0 001(801) BR-41S	1.6 19.0 001(801) BR-40S	1.8 19.0 001(801) BR-31S	2.1 19.0 001(801) BR-28S	2.3 19.0 001(801) BR-26S	2.5 19.0 001(801) BR-25S	1.2/10.0 21.0 494(802L) ENDO-012S	1.4/10.0 21.0 494(802L) ENDO-014S
黑色环极粗砂	001-014S	001-016S	001-018S	001-021S	001-023S	001-025S	494-012S	494-014S
圆球型 Ball , Round /BR							根管治疗车针 /Root canal therapy	

									
1.6/10.0 21.0 494(802L) ENDO-016S	1.8/10.0 21.0 494(802L) ENDO-018S	1.6(1.1)/1.8 19.0 010(805) SI-48S	3.2/5.0 19.0 038(811) EX-11S	3.5/7.0 20.0 039(811L) EX-12S	1.8(1.3)/5.0 19.0 237(808R) EX-20S	2.1(1.6)/5.0 19.0 237(808R) EX-21S	1.2/4.0 19.0 109(835) SF-31S	1.4/4.0 19.0 109(835) SF-42S	1.4/6.0 21.0 110(836) SF-21S
494-016S	494-018S	010-016S	038-032S	039-035S	237-018S	237-021S	109-012S	109-014S	110-014S
倒锥型 Inverted Cone /SI			双圆锥型 Double cone/EX		圆头倒锥型 Round End Inverted cone/EX		平头圆柱型 Flat End Cylinder/SF		

									
1.6/6.0 21.0 110(836) SF-15S	1.2/8.0 22.0 111(837) SF-11S	1.4/8.0 22.0 111(837) SF-12S	1.6/8.0 22.0 111(837) SF-13S	1.2/8.0 22.0 141(881) SR-11S	1.4/8.0 22.0 141(881) SR-12S	1.6/8.0 22.0 141(881) SR-13S	1.4/8.0 22.0 165(858) TC-21S	1.4/10.0 22.0 166(859) TC-10S	1.6/10.0 22.0 166(859) TC-11S
110-016S	111-012S	111-014S	111-016S	141-012S	141-014S	141-016S	165-014S	166-014S	166-016S
圆头柱型 Round End, Cylinder/SR							尖头锥型 Pointed cone/TC		

						
1.8/10.0 22.0 166(859) TC-16S	1.6(0.8)/6.5 21.0 171(846) TF-21S	1.4(0.7)/8.0 22.0 172(847) TF-09S	1.6(0.8)/8.0 22.0 172(847) TF-10S	1.4(0.7)/10.0 24.0 173(848) TF-11S	1.6(0.8)/10.0 22.0 173(848) TF-12S	1.8(0.9)/10.0 22.0 173(848) TF-13S
166-018S	171-016S	172-014S	172-016S	173-014S	173-016S	173-018S
平头锥型 Flat End Cone/TF						

					
1.6(0.8)/6.5 21.0 197(855) TR-21S	1.4(0.8)/8.0 22.0 198(856) TR-25S	1.6(0.9)/8.0 22.0 198(856) TR-24S	1.4(0.7)/10.0 22.0 199(850) TR-11S	1.6(0.8)/10.0 22.0 199(850) TR-12S	1.8(0.9)/9.0 22.0 199(850) TR-13S
197-016S	198-014S	198-016S	199-014S	199-016S	199-018S
圆头锥型 Round End Cone/CR/TR					

						
1.2/6.0 21.0 248(861) FO-41S	1.4/6.0 21.0 248(861) FO-20S	1.4/8.0 22.0 249(862) FO-21S	1.6/8.0 22.0 249(862) FO-22S	1.2/10.0 24.0 250(863) FO-11S	1.4/10.0 24.0 250(863) FO-12S	1.6/10.0 24.0 250(863) FO-13S
248-012S	248-014S	249-014S	249-016S	250-012S	250-014S	250-016S
尖头圆柱型 Ogival end Cylinder/FO						

									
1.8/3.6 19.0 257(368,830) FO-32S	2.1/5.0 19.0 257(368,830) ---	2.3/5.0 19.0 257(368,830) FO-34S	2.5/5.0 19.0 257(368,830) FO-24S	2.8/5.0 19.0 257(368,830) FO-25S	3.0/5.0 19.0 257(368,830) FO-26S	3.2/5.0 19.0 257(368,830) FO-27S	2.0/3.9 19.0 277(379,833) TR-22S	3.0/5.0 19.0 277(379,833) TR-28S	1.4/8.0 22.0 289(868,878) SO-21S
257-018S	257-021S	257-023S	257-025S	257-028S	257-030S	257-032S	277-020S	277-030S	289-014S
花蕾芽型 Bud/FO							蛋型 Egg Football/TR		鱼雷圆柱型 Torpedo/SO

								
Coarse Green	001(801) BR-46C	001(801) BR-41C	001(801) BR-40C	001(801) BR-31C	001(801) BR-28C	001(801) BR-26C	001(801) BR-25C	001(801) BR-23C
绿色环粗砂	001-012C	001-014C	001-016C	001-018C	001-021C	001-023C	001-025C	001-029C
圆球型 Ball, Round /BR								

								
494(802L) ENDO-012C	494(802L) ENDO-014C	494(802L) ENDO-016C	494(802L) ENDO-018C	010(805) SI-47C	010(805) SI-48C	038(811) EX-11C	039(811L) EX-12C	068(909) WR-13C
494-012C	494-014C	494-016C	494-018C	010-014C	010-016C	038-032C	039-035C	068-042C
根管治疗车针/Root canal therapy			倒锥型 Inverted Cone / SI			双圆锥型 Double cone/EX		轮型 Wheel/WR

								
109(835) SF-31C	109(835) SF-42C	109(835) ---	110(836) SF-20C	110(836) SF-21C	110(836) SF-15C	111(837) SF-11C	111(837) SF-12C	111(837) SF-13C
109-012C	109-014C	109-016C	110-012C	110-014C	110-016C	111-012C	111-014C	111-016C
平头圆柱型 Flat End Cylinder/SF								

					
130(885) CP-11C	130(885) CP-12C	130(885) CP-13C	131(886) ---	131(886) ---	131(886) ---
130-012C	130-014C	130-016C	131-012C	131-014C	131-016C
Ogival end Cylinder/CP					

					
140(880) SR-21C	141(881) SR-11C	141(881) SR-12C	141(881) SR-13C	142(882) ---	142(882) ---
140-014C	141-012C	141-014C	141-016C	142-014C	142-016C
圆头柱型 Round End, Cylinder/SR					

						
165(858) TC-21C	166(859) TC-10C	166(859) TC-11C	166(859) TC-16C	167(859L) ---	167(859L) ---	167(859L) ---
165-014C	166-014C	166-016C	166-018C	167-014C	167-016C	167-018C
尖头锥型 Pointed cone/TC						

								
170(845) TF-43C	171(846) TF-20C	172(847) TF-09C	172(847) TF-10C	172(847) TF-56C	172(847) TF-14C	173(848) TF-11C	173(848) TF-12C	173(848) TF-13C
170-014C	171-014C	172-014C	172-016C	172-018C	172-023C	173-014C	173-016C	173-018C
平头锥型 Flat End Cone/TF								

	1.4(0.7)/6.5 21.0 197(855) TR-20C		1.6(0.8)/6.5 21.0 197(855) TR-21C		1.8(1.0)/6.5 21.0 197(855) TR-62C		1.4(0.8)/8.0 22.0 198(856) TR-25C		1.6(0.9)/8.0 22.0 198(856) TR-24C		1.8(1.0)/8.0 22.0 198(856) TR-23C		1.4(0.7)/10.0 22.0 199(850) TR-11C		1.6(0.8)/10.0 22.0 199(850) TR-12C		1.8(0.9)/9.0 22.0 199(850) TR-13C		2.2(1.3)/9.0 22.0 199(850) TR-14C
197-014C 197-016C 197-018C 198-014C 198-016C 198-018C 199-014C 199-016C 199-018C 199-022C																			
圓頭錐型 Round End Cone/CR/TR																			


	1.4/8.0 22.0 249(862) FO-21C		1.6/8.0 22.0 249(862) FO-22C		1.2/10.0 24.0 250(863) FO-11C		1.4/10.0 24.0 250(863) FO-12C		1.6/10.0 24.0 250(863) FO-13C		1.8/9.5 21.0 190 FO-54C
249-014C 249-016C 250-012C 250-014C 250-016C 190-018C											
尖頭圓柱型 Ogival end Cylinder/FO											

	1.8/3.6 19.0 257(368,830) FO-32C		2.1/5.0 19.0 257(368,830) ---		2.3/5.0 19.0 257(368,830) FO-34C		2.5/5.0 19.0 257(368,830) FO-24C		2.8/5.0 19.0 257(368,830) FO-25C
257-018C 257-021C 257-023C 257-025C 257-028C									
花蕾芽型 Bud/FO									

	1.2/8.0 22.0 289(868,878) ---		1.4/8.0 22.0 289(868,878) SO-21C		1.6/8.0 22.0 289(868,878) ---		1.2/10.0 24.0 290(869,879) ---		1.4/10.0 24.0 290(869,879) ---		1.6/10.0 24.0 290(869,879) ---
289-012C 289-014C 289-016C 290-012C 290-014C 290-016C											
魚雷圓柱型 Torpedo/SO											

					
1.2 19.0 001(801) BR-46F	1.4 19.0 001(801) BR-41F	1.6 19.0 001(801) BR-40F	1.8 19.0 001(801) BR-31F	2.1 19.0 001(801) BR-28F	2.3 19.0 001(801) BR-26F
001-012F 001-014F 001-016F 001-018F 001-021F 001-023F					
紅色環紋細砂 圓球型 Ball, Round /BR					

								
1.4/6.0 21.0 110(836) SF-21F	1.2/8.0 22.0 111(837) SF-11F	1.4/8.0 22.0 111(837) SF-12F	1.6/8.0 22.0 111(837) SF-13F	1.4/6.0 21.0 140(880) SR-21F	1.2/8.0 22.0 141(881) SR-11F	1.4/8.0 22.0 141(881) SR-12F	1.6/8.0 22.0 141(881) SR-13F	1.2/0.2 21.0 150(839) EX-18F
110-014F 111-012F 111-014F 111-016F 140-014F 141-012F 141-014F 141-016F 150-012F								
平頭圓柱型 Flat End Cylinder/SF 圓頭柱型 Round End, Cylinder/SR End Cut								

							
1.2/6.0 21.0 164(852) TC-26F	1.2/8.0 22.0 165(858) ---	1.4/8.0 22.0 165(858) TC-21F	1.0/10.0 22.0 166(859) TC-27F	1.2/10.0 22.0 166(859) TC-28F	1.4/10.0 22.0 166(859) TC-10F	1.6/10.0 22.0 166(859) TC-11F	1.8/10.0 22.0 166(859) TC-16F
164-012F 165-012F 165-014F 166-010F 166-012F 166-014F 166-016F 166-018F							
尖頭錐型 Pointed cone/TC							

	1.4(0.8)/4.0 19.0 170(845) TF-43F		1.4(0.7)/6.5 21.0 171(846) TF-20F		1.6(0.8)/6.5 21.0 171(846) TF-21F		1.2(0.7)/8.0 22.0 172(847) TF-08F		1.4(0.7)/8.0 22.0 172(847) TF-09F		1.6(0.8)/8.0 22.0 172(847) TF-10F		1.4(0.7)/10.0 24.0 173(848) TF-11F		1.6(0.8)/10.0 22.0 173(848) TF-12F		1.8(0.9)/10.0 22.0 173(848) TF-13F
170-014F		171-014F		171-016F		172-012F		172-014F		172-016F		173-014F		173-016F		173-018F	
平头锥型 Flat End Cone/TF																	

	1.6(0.9)/4.0 19.0 196(849) CR-12F		2.0(1.0)/5.0 19.0 196(849) CR-11F		1.4(0.7)/6.5 21.0 197(855) TR-20F		1.6(0.8)/6.5 21.0 197(855) TR-21F		1.8(1.0)/6.5 21.0 197(855) TR-62F		1.4(0.8)/8.0 22.0 198(856) TR-25F		1.6(0.9)/8.0 22.0 198(856) TR-24F		1.8(1.0)/8.0 22.0 198(856) TR-23F		1.2(0.6)/10.0 22.0 199(850) ---		1.4(0.7)/10.0 22.0 199(850) TR-11F
196-016F		196-020F		197-014F		197-016F		197-018F		198-014F		198-016F		198-018F		199-012F		199-014F	
圆头锥型 Round End Cone/CR/TR																			


	1.6(0.8)/10.0 22.0 199(850) TR-12F		1.8(0.9)/9.0 22.0 199(850) TR-13F		2.2(1.3)/9.0 22.0 199(850) TR-14F		1.4(1.1)/4.0 19.0 234(808L) EX-31F		1.8(1.3)/5.0 19.0 237(808R) EX-20F		2.1(1.6)/5.0 19.0 237(808R) EX-21F
199-016F		199-018F		199-022F		234-014F		237-018F		237-021F	
圆头倒截锥型 Round End Inverted cone/EX											

	1.2/6.0 21.0 248(861) FO-41F		1.4/6.0 21.0 248(861) FO-20F		1.4/8.0 22.0 249(862) FO-21F		1.6/8.0 22.0 249(862) FO-22F		1.2/10.0 24.0 250(863) FO-11F		1.4/10.0 24.0 250(863) FO-12F		1.6/10.0 24.0 250(863) FO-13F
248-012F		248-014F		249-014F		249-016F		250-012F		250-014F		250-016F	
尖头圆柱型 Ogival end Cylinder/FO													

	1.6/3.2 19.0 257(368,830) FO-30F		1.8/3.6 19.0 257(368,830) FO-32F		2.1/5.0 19.0 257(368,830) ---		2.3/5.0 19.0 257(368,830) FO-34F		2.5/5.0 19.0 257(368,830) FO-24F		2.8/5.0 19.0 257(368,830) FO-25F		3.0/5.0 19.0 257(368,830) FO-26F		3.2/5.0 19.0 257(368,830) FO-27F
257-016F		257-018F		257-021F		257-023F		257-025F		257-028F		257-030F		257-032F	
花蕾芽型 Bud/FO															

	1.2/6.0 21.0 288(867,877) SO-20F		1.2/8.0 22.0 289(868,878) ---		1.4/8.0 22.0 289(868,878) SO-21F		1.6/8.0 22.0 289(868,878) ---
288-012F		289-012F		289-014F		289-016F	
魚雷圓柱型 Torpedo/SO							


	0.6 16.8		0.7/2.2 16.8		0.7/3.8 16.8		0.7/2.2 16.8		0.6/2.2 16.8		0.9/1.2 16.8		0.7/2.2 16.8		0.7/3.6 16.8		0.8/3.2 16.8		1.3/6.9 16.8
001 CD-50F		138 CD-51F		139 CD-52F		161 CD-53F		246 CD-54F		254 CD-55F		246 CD-56F		247 CD-57F		108 CD-58F		171 CD-59F	
001-S006F		138-S007F		139-S007F		161-S007F		246-S006F		254-S009F		246-S007F		247-S007F		108-S008F		171-S013F	
儿童车针 FG Short Burs for Kids / CD																			



黄色环 细砂


001-012EF 001-014EF 001-016EF 001-018EF 001-021EF 001-023EF

圆球型 Ball, Round /BR




164-012EF 165-012EF 165-014EF 166-010EF 166-012EF 166-014EF 166-016EF 166-018EF

尖头锥型 Pointed cone/TC




171-016EF 172-016EF 173-014EF 173-016EF 173-018EF

平头锥型 Flat End Cone/TF




196-016EF 196-020EF 197-012EF 197-014EF 197-016EF 197-018EF 198-012EF 198-014EF 198-016EF 198-018EF

圆头锥型 Round End Cone/CR/TR




199-012EF 199-014EF 199-016EF 199-018EF 234-014EF 237-018EF 237-021EF

圆头倒截锥型 Round End Inverted cone/EX



540-010EF 248-012EF 248-014EF 249-014EF 249-016EF 250-012EF 250-014EF 250-016EF


尖头圆柱型 Ogival end Cylinder/FO



257-016EF 257-018EF 257-021EF 257-023EF 257-025EF 257-028EF 257-030EF 257-032EF 288-012EF 289-014EF

花蕾芽型 Bud/FO

鱼雷圆柱型 Torpedo/SO

 25mm FG XL 316									
	1.0 25.0 001(801) BR-X45	1.2 25.0 001(801) BR-X46	1.4 25.0 001(801) BR-X41	1.6 25.0 001(801) BR-X40	1.8 25.0 001(801) BR-X31	1.2(0.9)/1.2 25.0 010(805) SI-X46	1.2/8.0 28.0 111(837) SF-X11	1.4/8.0 28.0 111(837) SF-X12	
FG Surgical 加长柄	001-X010M	001-X012M	001-X014M	001-X016M	001-X018M	010-X012M	111-X012M	111-X014M	
1.6/8.0 28.0 111(837) SF-X13	1.6/10.0 29.0 166(859) TC-X11	1.4(0.7)/10.0 29.0 173(848) TF-X11	1.6(0.8)/10.0 29.0 173(848) TF-X12	1.8(0.9)/10.0 29.0 173(848) TF-X13	1.6(0.8)/10.0 29.0 199(850) TR-X12	2.2(1.3)/11.0 29.0 200(850L) TR-X15	1.4/8.0 28.0 249(862) FO-X21	1.6/8.0 28.0 249(862) FO-X22	
111-X016M	166-X016M	173-X014M	173-X016M	173-X018M	199-X016M	200-X022M	249-X014M	249-X016M	
Super Coarse 1.6 25.0 001(801) BR-X40S	1.8 25.0 001(801) BR-X31S	Coarse 1.4 25.0 001(801) BR-X41C	1.6 25.0 001(801) BR-X40C	1.8 25.0 001(801) BR-X31C					
001-X016S	001-X018S	001-X014C	001-X016C	001-X018C					
 22.0mm RA, CA 204									
	1.2 22.0 001(801) BR-R46	1.4 22.0 001(801) BR-R41	1.6 22.0 001(801) BR-R40	1.8 22.0 001(801) BR-R31	2.3 22.0 001(801) BR-R26	1.0(0.8)/1.0 22.0 010(805) SI-R45	1.2(0.9)/1.2 22.0 010(805) SI-R46	1.4(1.0)/1.6 22.0 010(805) SI-R47	
RA, CA 弯机柄	001-R012M	001-R014M	001-R016M	001-R018M	001-R023M	010-R010M	010-R012M	010-R014M	
1.6(1.1)/1.8 22.0 010(805) SI-R48	1.0/4.0 22.0 109(835) SF-R41	1.2/4.0 22.0 109(835) SF-R31	1.4/4.0 22.0 109(835) SF-R42	1.2/8.0 24.0 111(837) SF-R11	1.4/8.0 24.0 111(837) SF-R12	1.6/8.0 24.0 111(837) SF-R13	1.2/8.0 24.0 141(881) SR-R11	1.4/8.0 24.0 141(881) SR-R12	1.4/8.0 24.0 165(858) TC-R21
010-R016M	109-R010M	109-R012M	109-R014M	111-R012M	111-R014M	111-R016M	141-R012M	141-R014M	165-R014M
1.6/10.0 27.0 166(859) TC-R11	1.0(0.6)/3.0 22.0 169(845) TF-R41	1.2(0.7)/4.0 22.0 170(845) TF-R42	1.8(1.0)/8.0 25.0 198(856) TR-R23	1.6(0.8)/10.0 27.0 199(850) TR-R12	1.8/3.6 22.0 257(368,830) FO-R32	1.6 22.0 001(801) BR-R40F	1.8 22.0 001(801) BR-R31F	2.3 22.0 001(801) BR-R26F	1.4/8.0 24.0 165(858) TC-R21F
166-R016M	169-R010M	170-R012M	198-R018M	199-R016M	257-R018M	001-R016F	001-R018F	001-R023F	165-R014F
1.8/3.6 22.0 257(368,830) FO-R32F	Extra Fine 1.6 22.0 001(801) BR-R40EF	1.8 22.0 001(801) BR-R31EF	2.3 22.0 001(801) BR-R26EF	去除肉芽组织 车针套装 Degranulation Kit					
257-R018F	001-R016EF	001-R018EF	001-R023EF	001-R010	001-R025	001-R030	001-R035		

用于制备烧结氧化锆的金刚石车针

Diamond burs for preparing sintered Zirconia

氧化锆是一种具有良好发展前景的假牙材料。氧化锆的特殊性能对制备提出了极高的要求。本系列金刚石钻头专门针对氧化锆陶瓷的极高硬度而设计。它们具有比标准金刚石钻头更长的使用寿命和更大的切削能力。它们是根管治疗中切割氧化锆冠的理想工具。

Zirconia is a material with a promising future in prosthetic dentistry. The special properties of zirconia place extreme demands on the preparation. This series of our diamond burs are tailored to the extreme hardness of zirconia ceramic. They have a long service life and much greater cutting capacity than standard diamond burs. They are ideal for sectioning zirconia crowns during endodontic treatment.

									
2.1/10.0 24.0	3.2/5.0 19.0	1.2/8.0 22.0	1.4/8.0 22.0	2.3(1.3)/9.0 22.0	1.6(0.8)/10.0 22.0	1.4/10.0 24.0	1.6(0.8)/8.0 22.0	1.8(0.9)/10.0 24.0	1.6(0.9)/8.0 22.0
494 (802L)	038 (811)	141 (881)	141 (881)	172 (847)	173 (848)	290 (869,879)	546 (847KR)	553 (848KR)	198 (856)
Z494-021M	Z038-032M	Z141-012M	Z141-014M	Z172-023M	Z173-016M	Z290-014M	Z546-016M	Z553-018M	Z198-016M
—	ZEX-11	ZSR-11	ZSR-12	ZTF-14	ZTF-12	—	—	—	ZTR-24
									
1.4 22.0	1.2/8.0 22.0	1.4/8.0 22.0	1.4/10.0 24.0	1.6(0.8)/8.0 22.0	1.8(0.9)/10.0 24.0	1.6(0.9)/8.0 22.0	2.5/5.0 19.0	1.6/3.2 19.0	1.0/11.5 24.0
697 (801L)	141 (881)	141 (881)	290 (869,879)	546 (847KR)	553 (848KR)	198 (856)	257 (368,830)	277 (379,833)	167 (859L)
Z697-014F	Z141-012F	Z141-014F	Z290-014F	Z546-016F	Z553-018F	Z198-016F	Z257-025F	Z277-016F	Z167-010F
ZBR-L41F	ZSR-11F	ZSR-12F	—	—	—	ZTR-24F	ZFO-24F	—	—
									
1.4 22.0	1.2/8.0 22.0	1.4/8.0 22.0	1.4/10.0 24.0	1.6(0.8)/8.0 22.0	1.8(0.9)/10.0 24.0	1.6(0.9)/8.0 22.0	2.5/5.0 19.0	1.6/3.2 19.0	1.0/11.5 24.0
697 (801L)	141 (881)	141 (881)	290 (869,879)	546 (847KR)	553 (848KR)	198 (856)	257 (368,830)	277 (379,833)	167 (859L)
Z697-014EF	Z141-012EF	Z141-014EF	Z290-014EF	Z546-016EF	Z553-018EF	Z198-016EF	Z257-025EF	Z277-016EF	Z167-010EF
ZBR-L41EF	ZSR-11EF	ZSR-12EF	—	—	—	ZTR-24EF	ZFO-24EF	—	—

快速切削全陶瓷，烧结氧化锆冠的黑色金刚石车针（超粗）



Fast reducing the All-Ceramics, Zirconia crown with the Black burs (Super coarse)

									
1.2/8.0 22.0	1.4/8.0 22.0	1.6(0.8)/10.0 22.0	1.6/10.0 22.0	1.8/10.0 22.0	1.8(0.9)/9.0 22.0	2.2(1.3)/9.0 22.0	1.6(0.8)/8.0 22.0	1.8(0.9)/10.0 24.0	2.3/5.0 19.0
141 (881)	141 (881)	173 (848)	166 (859)	166 (859)	199 (850)	199 (850)	546 (847KR)	553 (848KR)	257 (368,830)
Z141-012S	Z141-014S	Z173-016S	Z166-016S	Z166-018S	Z199-018S	Z199-022S	Z546-016S	Z553-018S	Z257-023S
ZSR-11S	ZSR-12S	ZTF-12S	ZTC-11S	ZTC-16S	ZTR-13S	ZTR-14S	—	—	ZFO-34S

牙体制备专用车针

1. 高品质专业备牙车针，有前后牙之分。
2. 制备肩台时不会出现羽状肩台，避免了因备牙而引起的修复体无法就位和颈缘崩裂现象。
3. 保证肩台宽度一致。

Especially for All-ceramic restoration. Pin guide not Diamond coated avoids injury of the gingival edge;
508T-014XC wearing-off the contact points mesial-distal;
001-014M determination of prep-bounds
198C-016XC 021XC 508T-016XC 021XC for circular wearing-off without damaging the prep-stage;
508T-016F;021F for precise grinding of the pre-stage with the help of a precision guiding-pin.

									
1.4/8.0 21.5	1.4 19.0	1.6/8.0 20.5	1.6/9.0 21.5	1.6/9.0 21.5	1.4/8.0 21.5	1.6 16.0	2.1/8.0 19.0	2.1/9.0 20.0	2.1/9.0 20.0
508T-014XC	001-014M	198C-016XC	508T-016XC	508T-016F	508T-014XC	001-016M	198C-021XC	508T-021XC	508T-021F



品名	FG 十支黄色塑料盒
Name	FG 10 pcs yellow plastic box
尺寸/Size	44mm*30.5mm
重量/Weight	6.6g



品名	FG 十支红色塑料盒
Name	FG 10 pcs red plastic box
尺寸/Size	44mm*30.5mm
重量/Weight	6.6g



品名	FG 十支蓝色塑料盒
Name	FG 10 pcs blue plastic box
尺寸/Size	44mm*30.5mm
重量/Weight	6.6g



品名	FG 十支绿色塑料盒
Name	FG 10 pcs green plastic box
尺寸/Size	44mm*30.5mm
重量/Weight	6.6g



品名	FG 十支黑色塑料盒
Name	FG 10 pcs black plastic box
尺寸/Size	44mm*30.5mm
重量/Weight	6.6g



品名	FG 五支黄色塑料盒
Name	FG 5 pcs yellow plastic box
尺寸/Size	35mm*29mm
重量/Weight	6.4g



品名	FG 五支红色塑料盒
Name	FG 5 pcs red plastic box
尺寸/Size	35mm*29mm
重量/Weight	6.4g



品名	FG 五支蓝色塑料盒
Name	FG 5 pcs blue plastic box
尺寸/Size	35mm*29mm
重量/Weight	6.4g



品名	FG 五支绿色塑料盒
Name	FG 5 pcs green plastic box
尺寸/Size	35mm*29mm
重量/Weight	6.4g



品名	FG 五支黑色塑料盒
Name	FG 5 pcs black plastic box
尺寸/Size	35mm*29mm
重量/Weight	6.4g



品名	FG 十支推拉盒
Name	FG 10 pcs plastic box
尺寸/Size	35mm*33mm
重量/Weight	6.8g



品名	HP 30支塑料盒
Name	HP 30 pcs plastic box
尺寸/Size	144mm*71mm
重量/Weight	77.2g



品名	HP 十支塑料盒
Name	HP 10 pcs plastic box
尺寸/Size	108mm*69mm
重量/Weight	30.5g



品名	HP 十支吸塑
Name	HP 10pcs blister card
尺寸/Size	61mm*165mm
重量/Weight	22.4g



品名	FG 五支吸塑
Name	FG 5pcs blister card
尺寸/Size	95.5mm*56mm
重量/Weight	4.8g



品名	FG 五支吸塑
Name	FG 5pcs blister card
尺寸/Size	81mm*47mm
重量/Weight	4.2g



品名	HP 五支吸塑
Name	HP 5pcs blister card
尺寸/Size	98.5mm*71mm
重量/Weight	10.5g



品名	FG 十支吸塑
Name	FG 10pcs blister card
尺寸/Size	80mm*69mm
重量/Weight	6.7g



品名	FG 十支吸塑
Name	FG 10pcs blister card
尺寸/Size	105.5mm*80mm
重量/Weight	6.9g



品名	FG 五连包
Name	FG Five packs
尺寸/Size	149.5mm*50mm
重量/Weight	3.2g

产品套装 Product Kit



备冠/前牙套装
Crown Preparation Anterior Teeth Kit



备冠/后牙套装
Crown Preparation Molar Teeth Kit



备洞套装
Cavity Preparation Kit



(高) 嵌体预备套装
Inlay Preparation Kit



烤瓷牙抛光套装
Ceramics/Composite Polishing Kit



儿童微创修复套装
Kiddies Kit



烤瓷制备套装
Ceramic Restoration Kit



备肩台套装
Porcelain Shouldered Abutment Polishing Kit



烤瓷制备/修复套装
Porcelain Preparation / Repair Kit



制备/修复套装
Preparation Repair Kit



护龈综合套装
Gingival Bevel Burs Kit



前牙全冠牙体预备套装
Crown Preparation Anterior Teeth Kit



磨牙全冠牙体预备套装
Crown Preparation Molar Teeth Kit



复合树脂金刚砂车针套装
Composite Resin Kit



微创窝洞制备套装 (粗砂)
Minimally Invasive Cavity Preparation Coarse Kit



微创窝洞制备套装 (细砂)
Minimally Invasive Cavity Preparation Fine Kit



树脂修形套装
Composite Repair Kit



陶瓷复合抛光套装
Porcelain Ceramics Composite Polishing Kit



瓷冠制备套装
Porcelain Crown Preparation Kit



常规粗砂套装
Regular Coarse Sand Kit



磨牙修复套装
Molars Restoration Kit



备肩台套装
Shoulder preparation Kit




备冠/前牙套装
Crown Preparation Anterior Teeth Kit



邻面去釉套装
Interproximal enamel kit



氧化锆金刚砂车针套装
Zirconia Diamond Kit



氧化锆金刚砂车针套装
Diamond Zirconia Finishing Kit



氧化锆金刚砂车针套装
Zirconia Kit



破除氧化锆金刚砂车针套装
Diamond Zirconia Cutting Kit



去除肉芽组织车针套装
Degranulation Kit

Diamond Ortho-strips

金剛石锯弓




Creating additional space by interproximal enamel reduction (IPR) :

When opening the interdental space and reducing, finishing and polishing enamel in orthodontics, care must be taken to avoid grooves and scratches. it's important to proceed step by step, following a coarse-to-fine sequence.



Item No.  1

Thickness (mm)



Separating

OS40

0.40 mm



Reducing

OS30

0.30 mm



Contouring

OS25


0.25 mm



Finishing

OS20

0.20 mm



Pre-polishing

OS15

0.15 mm



Polishing

OS10


0.10 mm



No-Diamond Coated
with serrated edge

OS05

0.05 mm



Diamond finishing strips

金刚石抛光条



Perforated Diamond finishing strips with and without serrated edge ;

Double side Coated / Single side Coated ;

with or without honeycomb design .

	45 μ m	Length x Width x Thickness : 148 x 3.75 x 0.16mm	No. : SWS37		Double side Coated
	45 μ m	Length x Width x Thickness : 148 x 3.75 x 0.13mm	No. : WS37		Single side Coated
	30 μ m	Length x Width x Thickness : 148 x 3.75 x 0.12mm	No. : SWS37F		Double side Coated
	30 μ m	Length x Width x Thickness : 148 x 3.75 x 0.10mm	No. : WS37F		Single side Coated
	15 μ m	Length x Width x Thickness : 148 x 3.75 x 0.10mm	No. : SWS37EF		Double side Coated
	15 μ m	Length x Width x Thickness : 148 x 3.75 x 0.08mm	No. : WS37EF		Single side Coated
	45 μ m	Length x Width x Thickness : 148 x 3.75 x 0.16mm	No. : SDS37		Double side Coated
	45 μ m	Length x Width x Thickness : 148 x 3.75 x 0.13mm	No. : DS37		Single side Coated
	30 μ m	Length x Width x Thickness : 148 x 3.75 x 0.12mm	No. : SDS37F		Double side Coated
	30 μ m	Length x Width x Thickness : 148 x 3.75 x 0.10mm	No. : DS37F		Single side Coated
	15 μ m	Length x Width x Thickness : 148 x 3.75 x 0.10mm	No. : SDS37EF		Double side Coated
	15 μ m	Length x Width x Thickness : 148 x 3.75 x 0.08mm	No. : DS37EF		Single side Coated



10 pcs of each per pack
or 10 pcs mixed 2 • + 4 • + 4 •



Tungsten Carbide Bur

钨钢车针



Advantages

1. Strictly comply with ISO Standard
2. Sharper and More Efficient
3. Durable and Valuable
4. FG, FG Short, FG Long, RA, RA Long etc



Ball, Plain Cut

ISO 500 ... 001 001	005	006	008	010	012	014	016	018	021	023	025
FG 314	FG 1/4	FG 1/2	FG 1	FG 2	FG 3	FG 4	FG 5	FG 6	FG 7	FG 8	FG 9
FG S 313			FG S1	FG S2	FG S3	FG S4	FG S5				
FG X 316			FG X1	FG X2	FG X3	FG X4	FG X5	FG X6	FG X7	FG X8	
RA 204		RA 1/2	RA 1	RA 2	RA 3	RA 4	RA 5	RA 6	RA 7	RA 8	
RA X 205				RA X2	RA X3	RA X4	RA X5	RA X6	RA X7	RA X8	



Inverted Cone

ISO 500 ... 010 001	008	010	012	014	016	018	021	023
FG 314	FG 34	FG 35	FG 36	FG 37	FG 38	FG 39		
RA 204	RA 34	RA 35	RA 36	RA 37	RA 38	RA 39	RA 40	RA 41



Cylinder, Plain Cut

ISO 500 ... 107 006	008	009	010	012	014	016			
Working Length/mm	4.0	4.0	4.0	4.0	4.0	4.0			
FG 314	FG 55	FG 56	FG 57	FG 58	FG 59	FG 60			
FG X 316			FG X57	FG X58					



Flat End Taper,
Plain Cut

ISO 500 ... 168 006	008	009	010	012	016	009	010	012	016
Working Length/mm	4.0	4.0	4.0	4.0	4.0	6.0	6.0	6.0	6.0
FG 314	FG 168	FG 169	FG 170	FG 171	FG 172	FG 169L	FG 170L	FG 171L	FG 172L
FG X 316			FG X170	FG X171					



Pear Plain Cut




ISO 500 ... 233 001	009	012	006	008	010	012	010	012
Working Length/mm	2.5	2.9	1.7	1.8	2.0	2.2	4.1	4.1
FG 314	FG 245	FG 246	FG 329	FG 330	FG 331	FG 332	FG 331L	FG 332L
FG S 313	FG S245							
RA 204	RA 245							



Cylinder, Cross Cut



ISO 500 ... 107 007	008	009	010	012	014	016	018	021	023
Working Length/mm	4.0	4.0	4.0	4.0	4.0	4.0	4.0	4.0	4.0
FG 314	FG 555	FG 556	FG 557	FG 558	FG 559	FG 560			
FG S 313		FG S556	FG S557						
FG X 316			FG X557	FG X558	FG X559				
RA 204		RA 556	RA 557	RA 558	RA 559	RA 560	RA 561	RA 562	RA 563
RA X 205			RA X557						



Cylinder, Cross Cut



ISO 500 ... 110 007	009	010	012	014
Working Length/mm	6.0	6.0	6.0	6.0
FG 314	FG 556L	FG 557L	FG 558L	FG 559L



Flat End Taper,
Cross Cut



ISO 500 ... 168 007	008	009	010	012	016	021	023	009	010	012	016
Working Length/mm	4.0	4.0	4.0	4.0	4.0	4.0	4.0	6.0	6.0	6.0	6.0
FG 314	FG 698	FG 699	FG 700	FG 701	FG 702			FG 699L	FG 700L	FG 701L	FG 702L
FG X 316			FG X700	FG X701	FG X702						
RA 204		RA 669	RA 700	RA 701	RA 702	RA 703	RA 704				
RA X 205				RA X701							



							
ISO 500 ... 137 007	008	009	010	012	014	016	
Working Length/mm	4.0	4.0	4.0	4.0	4.0	4.0	
FG 314	FG 1155	FG 1156	FG 1157	FG 1158	FG 1159	FG 1160	
FG S 313	FG S1157						

							
ISO 500 ... 194 006	009	010	012	016	ISO 500 ... 194 007	012	016
Working Length/mm	4.0	4.0	4.0	4.0	Working Length	4.0	4.0
FG 314	FG 1169	FG 1170	FG 1171	FG 1172	FG 314	FG 1701	FG 1702




						
ISO 500 ... 137 007	009	010	012	014		
Working Length/mm	4.0	4.0	4.0	4.0		
FG 314	FG 1556	FG 1557	FG 1558	FG 1559		
FG X 316	FG X1557	FG X1558				

							
ISO 500 ...	009	012	ISO 500 ... 468 ...	014	ISO 500 ... 150 001	010	012
Working Length/mm	3.0	5.5	Working Length/mm	4.0	Working Length/mm	0.2	0.2
FG 314	ET3	ET5	FG 314	FG 468	FG 314	FG 957	FG 958



								
ISO 500 ... 237.008	010	012	ISO 500 ... 137.008	010	012	ISO 500 ... 107 008	010	012
Working Length/mm	2.0	2.2	Working Length/mm	4.0	4.0	Working Length/mm	4.0	4.0
FG 314	FG 1931	FG 1932		FG 1957	FG 1958		FG 2057	FG 2058





Trimming & Finishing Carbide Burs




White Stone 白车石 Finish,Polish Filling Composite 去树脂,填充物及抛光。



Green Stone 绿车石 Grind Ceramics,Precious Metal 研磨烤瓷,贵金属等。






Diamond Stones 陶瓷金刚车石 Finish Zirconia, Lithium Disilicate, Ceramics . 研磨氧化锆、玻璃陶瓷、烤瓷。



Diamond Rubber Polishers 橡胶金刚石抛光磨头

Universally used to polish Zirconia, Lithium Disilicate, Ceramics, Composite etc.
广泛的用于抛光氧化锆、玻璃陶瓷、烤瓷、树脂等。



Dental Restorative Finisher & Polisher

牙科修复体修整与抛光磨头

For contouring, finishing and polishing cured composite restorations, and for polishing enamel surfaces after dental cleaning.

Aluminum Oxide Restorative Finisher

For contouring, finishing and smoothing cured composite restorations.

氧化铝复合树脂修整磨头

用于已固化复合树脂修复体的轮廓修整、表面精修和平滑。

Al₂O₃ • Polyurethane • Medium 200#
10 pieces / box



REF-100WB

10 × 1 mm



REF-060CB

6 × 7 mm



REF-050PB

5 × 11 mm

Diamond Restorative & Enamel Polisher - Medium

For polishing cured composite restorations and enamel surfaces after dental cleaning.

金刚石复合树脂及牙釉质中粗抛光磨头

用于已固化复合树脂修复体抛光及洁牙后的牙釉质表面抛光。

Diamond • Resin bond • Medium
5 pieces / box



REP-100DB

10 × 2 mm



REP-055TB

5.5 × 9.5 mm

Diamond Restorative & Enamel Polisher - Fine

For fine polishing of cured composite restorations and enamel surfaces after dental cleaning.

金刚石复合树脂及牙釉质精细抛光磨头

用于已固化复合树脂修复体精细抛光及洁牙后的牙釉质表面抛光。

Diamond • Resin bond • Fine
5 pieces / box



REP-100DW

10 × 2 mm



REP-055TW

5.5 × 9.5 mm

Cleaning Burs after debonding (脱胶后)洁牙磨头

Remove Calculus or Stains , No damage to the enamel.

 3 pieces mixed



2.3
9.0
Coarse



2.3
7.0



2.3
7.5



2.3
9.0
Medium



2.3
7.0



2.3
7.5

PC00C

PC02C

PC03C

PC00

PC02

PC03

 4 Medium



3.0
7.5



5.0
13.0



10.0
0.7



6.0
9.0

BR101

BR103

BR104

BR106

Silicon
矽粒子

 25



6x10 mm



10x0.7 mm



6x9 mm



SM103

SM104

SM106

53RS-2

Twin-Color Polishers for Composite

双色·树脂专用抛光磨头



Composite, Acrylic teeth's Polishing,
Polishing after remove Calculus or Stains , No damage to the enamel.
Polishing after remove the debonding

Item No.  12



Tr103C



Tr104C



Tr105C



Tr106C

Thickness (mm)

5 x 10.5mm

10 x 1 mm

10 x 2.0 (0.7) mm

6 x 8 mm

Item No.  12



Tr103M



Tr104M



Tr105M



Tr106M

Thickness (mm)

5 x 10.5mm

10 x 1 mm

10 x 2.0 (0.7) mm

6 x 8 mm

Item No.  12



Tr103F



Tr104F



Tr105F



Tr106F

Thickness (mm)

5 x 10.5mm

10 x 1 mm

10 x 2.0 (0.7) mm

6 x 8 mm

Silicon Polisher

For Ceramics

硅胶磨头
烤瓷抛光



White , Coarse

Cf101C

Cf103C

Cr101C

Cr103C

Cr104C

Cr105C

Cr106C

Cr108C

Silicon Polisher

For Ceramics

硅胶磨头
烤瓷抛光



Pink , Medium

Cf101M

Cf103M

Cr101M

Cr103M

Cr104M

Cr105M

Cr106M

Cr108M

Silicon Polisher

For Ceramics, Composite

硅胶磨头
烤瓷、树脂抛光



Blue , Fine

Cf101F

Cf103F

Cr101F

Cr103F

Cr104F

Cr105F

Cr106F

Cr108F

Silicon Polisher

For Metal, Alloy

硅胶磨头
金属、合金抛光



Black , Coarse

Af101C

Af103C

Ar101C

Ar103C

Ar104C

Ar105C

Ar106C

Ar108C

Silicon Polisher

For Metal, Alloy

硅胶磨头
金属、合金抛光



Brown , Medium

Af101M

Af103M

Ar101M

Ar103M

Ar104M

Ar105M

Ar106M

Ar108M

Silicon Polisher

For Metal, Alloy,
Composite

硅胶磨头
金属、合金、树脂抛光



Fine , Light Grey

Af101F

Af103F

Ar101F

Ar103F

Ar104F

Ar105F

Ar106F

Ar108F

Silicon Polisher

For Composite

硅胶磨头
树脂抛光



Fine , Yellow

Rf101F

Rf103F

Rr101F

Rr103F

Rr104F

Rr105F

Rr106F

Rr108F

Polishing Brushes

Nature Material, No hurt to teeth
Final polishing composite,
nature teeth or Ceramics
天然毛质、不伤牙
适用于树脂、真牙、烤瓷等



Felt
毛毡

R1Ns



Leather
羊皮

R1Ry



Goat Hair
羊毛

R1G



Thread
白棉线

R1T

Polishing Brushes

Abrasive Fiber
磨料丝 (含碳化硅
或氧化铝)

For Composite,
Ceramics

With Al Oxide 氧化铝



R1Al



R4Al



R5Al



R7Al

With Silicon Carbide 碳化硅



R1SiC



R4SiC



R5SiC



R7SiC



White Bristle

R4Dw



White Nylon

R4Cw



Clear Nylon

R4Cc



Blue Nylon

R4Cb



Pink Nylon

R4Cr



Purple Nylon

R4Cp



Yellow Nylon

R4Cy




White Goat Hair

R4G



White Bristle

R7Dw



Clear Nylon

R7Cc



White Goat Hair

R7G



White Bristle

R5Dw



White Nylon

R5Cw



Clear Nylon

R5Cc



Blue Nylon

R5Cb



Pink Nylon

R5Cr



Purple Nylon

R5Cp



Yellow Nylon

R5cy



White Goat Hair

R5G



4 Angles, Soft

R01W



6 Angles, Hard

R02W

Mandrels

夹 针



31RS



31RSL



35RS



52RS



53RS



26RS



51R



No. 201
for FG burs * 15



No. 211
for 10*FG + 5*RA



No. 221
for 5*FG + 10*RA



No. 241
for RA burs * 12



No. 202
for Root canal files * 18



No. 212
For HP burs * 18



No. 015A
for 18*FG



No. 015C
for 12*FG + 6*RA



No. 015E
for 18*RA



No. 015D
for 20*FG + 10*RA



No. 013E
for 16*FG + 9*RA



No. 022D
for 18*RA



No. 016B
for 8*HP



No. 023
132 * 107 * 65 mm
for FG, RA, HP, Files Mixed

Plastic Blocks, Sterilisable

塑料消毒盒

KIT 套装

Provide kinds of kits for clinic use , which offer the smartest operation ;
Whatever for Ceramics, Zirconia , Lithium Disilicate, Composite,Metal or Alloy ;
From Preparation to Final polishing ;
High Efficiency to shining, No need Glaze after polishing .



BR1901

BR-49, BR-45, BR-46, BR-41, BR-40,
BR-31, BR-28, BR-26, BR-25, BR-23 .

Diamond burs of Ball shape Assorted in FG shank
圆球状金刚石车针混装(高速)



BR1902

BR-X45, BR-X46, BR-X41, BR-X40, BR-X31,
BR-X40S, BR-X31S, BR-X41C, BR-X40C, BR-X31C


Diamond burs of Ball shape Assorted in FG LONG
圆球状金刚石车针混装(高速加长)



BR1903

BR-R46, BR-R41, BR-R40, BR-R31, BR-R26,
BR-R40F, BR-R31F, BR-R26F, BR-R40EF, BR-R26EF


Diamond burs of Ball shape Assorted in RA shank
圆球状金刚石车针混装(低速)



4100

508T-014S, BR-41, 198-016S, 508T-016S, 508T-016F ;
508T-014S, BR-540, 198-021S, 508T-021S, 508T-021F


Safe End Diamond Preparation
前后牙预备(安全端金刚石车针)



4101

TF-21, TF-10, TF-12, TC-10, WR-13
SR-11, TR-11, EX-18, TR-12F, TC-11F

Crown & Bridge Preparation for Anterior
前牙桥体/冠体预备



4102

EX-20S, SR-21, FO-34, TR-62C, TR-21
TR-24, SO-21, EX-18, TR-21F, EX-21F


Crown & Bridge Preparation for Posterior
后牙桥体/冠体预备



18041V

SI-45 x 2 ; SI-46 x 3,
SI-47 x 3, SI-48 x 2 .

Diamond burs of Inverted Cone Assorted in FG shank
常用倒锥金刚石车针混装(高速)



1805TA

TR-11 x 2, TR-12 x 2, TF-11 x 2 ;
TF-12, TC-10, TC-27, TC-28 x 1

Regular Diamond burs of Taper Assorted in FG shank
常用锥形金刚石车针混装 (高速)



1806

TC-11EF, 820-016EF, 820-031F,
820-016F, TC-10F

Interproximal Enamel Reduction (IPR) for Ortho
邻面去釉套装



4103

BR-45, BR-40, BC-32, EX-31
TF-20, SF-41, SF-21, SI-48

Cavity Preparation, Restoratives
备洞



4104

SI-48, DI-42, TF-31, TF-22, SF-21
EX-18, TF-21F, SF-21F, FO-22EF, EX-21EF

Inlay & Onlay Preparation
嵌体/高嵌体预备



4105

BR-49, SF-41, SR-41, TC-26, CR-41
EX-41, FO-40EF / FG1/4, FG1/2, FG 1

Micro Invasive Surgery
微创备牙



6000

BR-45, VP-20, VP-21, VP-22,
SR-41, TR-25, SO-20, TR-13F, CF101F

Porcelain Veneers' Preparation
瓷贴面备牙套装



6002

CD-50F, CD-51F, CD-52F, CD-53F, CD-54F
CD-55F, CD-56F, CD-57F, CD-58F, CD-59F

Children's Teeth Preparation
儿童牙预备



4106

BR-40EF, TC-11EF, TF-11EF, TF-12EF, TR-11EF,
EX-21EF, FO-40EF, FO-21EF, FO-22EF, FO-30EF

Ceramics / Composite Polishing
烤瓷牙/复合树脂抛光



FG008D

BR-49, BR-49, SI-49, SF-49, TF-49

Student or Beginner's Diamond FG 0.8mm
新手入门金刚石车针, 高速



FG009D

001-009M, 001-009M, 010-009M,
109-009M, 170-009M

Student or Beginner's Diamond FG 0.9mm
新手入门金刚石车针, 高速



RA008D

BR-R49, BR-R49, SI-R49, SF-R49, TF-R49

Student or Beginner's Diamond RA 0.8mm
新手入门金刚石车针, 慢速



RA009D

001-R009M, 001-R009M, 010-R009M,
109-R009M, 170-R009M


Student or Beginner's Diamond RA 0.9mm
新手入门金刚石车针, 慢速



TF7100

7004, 7404, 7214, 7714, 7903L

TC Trimming & Finishing burs, FG
钨钢精修车针, FG



1900 T

FG 1, FG 2, FG 3, FG 4, FG 5


Tungsten carbide burs of Ball shape Assorted in FG shank
圆球状钨钢车针混装(高速)



FG008T

FG 1, FG 1, FG 34, FG555, FG698

Tungsten Carbide of FG shank 0.8mm for Students
初学者用的钨钢车针 (头径0.8mm, 高速)



RA008T

RA 1, RA 1, RA 34, RA 555, RA 698

Tungsten Carbide of RA shank 0.8mm for Students
初学者用的钨钢车针 (头径0.8mm, 低速)



7109FW

FL2-Fw, CN1-Fw, RD1-Fw, RD2-Fw, TC1-Fw,

FG White Stone Burs, to Finish Composite
白车石, 高速



7109RW

FL2-Rw, CN1-Rw, RD1-Rw, RD2-Rw, TC1-Rw,

RA White Stone Burs, to Finish Composite
白车石, 低速



7109FG

FL2-Fg, CN1-Fg, RD1-Fg, RD2-Fg, TC1-Fg,

FG Green Stone burs, Finish Ceramics
绿车石, 高速



7109RG

FL2-Rg, CN1-Rg, RD1-Rg, RD2-Rg, TC1-Rg,

RA Green Stone burs, Finish Ceramics
绿车石, 低速



1906

Cr103C x 2, Cr104C, Cr106C, Cr103M, Cr104M, R4AL, R5AL; R4Cb, R5Cb, R4Cr, R5Cr

Composite Finish and Polish RA
树脂成型和抛光,弯机柄



1907

Cr103C, Cr104C, Cr106C, Cr103M, Cr104M, Cr106M, R4AL, R1AL, R5AL; R1Ry, R1G, R1T

Composite Finish and Polish RA
树脂成型和抛光,弯机柄



QR2101

2 step: Sr003F: Pro-Polishing
Sr003UF: Polishing

Composite Polish, RA
树脂抛光,弯机柄



QR2102

2 step: Sr304F: Pro-Polishing
Sr304UF, Sr003UF: Polishing

Composite Polish, RA
树脂抛光,弯机柄



QR2103

3 step: Sr304CM: Smooth;
Sr003F: Pro-Polishing; Sr003UF: Polishing

Composite Polish, RA
树脂抛光,弯机柄



QR2104

3 step: Sr003M: Smooth;
Sr003F: Pro-Polishing; Sr003UF x 2: Polishing

Composite Polish, RA
树脂抛光,弯机柄



CP 6011

2 Step: Sr 204 M x 3: Pro-polish;
Sr 204 F x 3: Polish

Composite, Ceramic, Zirconia Polishers . RA
树脂、烤瓷、氧化锆抛光



CP 6001

3 Step: Sr 204 C x 2: Smooth;
Sr 204 M x 2: Pro-polish; Fine: Sr 204 F x 2: Polish

Composite, Ceramic, Zirconia Polishers . RA
树脂、烤瓷、氧化锆抛光



7115

Smooth: DCR101, DCR102, DCR122, DCR105
Pro-polish: Sr101F, Sr103F, Sr105F, Sr106F
Polish: Sr101UF, Sr103UF, Sr105UF, Sr106UF


Diamond stones & polishers for Zirconia / Lithium
Disilicate, 氧化锆, 玻璃陶瓷研磨抛光



1904

X151, X1151, EC25, EC28, FG X6,
EX-28, BR-X31, TR-X12, TR-X15, TF-X13

Extraction Burs
拔阻生齿车针



1905

FG X2, FG X4, FG X6, FG X558, FG X701,
BR-X41, BR-X40, BR-X31, TF-X12, TR-X12

Oral Surgery burs FG Long
手术车针(高速加长)



4110V

Smooth: Rw01, Rw02 x 2
Pro-polish: Cr101M, Cr103M, Cr104M, Cr106M,
Polish: Cr101F, Cr103F, Cr104F, Cr106F

Composite Finish and Polish RA
树脂成型和抛光,弯机柄



4112

Smooth: Cr101C, Cr103C, Cr104C, Cr106C
Pro-polish: Cr101M, Cr103M, Cr104M, Cr106M
Polish: Cr101F, Cr103F, Cr104F, Cr106F

Porcelain / Natural Teeth Polishing RA
烤瓷牙, 自然牙抛光, 弯机柄



4114

Smooth: Ar101C, Ar103C, Ar104C, Ar106C
Pro-polish: Ar101M, Ar103M, Ar104M, Ar106M
Polish: Ar101F, Ar103F, Ar104F, Ar106F

Amalgam / Metal Polishing RA
银汞合金, 金属抛光, 弯机柄



4115

Smooth: Sr103M, Sr105M, Sr106M; Pro-polish: Sr103F, Sr105F, Sr106F
Polish: Sr101EF, Sr103EF, Sr105EF, Sr106EF,
Final Polish: Sr101UF, Sr103UF

Diamond Polishers for Zirconia / Lithium Disilicate
氧化锆, 玻璃陶瓷抛光, 弯机柄

Crown & Bridge Preparation for Anterior
前牙桥体/冠体预备



4101A

TF-21, TF-10, TF-12 x 2, WR-13, TC-21,
TR-11, TC-10/EX-18, SR-11/ TR-12F, TC-11F




Crown & Bridge Preparation for Posterior
后牙桥体/冠体预备



4102A

EX-11, EX-20S, TR-22S, BR-31C, TR-62C, FO-34/
TR-24, TR-21/EX-18, SO-21/EX-21F, TR-21F




Cavity Preparation, Restoratives
备洞



4103A

BR-45, BR-40, BR-46, BR-31, BR-41/BC-32,
SI-47, SI-48/SF-41, SF-21/EX-31, TF-20




Inlay/Onlay Preparation
嵌体/高嵌体预备



4104A

DI-42, SI-48, TF-22, TF-31 / TF-21F, SF-21F /
TF-21, SF-21, EX-18, BC-32 / EX-21EF, FO-22EF




Porcelain Veneer's Preparation
瓷贴面备牙




6000A

SR-41, BR-45 / VP-20 x 2, VP-21 x 2, VP-22 x 2/
SO-20, TR-25, TR-13F, TC-28EF



Trimming & Finishing Tungsten
Carbide burs FG 钨钢精修车针



6003A

7004, 7005, 7006, 7404, 7406, 7114,
7214, 7714, 7901, 7902, 7903, 7903L




Zirconia, Lithium Disilicate Polishing kit, RA
氧化锆, 玻璃陶瓷抛光(低速)




7115A

Smooth : DCR101 , DCR102, DCR122, DCR105 ; Pro-polish: Sr101F, Sr103F, Sr105F, Sr106F
Polish: Sr101UF, Sr103UF, Sr105UF, Sr106UF




Zirconia, Lithium Disilicate High Gloss polishing kit, RA
氧化锆抛光(低速)



4115A

Smooth : Sr103M, Sr105M, Sr106M ; Pro-polish : Sr103F, Sr105F, Sr106F
Polish: Sr101EF, Sr103EF, Sr105EF, Sr106EF, Final Polish: Sr101UF, Sr103UF



Composite or filling material Polisher RA
美学树脂、填充材料抛光



CP6001A

3 Step: Sr 204 C x 2 : Smooth;
Sr 204 M x 2 : Pro-polish ; Fine : Sr 204 F x 2 : Polish



Composite Polishing kit , RA
树脂抛光套装(低速)



4110A2

Smooth: FL2-Rwx 2 , CN1-Rw x 2 ; Pro-polish: Cr101M, Cr103M, Cr104M, Cr106M,
Polish: Cr101F, Cr103F, Cr104F, Cr106F



Ceramics / Porcelain Teeth Polish , RA
烤瓷牙抛光(低速)



4112A

Smooth: Cr101C, Cr103C, Cr104C, Cr106C
Pro-polish : Cr101M, Cr103M, Cr104M, Cr106M
Polish: Cr101F, Cr103F, Cr104F, Cr106F



Amalgam, Alloy, Metal Polish , RA
银汞合金抛光(低速)




4114A

Smooth : Ar101C, Ar103C, Ar104C, Ar106C
Pro-polish : Ar101M, Ar103M, Ar104M, Ar106M
Polish: Ar101F, Ar103F, Ar104F, Ar106F



Micro Invasive Surgery Burs
微创手术用车针



4105A

BR-49 x 2, SR-41 x 2, SF-41, EX-41, TC-26,
CR-41, FG1/4, FG1/2, FG 1, FO-40EF

Extra Fine Diamond burs for polishing
特细粒度金刚砂抛光车针,高速



4106A

CR-11EF, BR-40EF, TC-27EF, TC-11EF, TF-12EF, TF-11EF,
FO-40EF, TR-11EF, FO-22EF, FO-21EF, FO-30EF, EX-21EF

Composite, Ceramics Finishing kit, FG
树脂成型(高速)



0302A

FL2-Fw x 3 ; CN1-Fw x 3
FL2-Fg x 3 ; CN1-Fg x 3

Composite, Ceramics Finishing kit, RA
树脂成型(低速)



0306A

FL2-Rw x 3 ; CN1-Rw x 3
FL2-Rg x 3 ; CN1-Rg x 3

White Stone Burs kit, FG
白车石(高速)



302F1

FL2-Fw x 2, CN1-Fw x 2, RD1-Fw x 2
RD2-Fw x 2, TC1-Fw x 2, IC2-Fw x 2

White Stone Burs kit, RA
白车石(低速)



306R1

FL2-Rw x 2, FL3-Rw x 2, CN1-Rw x 2,
RD1-Rw x 2, RD2-Rw x 2, TC1-Rw x 2.

HP Kit to adjust Zirconia
氧化锆修整椅旁套装(HP)



Z31

ZF33, DY001EF to Smoothing
D203F for Pre-polishing ; D203EF for Polishing

HP Kit to adjust Lithium Disilicate
玻璃陶瓷修整椅旁套装(HP)



L31

ZF33, DY001EF to Smoothing
D203F for Pre-polishing ; D203UF for Polishing

Samples Collection For Clinic Items

样品册



Flip-Cart Plastic Pages ;
Smart, Easily fill all burs ;
Economic and Recyclable ;
Nice-looking for showing to Customers;
Convenient to take to anywhere .



No. 170B5

Empty Collection Wallet
样品本(不带车针)



4 Colors Available ,
100 holes per page,
5 pages per set,
Smart to FG or RA Burs



No. 1701D

Diamond Burs Collection Wallet
金刚石车针样品本



More than 300 specification ,
Diamond burs ,
FG, FG XL, FG 2XL,
FG Short, RA, RA Long



No. 1702T

Tungsten Carbide Collection Wallet
钨钢车针样品本



Nearly 200 specification ,
TC burs
FG, FG XL, FG Short, RA, RA Long
T & F Gold, Endo z etc.



No. 1703P

Polishing burs Collection Wallet
抛光器械样品本



Nearly 150 specifications ,
Silicon Rubber burs,
Diamond Abrasives & polishers,
Polishing brushes etc




Four Colors Available

Brown




Black

Blue

Yellow



BurDental 博登科技



博登科技有限公司

B&D Technologies Co.,Ltd

地址: 广东省广州市番禺区洛浦街恒达产业园A-522

Address : A-522, Heng Da Industrial Park, Luo Pu Street, Pan Yu District, 510000,
GuangZhou, GuangDong Province, China

网站(Web): www.burdental.com ; 邮箱(Email): info@burdental.com

手机(Mobile, WhatsApp) : 0086 1866 6599 173