

KIT 套装

Provide kinds of kits for clinic use , which offer the smartest operation ;
 Whatever for Ceramics, Zirconia , Lithium Disilicate, Composite, Metal or Alloy ;
 From Preparation to Final polishing ;
 High Efficiency to shining, No need Glaze after polishing .



BR1901

BR-49, BR-45, BR-46, BR-41, BR-40,
 BR-31, BR-28, BR-26, BR-25, BR-23 .

Diamond burs of Ball shape Assorted in FG shank
 圆球状金刚石车针混装(高速)



BR1902

BR-X45, BR-X46, BR-X41, BR-X40, BR-X31,
 BR-X40S, BR-X31S, BR-X41C, BR-X40C, BR-X31C


Diamond burs of Ball shape Assorted in FG LONG
 圆球状金刚石车针混装(高速加长)



BR1903

BR-R46, BR-R41, BR-R40, BR-R31, BR-R26,
 BR-R40F, BR-R31F, BR-R26F, BR-R40EF, BR-R26EF

Diamond burs of Ball shape Assorted in RA shank
 圆球状金刚石车针混装(低速)



4100

508T-014S, BR-41, 198-016S, 508T-016S, 508T-016F ;
 508T-014S, BR-540, 198-021S, 508T-021S, 508T-021F


Safe End Diamond Preparation
 前后牙预备(安全端金刚石车针)



4101

TF-21, TF-10, TF-12, TC-10, WR-13
 SR-11, TR-11, EX-18, TR-12F, TC-11F

Crown & Bridge Preparation for Anterior
 前牙桥体/冠体预备



4102

EX-20S, SR-21, FO-34, TR-62C, TR-21
 TR-24, SO-21, EX-18, TR-21F, EX-21F


Crown & Bridge Preparation for Posterior
 后牙桥体/冠体预备



18041V

SI-45 x 2 ; SI-46 x 3,
 SI-47 x 3, SI-48 x 2 .

Diamond burs of Inverted Cone Assorted in FG shank
 常用倒锥金刚石车针混装(高速)



1805TA

TR-11 x 2, TR-12 x 2, TF-11 x 2 ;
 TF-12, TC-10, TC-27, TC-28 x 1


Regular Diamond burs of Taper Assorted in FG shank
 常用锥形金刚石车针混装 (高速)



1806

TC-11EF, 820-016EF, 820-031F,
 820-016F, TC-10F

Interproximal Enamel Reduction (IPR) for Ortho
 邻面去釉套装



4103

BR-45, BR-40, BC-32, EX-31
 TF-20, SF-41, SF-21, SI-48

Cavity Preparation, Restoratives
 备洞



4104

SI-48, DI-42, TF-31, TF-22, SF-21
 EX-18, TF-21F, SF-21F, FO-22EF, EX-21EF

Inlay & Onlay Preparation
 嵌体/高嵌体预备



4105

BR-49, SF-41, SR-41, TC-26, CR-41
 EX-41, FO-40EF / FG1/4, FG1/2, FG 1

Micro Invasive Surgery
 微创备牙



6000

BR-45, VP-20, VP-21, VP-22,
SR-41, TR-25, SO-20, TR-13F, CF101F

Porcelain Veneers' Preparation
瓷贴面备牙套装



6002

CD-50F, CD-51F, CD-52F, CD-53F, CD-54F
CD-55F, CD-56F, CD-57F, CD-58F, CD-59F

Children's Teeth Preparation
儿童牙预备



4106

BR-40EF, TC-11EF, TF-11EF, TF-12EF, TR-11EF,
EX-21EF, FO-40EF, FO-21EF, FO-22EF, FO-30EF


Ceramics / Composite Polishing
烤瓷牙/复合树脂抛光



FG008D

BR-49, BR-49, SI-49, SF-49, TF-49

Student or Beginner's Diamond FG 0.8mm
新手入门金刚石车针,高速



FG009D

001-009M, 001-009M, 010-009M,
109-009M, 170-009M

Student or Beginner's Diamond FG 0.9mm
新手入门金刚石车针,高速



RA008D

BR-R49, BR-R49, SI-R49, SF-R49, TF-R49

Student or Beginner's Diamond RA 0.8mm
新手入门金刚石车针,慢速



RA009D

001-R009M, 001-R009M, 010-R009M,
109-R009M, 170-R009M


Student or Beginner's Diamond RA 0.9mm
新手入门金刚石车针,慢速



TF7100

7004, 7404, 7214, 7714, 7903L

TC Trimming & Finishing burs, F6
钨钢精修车针, FG



1900 T

FG 1, FG 2, FG 3, FG 4, FG 5

Tungsten carbide burs of Ball shape Assorted in FG shank
圆球状钨钢车针混装(高速)



FG008T

FG 1, FG 1, FG 34, FG555, FG698

Tungsten Carbide of FG shank 0.8mm for Students
初学者用的钨钢车针 (头径0.8mm, 高速)



RA008T

RA 1, RA 1, RA 34, RA 555, RA 698

Tungsten Carbide of RA shank 0.8mm for Students
初学者用的钨钢车针 (头径0.8mm, 低速)



7109FW

FL2-Fw, CN1-Fw, RD1-Fw, RD2-Fw, TC1-Fw,

FG White Stone Burs, to Finish Composite
白车石, 高速



7109RW

FL2-Rw, CN1-Rw, RD1-Rw, RD2-Rw, TC1-Rw,

RA White Stone Burs, to Finish Composite
白车石, 低速



7109FG

FL2-Fg, CN1-Fg, RD1-Fg, RD2-Fg, TC1-Fg,

FG Green Stone burs, Finish Ceramics
绿车石, 高速



7109RG

FL2-Rg, CN1-Rg, RD1-Rg, RD2-Rg, TC1-Rg,

RA Green Stone burs, Finish Ceramics
绿车石, 低速



1906

Cr103C,Cr 104C, Cr106C;Cr103M, Cr104M, Cr106M, R4AL, R1AL,R5AL;R1Ry, R1G, R1T

Composite Finish and Polish RA
树脂成型和抛光,弯机柄



1907

Cr103C x 2,Cr 104C, Cr106C; , Cr103M, Cr104M; R4AL, R5AL; R4Cb, R5Cb, R4Cr, R5Cr

Composite Finish and Polish RA
树脂成型和抛光,弯机柄



QR2101

2 step : Sr003F : Pro-Polishing
Sr003UF: Polishing

Composite Polish, RA
树脂抛光,弯机柄



QR2102

2 step : Sr304F : Pro-Polishing
Sr304UF, Sr003UF: Polishing

Composite Polish, RA
树脂抛光,弯机柄



QR2103

3 step : Sr304CM: Smooth;
Sr003F : Pro-Polishing; Sr003UF: Polishing

Composite Polish, RA
树脂抛光,弯机柄



QR2104

3 step : Sr003M:Smooth;
Sr003F: Pro-Polishing; Sr003UF x 2: Polishing

Composite Polish, RA
树脂抛光,弯机柄



CP 6011

2 Step : Sr 204 M x 3: Pro-polish;
Sr 204 F x 3: Polish

Composite, Ceramic, Zirconia Polishers . RA
树脂、烤瓷、氧化锆抛光



CP 6001

3 Step: Sr 204 C x 2 : Smooth;
Sr 204 M x 2 : Pro-polish; Fine : Sr 204 F x 2: Polish

Composite, Ceramic, Zirconia Polishers . RA
树脂、烤瓷、氧化锆抛光



7115

Smooth : DCR101 , DCR102, DCR122,DCR105
Pro-polish: Sr101F, Sr103F, Sr105F, Sr106F
Polish: Sr101UF, Sr103UF, Sr105UF, Sr106UF

Diamond stones & polishers for Zirconia / Lithium Disilicate, 氧化锆,玻璃陶瓷研磨抛光



1904

X151, XL151, EC25, EC28, FG X6,
EX-28, BR-X31, TR-X12, TR-X15, TF-X13

Extraction Burs
拔阻生齿车针



1905

FG X2, FG X4, FG X6, FG X558, FG X701,
BR-X41, BR-X40, BR-X31, TF-X12, TR-X12

Oral Surgery burs FG Long
手术车针(高速加长)



4110V

Smooth: Rv01, Rv02 x 2
Pro-polish: Cr101M, Cr103M, Cr104M, Cr106M,
Polish: Cr101F, Cr103F, Cr104F, Cr106F

Composite Finish and Polish RA
树脂成型和抛光,弯机柄



4112

Smooth: Cr101C, Cr103C, Cr104C, Cr106C
Pro-polish: Cr101M, Cr103M, Cr104M, Cr106M
Polish: Cr101F, Cr103F, Cr104F, Cr106F

Porcelain / Natural Teeth Polishing RA
烤瓷牙,自然牙抛光,弯机柄



4114

Smooth: Ar101C, Ar103C, Ar104C, Ar106C
Pro-polish: Ar101M, Ar103M, Ar104M, Ar106M
Polish: Ar101F, Ar103F, Ar104F, Ar106F

Amalgam / Metal Polishing RA
银汞合金,金属抛光,弯机柄



4115

Smooth : Sr103M, Sr105M, Sr106M ; Pro-polish :Sr103F, Sr105F, Sr106F
Polish: Sr101EF, Sr103EF, Sr105EF, Sr106EF,
Final Polish: Sr101UF, Sr103UF


Diamond Polishers for Zirconia /Lithium Disilicate
氧化锆,玻璃陶瓷抛光,弯机柄

Crown & Bridge Preparation for Anterior
前牙桥体/冠体预备



4101A

TF-21, TF-10, TF-12 x 2, WR-13, TC-21,
TR-11, TC-10/EX-18, SR-11/ TR-12F, TC-11F




Crown & Bridge Preparation for Posterior
后牙桥体/冠体预备



4102A

EX-11, EX-20S, TR-22S, BR-31C, TR-62C, FO-34/
TR-24, TR-21/ EX-18, SO-21/EX-21F, TR-21F




Cavity Preparation, Restoratives
备洞



4103A

BR-45, BR-40, BR-46, BR-31, BR-41/ BC-32,
SI-47, SI-48/ SF-41, SF-21/ EX-31, TF-20



Inlay/Onlay Preparation
嵌体/高嵌体预备



4104A

DI-42, SI-48, TF-22, TF-31 / TF-21F, SF-21F /
TF-21, SF-21, EX-18, BC-32 / EX-21EF, FO-22EF



Porcelain Veneer's Preparation
瓷贴面备牙




6000A

SR-41, BR-45 / VP-20 x 2, VP-21 x 2, VP-22 x 2/
SO-20, TR-25, TR-13F, TC-28EF



Trimming & Finishing Tungsten
Carbide burs FG 钨钢精修车针



6003A

7004, 7005, 7006, 7404, 7406, 7114,
7214, 7714, 7901, 7902, 7903, 7903L



Zirconia, Lithium Disilicate Polishing kit, RA
氧化锆, 玻璃陶瓷抛光(低速)



7115A

Smooth : DCR101 , DCR102, DCR122, DCR105 ; Pro-polish: Sr101F, Sr103F, Sr105F, Sr106F
Polish: Sr101UF, Sr103UF, Sr105UF, Sr106UF



Zirconia, Lithium Disilicate High Gloss polishing kit, RA
氧化锆抛光(低速)



4115A

Smooth : Sr103M, Sr105M, Sr106M ; Pro-polish : Sr103F, Sr105F, Sr106F
Polish: Sr101EF, Sr103EF, Sr105EF, Sr106EF, Final Polish: Sr101UF, Sr103UF



Composite or filling material Polisher RA
美学树脂、填充材料抛光



CP6001A

3 Step: Sr 204 C x 2 : Smooth;
Sr 204 M x 2 : Pro-polish ; Fine : Sr 204 F x 2 : Polish



Composite Polishing kit , RA
树脂抛光套装(低速)



4110A2

Smooth: FL2-Rwx 2 , CN1-Rw x 2 ; Pro-polish: Cr101M, Cr103M, Cr104M, Cr106M,
Polish: Cr101F, Cr103F, Cr104F, Cr106F



Ceramics / Porcelain Teeth Polish , RA
烤瓷牙抛光(低速)



4112A

Smooth: Cr101C, Cr103C, Cr104C, Cr106C
Pro-polish : Cr101M, Cr103M, Cr104M, Cr106M
Polish: Cr101F, Cr103F, Cr104F, Cr106F



Amalgam, Alloy, Metal Polish , RA
银汞合金抛光(低速)




4114A

Smooth : Ar101C, Ar103C, Ar104C, Ar106C
Pro-polish : Ar101M, Ar103M, Ar104M, Ar106M
Polish : Ar101F, Ar103F, Ar104F, Ar106F



Micro Invasive Surgery Burs
微创手术用车针



4105A

BR-49 x 2, SR-41 x 2, SF-41, EX-41, TC-26,
CR-41, FG1/4, FG1/2, FG 1, FO-40EF

Extra Fine Diamond burs for polishing
特细粒度金刚砂抛光车针,高速



4106A

CR-11EF, BR-40EF, TC-27EF, TC-11EF, TF-12EF, TF-11EF,
FO-40EF, TR-11EF, FO-22EF, FO-21EF, FO-30EF, EX-21EF

Composite, Ceramics Finishing kit, FG
树脂成型(高速)



0302A

FL2-Fw x 3; CN1-Fw x 3
FL2-Fg x 3; CN1-Fg x 3

Composite, Ceramics Finishing kit, RA
树脂成型(低速)



0306A

FL2-Rw x 3; CN1-Rw x 3
FL2-Rg x 3; CN1-Rg x 3

White Stone Burs kit, FG
白车石(高速)



302F1

FL2-Fw x 2, CN1-Fw x 2, RD1-Fw x 2
RD2-Fw x 2, TC1-Fw x 2, IC2-Fw x 2

White Stone Burs kit, RA
白车石(低速)



306R1

FL2-Rw x 2, FL3-Rw x 2, CN1-Rw x 2,
RD1-Rw x 2, RD2-Rw x 2, TC1-Rw x 2.

HP Kit to adjust Zirconia
氧化锆修整椅旁套装(HP)



Z31

ZF33, DY001EF to Smoothing
D203F for Pre-polishing; D203EF for Polishing

HP Kit to adjust Lithium Disilicate
玻璃陶瓷修整椅旁套装(HP)



L31

ZF33, DY001EF to Smoothing
D203F for Pre-polishing; D203UF for Polishing